


Good-Lite Product Catalog

YOUR SOURCE ^{for} VISION



Cognitive Test • Grating Acuity • Wall Charts • Handheld Charts • EyeSpy 20/20
18 Charts • ETDRS Charts • Low Vision • Intermediate Vision • Near Vision • Re
Cards • ETDRS Cabinet • Spot Vision Screener • ESV1200 - 9x14 • ESV1018 - 10
nsta-Line • CSV-1000 • Super Pinhole • Accessories • Screening Software • Test
software • Fixation • Occluders • Fun Frames • Prisms • Hyperopia • Retinoscopy
Models • Stereo Test • Worth 4-Dot • Color Books & Tests • D15 Style Color Te
saturated Color Test • Projector Slides • Continuous Text • Adult Charts • Child
charts • Pediatric Charts • Vectographic Slides • Amsler Grid • Campimeter • Pris
ar • Loose Prism • Spanish • Vectograph • Maddox • Cylinders • Phoropter • Tri
s/Frames • Filters • Risley • Eye Patches • Light Numerical • Contr
sitivity • Visual Field • Cognitive Vision • Adaptation • Color Vision • Optokinet

GOOD-LITE



Celebrating More Than 80 Years of Good-Lite

The Good-Lite Company is celebrating more than 80 years as the industry leader in illuminated cabinets; evidence-based eye charts, such as the LEA Symbols, and many more vision screening and testing products available in this catalog and online.

The combined efforts of Robert Good, MD, Professor of Ear, Nose, and Throat at the Chicago College of Medicine and Surgery, and his son, Palmer Good, MD, ophthalmologist, created the Good-Lite Company in 1930 with the development of the first Good-Lite medical headlight.

Later, Palmer added vision testing equipment to the product line and led the Good-Lite Company to its honored position among today's vision screening and testing providers.

Good-Lite's numerous achievements over the last 80 years are noteworthy. Palmer Good worked closely with Louise Sloan, PhD, of the Wilmer Eye Institute at Johns Hopkins to develop the Sloan Letters for distance and near testing. Good-Lite is proud to have produced the first commercial chart using Sloan Letters, copyrighted in 1961 by Palmer Good. Sloan letters continue to be the optotype of choice for many medical professionals, and were the precursor to ETDRS charts in worldwide use in clinical trials.

In the 1970s, Palmer worked with Otto Lippmann, MD, to develop HOTV optotypes using Sloan letters. Lippmann's HOTV optotypes — a modified version of a Stycar matching letters test based on Snellen letters — are useful for screening and testing vision of children as young as age 3 years. The American Academy of Pediatrics, et al., includes HOTV optotypes on its list of recommended tests for preschool vision screening.

In 2002, Good-Lite proudly received exclusive, worldwide manufacturing rights to the complete LEA Test System developed by Lea Hyvärinen, MD, PhD, pediatric ophthalmologist, of Finland. Dr. Hyvärinen started the first vision rehabilitation service in Finland in 1976, designed the LEA Symbols in 1976, and first described her 4 symbols—apple, house, square, and circle—in a peer-reviewed journal in 1980. This System includes the LEA Symbols, the only pediatric optotypes that blur equally at threshold to prevent guessing, and the LEA Numbers. The core tests in this unique system assess near and distance visual acuity, contrast sensitivity, visual field deficits, cognitive vision, adaptation, and color deficiency. They are useful for screening and testing children of all ages, as well as for conducting functional vision assessments

in children and adults with developmental delays/disabilities. The American Academy of Pediatrics, et al., includes LEA Symbols on its list of recommended tests for preschool vision screening.

Good-Lite values its key role in helping vision screeners and eye care professionals detect and treat vision disorders. We take pride in meeting the diverse needs of all those who assess vision, including the many clinical trials investigators who look to Good-Lite for state-of-the-art, evidence-based, standardized equipment. As we celebrate our 80+ years, Good-Lite remains family owned and operated and is currently in the capable hands of the third and fourth generations of the Good family.

In celebration of more than 80 years of success, we thank you, our valued customers, for your continued confidence in us. We appreciate your business and you have our assurance that we will continue to provide you with dependable, cost-efficient products and professional services in the years to come.

Good-Lite's Vision Experts



Lea Hyvärinen, MD

Lea Hyvärinen, MD, PhD, Professor h.c., Rehabilitation Sciences, University of Dortmund, Germany; Senior Lecturer, Developmental Neuropsychology, University of Helsinki, Finland, started to develop her LEA Vision Test System in 1976 in Helsinki and has worked on it with Chris Greening in 17 years. The Test System now contains tests for numerous clinical test situations and vision screening of children and adults with different communication needs.

The original tests on the printed near vision card were: a standard test, a test with wide spacing and a tightly crowded test. There were two tests for distance, a hand-made grating acuity test and for low contrast, photographic tests at the same contrast levels as today. In early 1990s came the LEA Puzzle, first black & white, later colorful; Panel 16 color vision test; LEA GRATINGS; Cone Adaptation Test; Rectangle Game; Mailbox Game; Heidi Expressions Game and the LEA Numbers. The additions of the recent years are LEA Grating Acuity Test, LEA Flicker Wand, the 40M and 60M cards for presentation as a line test or singly and the LEA Low Contrast Grating Acuity Test.

Dr. Hyvärinen's contributions to vision include training courses and lectures on functional vision testing for children with disabilities. The next pages depict the Core Tests and their numerous variants that are still available for numerous screening projects and programs.



P. Kay Nottingham Chaplin, EdD Director, Vision and Eye Health Initiatives

P. Kay Nottingham Chaplin, EdD, helped Geoffrey E. Bradford, MD, pediatric ophthalmologist at West Virginia University (WVU) Eye Institute and member of the American Academy of Pediatrics' Section on Ophthalmology Executive Committee, to create the Vision Initiative for Children, a program that trained and equipped individuals to screen the vision of preschoolers.

Dr. Chaplin is a member of the national Expert Advisory Committee and Education Subcommittee to the National Center for Children's Vision and Eye Health, an initiative housed at Prevent Blindness America. Her role is to help implement recommendations of the National Expert Panel for a universal preschool vision screening strategy, to develop a vision screening training and certification program, and to contribute to stakeholder education.

Dr. Chaplin has lectured, trained, and consulted at more than 100 national, state, and local venues, including conferences for the:

- National Association of School Nurses,
- National Head Start Association,
- National Association of Pediatric Nurse Practitioners,
- Society for Physician Assistants in Pediatrics, and the
- National Assembly on School-Based Health Care.

GOOD-LITE

CATALOG CONTENTS

Follow **Lea Hyvärinen, MD, PhD**, online:



Dr. Lea Hyvarinen and
Lea-Test Ltd



www.drleahyvarinen.com



www.lea-test.fi/

lea SYMBOLS®

lea-numbers®



These products and
much, much more can be
found on our website at:

www.good-lite.com



LEA Core Vision Tests — Page 2

The Core Tests based on LEA Symbols and LEA Numbers designed for screening and assessment of visual functioning.

Color Vision Tests — Page 9

Panel 16	Magnetic Color Test	HMC Anomaloscope
Color Plate Books	HRR Color Tests	Farnsworth 100 Hue

Computerized Vision — Page 12

Vision Testing Software EyeSpy 20/20™	Vision Screening Software iPad Apps
--	--

Distance Visual Acuity Testing — Page 15

Wall & Handheld Charts	Pediatric Vision	LEA Symbols
LEA Numbers	Sloan Letters	HOTV

Distance Visual Acuity Screening — Page 17

Spot Handheld Vision Screener	AAPOS Vision Screening Kit	AAPOS Supplemental Vision Screening Kit
-------------------------------	----------------------------	---

Near Vision and Reading Visual Acuity — Page 18

LEA Symbols	LEA Numbers	Lighthouse Near Vision
Sloan Letters	Reading Cards	Spanish Reading Cards

Low Vision and Low Contrast — Page 25

Rehabilitation Charts	Low Vision Devices	Magnifiers & Domes
Low Contrast Charts	Low Contrast Flip Charts	ELCT Contrast Chart

Illuminated Cabinets and Charts — Page 33

ESV3000 & ESC2000 ETDRS	ESV1200 & 1500	Insta-Line Quantum
ESV1018	CSV-1000	Super Pinhole

Stereoacuity Screening — Page 48

Stereoacuity Screening	Random Dot E	Disparity Targets
PASS Smile Tests	Polarized Vectographs	Worth 4-Dot

Examination Tools — Page 52

Prisms	Neutral Density	Trial Frames and Lenses
+/- Flippers	Occluders	Fixation

General Optical Aids — Page 63

Marsden Ball	Eye Models	Tangent Screen
Slit Lamp & Lenses	Alger Brushes	Anti-Suppression

3 Easy Ways to Order from Good-Lite:



Order online at
www.good-lite.com



Order by Phone
800-362-3860



Order by FAX
888-362-2576

Visual Acuity Tests: lea SYMBOLS®

Visual acuity tests are used more than any other tests in ophthalmology, optometry, and in vision screening in preschool, school, and occupational health care. Visual acuity can be measured with several different tests that depict different qualities of vision. The basic test is the line test, which has at least five optotypes (test symbols) on the test lines with spaces between the optotypes

that are equal to the width of the optotypes of that line. The distance between the test lines is equal to the height of the lower line (Visual Acuity Measurement Standard (1984) Consilium Ophthalmologicum Universale, Visual Function Committee).

In the LEA Test System series there are line tests, tests with more crowded symbols, and single symbol tests to allow assessment of

function of the visual system in these three functionally different situations. The number of tests has grown over the years because the needs in screening and assessment of children and adults vary at different ages and different functional levels. For more information, go to the Instruction Section of <http://www.lea-test.fi>.



LEA 3-D Puzzle

The LEA 3-D Puzzle is designed for supporting the development of the concept of "same", first with colors, then with black and white concrete forms that are prerequisites for the measurement of recognition acuity with the LEA Symbols tests at an early age. During the play you can also observe eye-hand coordination and awareness of directions, which are often problems for children with brain damage, even mild brain damage. Package includes the booklet: *Assessing Vision Development Through Picture and Shapes* and instructions. Tray is 6.75" x 6.75" (17.2 cm x 17.2 cm). The 4 symbol puzzle pieces are approximately 2" x 2" (5.1 cm x 5.1 cm) each.

..... 251600



LEA Symbols Playing Cards

These specially designed playing cards make it easy to measure near visual acuity in young children, while teaching concepts of similar/different, big/small, bigger/smaller. The 4 packs contain 16 cards each with symbols of varying sizes: 16M and 10M, 6.3M and 4M, 2.5M and 1.6M, and 1M and 0.63M. Includes training cards and instructions. The paper cards measure 1.6" x 2.75" (4 cm x 7 cm).

..... 252500



LEA Symbols Domino Cards

Assess the near visual acuity of 2- and 3-year-old children, or use the cards at any age (older than 2+ years) to train an amblyopic eye. Set contains 3 packs of 16 domino-like cards with a large symbol at one end and a small symbol at the opposite end. Symbol sizes are 3.2M and 0.8M, 2.0M and 0.5M, and 1.25M and 0.3M. Each symbol on the card is printed with its M unit, decimal value, 6 meter, and 20 foot equivalency for visual acuity value at 40 cm testing distance. Includes training cards and instructions. Paper cards measure 1.6" x 2.75" (4 cm x 7 cm).

..... 251500



LEA Symbols Near Vision Screener (16 inches/40 centimeters)

This test is designed for amblyopia screening and testing children who have difficulties when tested with the standard near vision card (#250800). The background is gray with the exception of a white rectangle spaced at 100% surrounding each line of symbols. Line sizes from 20/400 to 20/16 (6/120 to 6/4.8) equivalent, 0.05 to 1.25. Printed on two sides. Response key is printed on the test card. Includes training cards and instructions. Screener measures 8" x 10" (20.3 cm x 25.4 cm).


..... 252000



LEA Symbols Near Vision Card (16 inches/40 centimeters)

This test measures near visual acuity with proportionally spaced (logMAR) lines on the front and lines with 25% and 50% spacing on the back. The tighter spacing often reveals difficulties in seeing tightly spaced details, which should be known when a child is learning to read. Line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Response key printed on test card. Includes training cards and instructions. Card measures 8" x 10" (20.3 cm x 25.4 cm).

..... 250800



LEA Symbols Pocket-Sized Near Vision Card (16 inches/40 centimeters)

This test is designed for testing patients' visual acuity during rounds and in similar situations with a need for quick testing. The test has proportionally spaced (logMAR) lines on one side and tightly crowded lines (50% spacing) on the reverse side. Line sizes from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. 3.75" x 6.25" (9.5 cm x 15.9 cm).

..... 250950



LEA Symbols Distance Vision Screener (10 feet/3 meters)

This test is designed for amblyopia screening and testing children who have difficulties when tested with the standard distance charts (#250250 or 250150). The background is gray with the exception of a white rectangle spaced at 50% surrounding each line of symbols. Line sizes from 20/80 to 20/16 (6/24 to 6/4.8) equivalent, 0.25 to 1.25. Includes response key and training cards. Pages are offset for easy flipping. Instructions are printed on the inside cover. 10 pages, 3.25" x 9" (8.25 cm x 22.9 cm).

..... 259900



LEA Symbols 10-Line Distance Chart (10 feet/3 meters)

This test is the most commonly used test in the health care of children at the age of 3 to 4 years and in amblyopia screening in many countries.

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/16 (6/60 to 6/4.8) equivalent, 0.25 to 1.25. The 20/160 (6/48) and 20/100 (6/30) lines have been omitted. The reverse side is black and can be used as the test surface for the color vision test (#260260). Has a hole for hanging the chart on the wall. Includes response key, training cards, and instructions. 16" x 16" (40.6 cm x 40.6 cm). Folds to 16" x 8" (40.6 cm x 20.3 cm).

..... 250250



LEA Symbols Crowded Symbol Book (10 feet/3 meters)

The Crowded Symbol Book is designed for testing children and adults who have difficulty in reading lines of optotypes. In this test, the child or adult identifies only the middle optotype. The four symbols surrounding the central symbol similarly to the way optotypes are crowded on line acuity charts. As a result of this crowding effect, this test often permits measurement of line acuity in young children much earlier than the line visual acuity charts and, thus, helps in early detection of amblyopia and impaired vision. Line sizes from 20/63 to 20/16 (6/19 to 6/4.8) equivalent, 0.63 to 1.25. Includes response key, training cards, and instructions. 16 offset pages, 5" x 5" (12.7 cm x 12.7 cm).

..... 250700



LEA Symbols Single Symbol Book (10 feet/3 meters)

The Single Symbol Book is designed for testing children and adults who have difficulties in perception of symbols in lines. The visual acuity values measured with single symbols are often several lines better than the line acuity values in persons with impaired vision and, thus, overestimate the quality of visual functioning. Single optotype line sizes from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Includes isolating mask, response key, training cards, and instructions. 13 offset pages for easy flipping. 5" x 5" (12.7 cm x 12.7 cm).


..... 250600



LEA Symbols Flash Cards (10 feet/3 meters)

This set of flash cards was designed for testing distance visual acuity in children and adults who have difficulty with visual acuity charts, the Crowded Symbol Book Test, or the Single Symbol Book Test. Line sizes from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Includes response key, training cards, and instructions. 30 cards per set, 4" x 5" (10 cm x 12.7 cm) to 5" x 5" (12.7 cm x 12.7 cm).

..... 252700



LEA Symbols Cards, 40M and 60M

The 5 LEA Symbols Cards, 40M and 60M, can be presented as single optotypes or arranged in a row as a line test. When the cards are next to each other, the space between the symbols is equal to the width of the symbols. To determine the visual acuity use the following formula:

$$VA = \frac{\text{Viewing Distance Used (meters)}}{M\text{-value}}$$

60M cards are approx. 8.6" x 8.6" (21.9 cm x 21.9 cm). 40M cards are approx. 6" x 6" (15.25 cm x 15.25 cm).

..... 256100



Leo Learns By Doing (DVD, 17 minutes)

When a blind child is born into a family, the situation is new and confusing. Parents have thousands of questions,

but very few answers, such as:

- How am I able to do everything possible for my child?
- How can our everyday life become ordinary, good, and safe, yet active?
- How can I best support my child?
- Am I overprotective or encouraging as a parent?
- How will my blind child receive as many experiences as the seeing child?

These questions and more are answered in this DVD.


..... 790000



LEA Symbols 15-Line Distance Chart (10 feet/3 meters)

This test is designed for vision screening at schools and assessment of vision of children with special needs beginning with age of 3 years. Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. The reverse side is black and can be used as the test surface for the color vision test (260260). Has a hole for hanging the chart on the wall. Includes response key, training cards, and instructions. 17" x 21" (43 cm x 53.4 cm). Folds to 17" x 10.5" (43 cm x 26.7 cm).

Set **250150**




LEA Symbols 13-Line Translucent Distance Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/8 (6/38 to 6/2.4) equivalent, 0.16 to 2.50 with repetitive left eye/right eye groups. This chart can either be wall mounted or used in an illuminated cabinet found on page 26. Includes response key, training cards, and instructions.

Set **252400**

Chart Only **252401**




LEA Symbols 13-Line Translucent Distance Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/8 (6/38 to 6/2.4) equivalent, 0.16 to 2.50. This chart can either be wall mounted or used in an illuminated cabinet found on page 26. Includes response key, training cards, and instructions.

Set **250400**

Chart Only **250401**




LEA Symbols 12-Line Translucent Distance Chart (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/100 to 20/8 (6/30 to 6/2.4) equivalent, 0.20 to 2.50. This chart can either be wall mounted or used in an illuminated cabinet found on page 26. Includes response key, training cards, and instructions.

Set **255100**

Chart Only **255101**



LEA Symbols 15-Line Translucent ETDRS-Style Distance Chart for ESV3000 Cabinet (10 feet/3 meters)

LEA Symbols distance visual acuity test for the ETDRS Illuminated Cabinet. Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Includes response key, training cards, and instructions. 24.5" x 25.5" (62 cm x 65 cm).

..... **250300**



LEA Symbols Chart for Vision Rehabilitation (40 inches/1 meter)

This test is designed for testing children and adults with visual impairment at a distance of 1 meter. The largest symbols are 50M in size, corresponding to 1m/50M = to 1/50 or decimal visual acuity 0.02, equivalent values 20/1000, 6/300. The smallest symbols are 1M in size corresponding to 1m/1M = to 1/1 or decimal visual acuity 1.00, equivalent values 20/20, 6/6. Designed courtesy of August Colenbrander, MD. Chart measures 16" x 21.5" (40.7 cm x 54.6 cm). Folds to 16" x 10.75" (40.7 cm x 27.3 cm).

..... **258000**

Visual Acuity Tests:



LEA Numbers were developed for international studies in occupational health and school health services. Letters like the Sloan letters are used in the

Western world and in the old colonies of western empires but not in a large part of the

Lea-numbers®

rest of the world. The Arabic numbers are known world over, even in countries that use other types of numbers. The four numbers 5, 6, 8, and 9 were chosen because they could be designed to blur equally.

The numbers 5, 6, 8, and 9 were designed by Dr. Hyvärinen in 1993 and calibrated with the Landolt C (Juvi Mustonen, OD; Risto Näsänen, PhD; and Lea Hyvärinen, MD, Vision Laboratory, Department of Physiology, University of Helsinki, Finland 1994).



LEA Numbers 15-Line Translucent ETDRS-Style Distance Chart for ESV3000 Cabinet (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50.

..... 270200



LEA Numbers Near Vision Card (16 inches/40 centimeters)

Designed to measure visual acuity with a standard line test and with three tighter spaced tests; 50%, 25%, and 12.5%, which often reveal difficulties in seeing tightly crowded details. Line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Card measures 8" x 10" (20.3 cm x 25.4 cm).

..... 270900



LEA Numbers Near Vision Card (24 inches/60 centimeters)

Designed to measure visual acuity with a

standard line test and with three tighter spaced tests; 50%, 25%, and 12.5%, which often reveal difficulties in seeing tightly crowded details. Line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Card measures 8" x 10" (20.3 cm x 25.4 cm).

..... 274100



LEA Numbers Pocket-Sized Near Card (16 inches/40 centimeters)

This test is designed for testing patients' visual

acuity during rounds and in similar situations with a need for quick testing. The test has proportionally spaced (logMAR) lines on one side and tightly crowded lines (50% spacing) on the reverse side. Line sizes from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. 3.75" x 6.25" (9.5 cm x 15.9 cm).

..... 271000



LEA Numbers 15-Line Distance Chart (10 feet/3 meters)

This test is designed for vision screening at schools and assessment of vision of children with special needs beginning with age of 3 years. Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Has hole for hanging on the wall. 18" x 20" (45.7 cm x 50.8 cm). Folds to 18" x 10" (45.7 cm x 25.4 cm).

..... 271100



A LEA Numbers 13-Line Translucent Distance Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/8 (6/38 to 6/2.4) equivalent, 0.16 to 2.50.

..... 271200

B LEA Numbers 12-Line Translucent Distance Chart (13 feet, 4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/100 to 20/8 (6/30 to 6/2.4) equivalent, 0.20 to 2.50.

..... 271300



LEA Numbers Chart for Vision Rehabilitation (40 inches/1 meter)

This test is designed for testing children and adults with visual impairment at a distance of 1 meter. The largest symbols are 50M in size, corresponding to 1m/50M = to 1/50 or decimal visual acuity 0.02, equivalent values 20/1000, 6/300. The smallest symbols are 1M in size corresponding to 1m/1M = to 1/1 or decimal visual acuity 1.00, equivalent values 20/20, 6/6. Designed courtesy of August Colenbrander, MD. Chart measures 16" x 21.5" (40.7 cm x 54.6 cm). Folds to 16" x 10.75" (40.7 cm x 27.3 cm).

..... 272100

PLEASE NOTE: 9 x 14 LEA Numbers Charts can either be wall mounted or used in an illuminated cabinet found on page 33.



LEA Numbers Single Number Book (10 feet/3 meters)

The Single Number Book is ideal for situations in which the child/adult cannot perform when tested with a line test. This book can help determine the best resolution of the visual pathway in a child/adult with amblyopia or impaired vision. Single optotype line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Includes isolating mask and instructions. 13 offset pages for easy flipping. 5" x 5" (12.7 cm x 12.7 cm).

..... **274000**



LEA Numbers Cards, 40M and 60M

These five LEA Numbers Cards can be used for presentation as single optotypes or arranged in a row as a line test. When the cards are next to each other, the space between the numbers is equal to the width of the symbols. To determine the visual acuity use the following formula:

$$VA = \frac{\text{Viewing Distance Used (meters)}}{M\text{-value}}$$

60M cards are approx. 6" x 9" (15.25 cm x 22.9 cm). 40M cards are approx. 4" x 6" (10.2 cm x 15.25 cm).

..... **272400**

Contrast Sensitivity

The contrast sensitivity tests that Dr. Hyvärinen has developed for clinical examinations and research projects have been used in the evaluation of vision of thousands of patients in numerous countries revealing changes in visual functioning that

were not detected with the high contrast visual acuity tests. Measurement of contrast sensitivity should be a part of routine examinations because vision at low and intermediate contrast levels is important in communication and orientation in space.




LEA Symbols Low Contrast Test, 10M

Test for quick measurement of low contrast visual acuity at 25%, 10%, 5%, 2.5%, and 1.25% contrast by measuring the distance where the symbols are seen. 10M symbols correspond to visual acuity 0.10 at 1 m, 0.5 m, and 0.3 at 3 m. VA values are, thus, easy to calculate. Test includes a table for conversion to Snellen equivalent VA values, a response key, training cards, recording forms, and instructions. 6 pages, 3.25" x 9" (8.25 cm x 22.9 cm).

..... **251100**

Mixed Order **251200**



LEA Symbols Translucent Low Contrast Charts for Illuminated Cabinet

A low contrast test is a special way to measure visual acuity to determine how well eyes function in low light situations and/or to determine how well eyes distinguish objects from similarly colored or shaded backgrounds. Testing is identical to the measurement of visual acuity at high contrast level. The 2.5% test is the most practical test in clinical use. Available in 5 contrast levels: 1.25%, 2.5%, 5%, 10%, and 25%. Includes response key, training cards, recording forms, and instructions.

- A** 1.25% **253700**
- B** 2.5% **253800**
- C** 5%..... **253900**
- D** 10%..... **254000**
- E** 25% **254100**



LEA Numbers Low Contrast Test, 10M

Test to easily measure, record, and detect changes in the transfer of visual information when the change only affects visual acuity at low contrast levels. Chart presents numbers at the following contrast levels: black, 25%, 10%, 5%, 2.5%, and 1.25%. Pages are offset for easy flipping. Each chart includes recording forms and instructions. 6 pages, 3.25" x 9" (8.25 cm x 22.9 cm).

..... **270400**

Mixed Order **270500**



Hiding Heidi Low Contrast Face Test

See page 8 for more details.

Four cards printed on both sides in the following contrast levels: black, 25%, 10%, 5%, 2.5%, and 1.25%. Cards are 9" x 9" (23 cm x 23 cm). Includes instructions.

..... **253500**



LEA Low Contrast Grating Acuity Test

Measurement of contrast sensitivity curve is based on 4 gratings at 2 contrast levels, 10%, and 2.5% contrast using 2 grating frequencies: 0.5 cpcm and 8 cpcm (10% level) or 0.5 cpcm and 4 cpcm (2.5% level). Testing is based on finding the distance where the person is at threshold in seeing the direction of the grating lines. Includes nomograms, instructions, and recording forms. Measures 8" (20 cm) in diameter.

..... **253400**

Cognitive Vision



Visual information undergoes several neural processes before it is perceived as images. When vision changes, the changes may occur either in the image quality, due to disorders of the eyes and/or visual pathways, or picture perception, due to loss of function in some of the specific analytic functions in the visual cortices. The number of visually impaired children with brain damage-related visual impairment is increasing in all countries.

Tests for cognitive vision were designed for early diagnosis of some usual cognitive problems.

It is important that all cognitive visual functions are assessed, from seeing the orientation of lines and size of objects to different recognition functions and more complicated perceptual tasks. For example, a child who has lost recognition of facial features only may be misunderstood when he/she does not greet people whom the child knows well. Likewise, children who do not perceive facial expressions may be taken as autistic when they cannot function in a group of toddlers whose communication is through body language. For more information, go to the Instructions Link:

www.lea-test.fi

LEA Numbers Low Contrast Charts for Illuminated Cabinet

A low contrast test is a special way to measure acuity to determine how well eyes function in low light situations and/or to determine how well eyes distinguish objects from similarly colored or shaded backgrounds. Testing is identical to the measurement of visual acuity at high contrast level. The 2.5% test is the most practical test in clinical use. Available in 5 contrast levels: 1.25%, 2.5%, 5%, 10%, and 25%. Includes response key, training cards, recording forms, and instructions.

- A** 1.25% **271500**
- B** 2.5% **271600**
- C** 5% **271700**
- D** 10% **271800**
- E** 25% **271900**



LEA 3-D Puzzle

The LEA 3-D Puzzle is designed for training and assessing typically developing toddlers and other young children. The puzzle also serves as an assessment tool for individuals with brain damage, or for older children and adults with low cognitive abilities. Includes instructions and the booklet: *Assessing Vision Development Through Pictures and Shapes*. Tray is 6.75" x 6.75" (17.2 cm x 17.2 cm). The 4 symbol puzzle pieces are approximately 2" x 2" (5.1 cm x 5.1 cm) each.

..... **251600**


Visual Field



LEA Flicker Wand

The LEA Flicker Wand is 26.5" (67.3 cm) long with a 20" (50.8 cm) flexible wand with a diode at the end that can be used as either a flickering or non-flickering stimulus at 4, 40, and 400 cd/m². The curved wand is moved from behind the person forward in all four quadrants of the visual field to detect losses in the periphery of the visual field.

..... **280000**



Small and Medium Heidi Fixation Targets

The small 2" (5 cm) diameter "smiling face" fixation paddle was developed based on studies by Robert Fanz. Features include grating on reverse side, small pictures of a bunny and a giraffe for demanding fixation, and pupil size and corneal diameter measurement on the handle. The medium 5" (13 cm) diameter fixation paddle is for assessing children with severe disabilities to determine the distance at which the children respond to a high-contrast face figure.

- A** Medium Target **253100**
- B** Small Target **253000**



LEA Grating Acuity Test, a Resolving/Discrimination Test

This test at full contrast requires discrimination of

the direction of long lines, which is a more demanding task than the detection acuity measured with LEA GRATINGS. The test gratings are presented to children and adults at different distances to determine which distance the direction of lines can be perceived. Knowing the distance and which grating was used, the grating acuity values can be read from the nomograms. Includes, nomograms, instructions, and storage case. 4 disks 8" (20 cm) in diameter.

..... **251300**



LEA GRATINGS, a Preferential Looking Test
LEA GRATINGS are for preferential looking test situations with infants or children and adults with disabilities to determine detection acuity. Grating levels printed on each handle are: 0.25, 0.5, 1.0, 2.0, 4.0, and 8.0 cpcm (cycles per centimeter of surface). Includes instructions and storage case. 4 paddles 8" (20 cm) in diameter.

..... **253300**



Hiding Heidi Low Contrast Face Test

Visual communication is the most important communication method during the first year of life. The ability to detect objects of low contrast is an important component of the visual system. For example, facial expressions are mediated by faint shadows and changes in the contours of the mouth and eyes. Determining the levels of contrast that an infant can detect helps provide baseline data for evaluating future changes. For example, deviations from usual behavior may indicate disorders that leave vision at high-contrast levels unaffected.

Four cards printed on both sides in the following contrast levels: black, 25%, 10%, 5%, 2.5%, and 1.25%. Cards are 9" x 9" (23 cm x 23 cm). Includes instructions.

..... 253500



Heidi Expressions Test Game

Some children with visual impairments may achieve near normal results in routine vision tests (large visual field and normal or near normal visual acuity), but cannot interpret facial expressions or recognize people's faces. The Heidi Expressions Test Game improves early evaluation of vision for communication. 18 cards per set. Cards are 4" x 4" (10.16 cm x 10.16 cm).

..... 254505



LEA Rectangles Game

This educational toy is a modified version of Effron's rectangles and is used to assess size perception. The Rectangles Game contains two groups of rectangles. The total surface area of each rectangle is the same. A child with normal eyesight can distinguish between the longer and the shorter rectangles—often before the age of 3 years. If the child has difficulty making this distinction or does not understand the instructions, the LEA Rectangles Game can help determine whether or not the child has a specific visual acuity problem. Includes instructions and 10 colored rectangles (5 black, 5 gray) in 5 sizes with a total area of 4" (10.16 cm) each.

..... 254600



LEA Mailbox Game

Used to assess visual perception of line orientation. Tests two components:
1) information for

the hand movements in the parietal lobe and 2) picture perception in the inferior temporal lobe. To observe the first function, the child is asked to drop a card through the slot of the LEA Mailbox Game. The child should demonstrate perception of the orientation of the slot by turning the card to the correct position before moving it half the distance to the slot (dorsal stream function). A ventral stream function can be observed through interactive play using the LEA Mailbox Game. If a child cannot visually perceive the orientations of lines, drawing them may be difficult, and drawing angles will make no sense. Includes instructions.

..... 254700



"WHAT and HOW Does this Child See?"

Is a question that is asked each day by many professionals in early intervention and education about children with disabilities. This book attempts to help you structure your questions so that you can look for answers together with other members of your team. This book is a collection of questions and answers from lectures presented by **Lea Hyvärinen, MD, PhD, FAAP**, and **Namita Jacob, PhD**, in a number of countries.

..... 791000

Adaptation Test



Cone Adaptation Test

The test consists of 15, 2" x 2" (5 cm x 5 cm) red, blue, and white plastic chips designed to help parents, teachers, and doctors assess children's visual difficulties in twilight. The test situation is not a formal test of night vision, but gives useful information on visual adaptation to lighting changes.

..... 252900

Color Vision Test



Panel 16 Quantitative Color Vision Test

The Panel 16 Color Vision Test includes 16 test caps and 17 restriction rings with one of the caps being a "pilot". The diameter of the stimulus area is 1.3" (3.3 cm), and can be reduced by using dark gray restriction rings with an opening of 0.47" (1.2 cm) in diameter. Acrylic case includes a cloth storage pouch. Travel tube comes with black cloth. Both come with recording forms and instructions. Single Set

A With Acrylic Case 260200

B With Travel Tube 260260

Double Set (Shown on page 10)

With Acrylic Case 260100

With Travel Tube 260160

Color Vision Tests




Good-Lite ColorCheck Complete Vision Screener

This test offers both LEA Numbers and Pediatric LEA Symbols screening for protan/deutan (red/green) and plates screening for tritan (blue/yellow). Color plates are larger making the screening of children, adults and low vision patients much easier.

Includes user-friendly instructions, tabbed sections, and protective coating on each test page.

..... 730050



ColorCheck Handheld Color Vision Screening Tablet

Color vision screening is now as easy as using your favorite app. The M&S Smart System Tablet (Google Nexus 7) comes pre-installed with ColorCheck color vision application to easily and quickly screen both adults and children. ColorCheck uses Pediatric LEA Symbols and Numbers to check for protan/deutan/tritan (red/green/blue) color deficiencies. The protocol is completely randomized ensuring accurate results.

..... 730040





Richmond HRR Standard Pseudoisochromatic Test

Quickly screens for normal color vision and congenital and acquired color vision defects, classifying by type—protan/deutan (red/green) and tritan (blue/yellow)—and extent of defect. This test is ideal for pediatric testing, because it uses simple shapes (O, X, Δ) that are easily recognized by children. Special non-glare laminate keeps colors from becoming distorted and prevents fingerprint smudges. Reversible occluding glasses are included.

Unlaminated 730006

Laminated 730005



Richmond HRR/Amsler Combination Color Book

HRR and Amsler Grid combined into one easy to use book. The book contains 7 different amsler charts including Standard Chart, Diagonal, Red on Black Standard and Minute Juxta-Central.

Laminated 730020

Unlaminated 730021



HMC Oculus Anomaloscope

For precision diagnosis of color vision in the red/green and blue/green areas with integrated automatic neutral adaptation. The unit is normally controlled using the MS-Windows® based software and included cable from a PC or laptop. Can also be controlled using the optional keypad (P/N 4667) with LCD display.

- Meets requirements of Nagel and DIN-Norm 6160
- Adjustable automatic neutral adaptation (standard 'C' lighting)
- Standard test field of 2°
- Based on additive mixing of colors

Pentium 2 or newer and Windows XP Only. Interface cable with RS-232 9 pin (RS 232 / V24, Dub D-jack, 9-pole) serial port.

..... 466600



LCD Control Panel for Anomaloscope

LCD control panel provides an alternative way to control the Oculus Anomaloscope.

..... 466700

Arrangement Color Vision

Panel 16 Quantitative Color Vision Test

The Panel 16 Color Vision Test sets display the same hues as the Farnsworth D-15. The diameter of the stimulus area is 1.3" (3.3 cm), and can be reduced by using dark gray restriction rings with an opening of 0.47" (1.2 cm) in diameter. Travel tube comes with black cloth.



Single Set

Includes 16 test caps, 16 restriction rings, acrylic case, cloth storage pouch, recording forms, and instructions.

- A With Case 260200
- B With Travel Tube..... 260260



Double Set

Includes 2 sets of 16 test caps (32 total), 32 restriction rings, acrylic case, cloth storage pouch, recording forms, and instructions.

- C With Case 260100
- D With Travel Tube 260160



Farnsworth D-15 Color Test Kit

Quickly and easily screen for red-green color deficiencies by asking subjects to separate 15 color disks. The Farnsworth D-15 allows basic color deficiency screening. The diameter of the stimulus area is .47" (1.2 cm).

..... 730022



Lanthony Desaturated Color Test

The Lanthony Desaturated Color test is similar to the Farnsworth D-15, but the 15 disks are desaturated (lighter in hue). This allows the screener to determine whether or not the subject has the ability to discriminate between subtle color variations. The test can distinguish between normal color perception and mild deficiency in protan/deutan (red-green) or tritan (blue-yellow). The diameter of the stimulus area is .47" (1.2 cm).

..... 730023



Richmond Magnetic Color Test

Test color vision by moving a magnet to arrange the colored disks. This complete Farnsworth D-15 test is enclosed in a sealed plastic box to prevent the disks from becoming lost or dirty. The diameter of the stimulus area is .47" (1.2 cm).

..... 730004




Richmond Desaturated Magnetic Color Test

Assess fine color discrimination with the use of this desaturated Farnsworth D-15 test. The color disks come in an enclosed plastic box and are arranged using a magnet. The colors used in this set are desaturated (lighter in hue), so that finer discrimination is required to arrange the disks. The diameter of the stimulus area is .47" (1.2 cm).

..... 730025

Arrangement Color Vision Size Comparison



Good-Lite offers arrangement color tests with color chips in 2 sizes. The Panel 16 Quantitative Color Test measures 32 mm in diameter with a "Reduction Ring" that restricts the target to the standard size of 12 mm.



Farnsworth 100 Hue Color Test

Farnsworth 100 Hue Test comes boxed in clear acrylic containers which make scoring easy and reduces the potential for test disruption by holding the patient's disc sequence intact. The Farnsworth 100 Hue test includes 93 color discs; 4 clear acrylic boxes; a tray; a durable, laminated score template for copying; a supply of non-latex, lightly powdered gloves for the patient's hands to maintain cleanliness and accuracy of the color discs; comprehensive instructions; and a CD disk with a scoring template for use with MS Excel. The scoring template calculates a numerical score and provides a graphical display or print-out of the patient's score. Dimensions: 56 cm x 17 cm x 2.4 cm.

..... **730024**



Color Test Daylight Illuminator with Platform

The Color Test Daylight Illuminator with Platform is specifically designed to provide the optimal light when using Panel 16, Farnsworth 100 Hue, Farnsworth D-15 and Lanthony Desaturated D-15 color tests. The Daylight bulb gives off 6280 degree K light which closely matches the standards set by the International Commission of Illumination and had been approved by the U.S. Food and Drug Administration and the U.S. Defense Personnel Support Center for the armed forces. The removable platform measures 10" x 18" (26 cm x 46 cm) and includes a black rubber dull mat.

Daylight Illuminator comes with Universal Power Supply adapting to any combination of 50 or 60 Hz and 110 or 220 volt.

..... **730027**



Color Test Daylight Illuminator with Easel

The redesigned easel illuminator provides the proper angle for Pseudoisochromatic Plates and Cap Style tests. The easel helps avoid damaging fingerprints on the plates because the plates are presented standing on their own. The easel is removable for installation of a flat tray in order to use the Illuminator with arrangement tests. Includes convenient hand grips for ease of transport, light switch, and a 4-foot cord.

..... **730003**



Arrangement Test Tray

Raised platform tray for the Magnetic Farnsworth or Roth 28 tests is 10" x 18" x 1.75" (26 cm x 46 cm x 3.5 cm) and includes a black dull rubber mat.

..... **730028**



Color Book Easel Tray

Raised easel tray for pseudoisochromatic book when viewed using the Color Test Daylight Illuminator. Places the color book at the proper viewing angle for testing.

..... **730018**

Distance Computerized Vision Testing System (VTS)



VTS Complete Computerized Vision Testing System

The Good-Lite Company brings our extensive 80+ years of vision testing experience to the Vision Testing System (VTS), a seamless transition from projector charts to state-of-the-art computerized vision testing system!

VTS Testing System Features

Multiple Optotypes in One System

Now easily test all your patients; children as well as adults. The VTS system contains LEA Symbols, LEA Numbers, HOTV, Sloan Letters, and more. Display characters in a variety of formats including logMAR, vertical/horizontal single lines, crowded or uncrowded single optotype, and low contrast fonts. Simply change the optotype or display format with the press of a button on the included remote control.

Secondary Calibration Distance

Seamlessly switch between 2 calibration distances on your VTS System to accommodate a clinical trial requirement, clinical trial requirement, patient with low vision, or a patient who uses a wheelchair.

LEA Symbols Pediatric Testing

Unlike optotypes that look like LEA Symbols, the LEA Symbols have been validated, calibrated against the Landolt C, and are the only pediatric optotypes to blur equally at threshold to prevent guessing.

LEA Symbols are on the recommended list of eye charts in a consensus policy statement of the AAP, AOA, AAPOS, and AACO.

Optotype Indicator

The advanced optotype indicator allows the clinician to easily direct both children and adults as to which optotypes to identify.

VTS Includes:

- | | |
|---------------|---------------------------|
| Sloan Letters | Low Contrast Testing |
| LEA Symbols | Crowding/Interaction Bars |
| LEA Numbers | Worth 4-Dot |
| Landolt "C" | Phoria Fixation Disparity |
| HOTV | Three Lines Same Size |
| Tumbling "E" | Single Vertical Line |
| Red/Green | Color Vision Testing |

Performance Attributes:


- Optotype Indicator
- Unlimited User-Defined Protocols
- Cycling Optotype/Multi-Lines/Letters
- Complete Randomization
- Streamlined Video Select and Display
- Secondary Calibration Distance
- Quick Copy Video Control
- ETDRS Selectable
- Fixation Targets/Stimulus

VTS Software with Computer* 300100

VTS Software Only 300105

* Please contact Good-Lite for details on computer


VTS Features



LEA Symbols with Indicator




High Contrast Characters



Single Number with 50% Crowding



Red/Green Test



Worth 4-Dot Testing



Color Vision Testing

Distance Vision Screening is now Fun and Easy!



VISIONquest 20/20

Because Every Child Deserves To See



Distance

Color

Stereoacuity



EyeSpy 20/20™ Vision Screening Software

EyeSpy 20/20™ employs a revolutionary proprietary approach to vision screening technology - a software program created by the Founders of VisionQuest 20/20. A child's vision is assessed in a standardized, cost-effective manner while the child "plays" a 2- to 3-minute video game.

EyeSpy 20/20™ has been independently validated by the Medical University of South Carolina's Storm Eye Institute for scientific validity and reliability. The results of the double-blind study suggest that EyeSpy 20/20™ provides results comparable to those achieved by a certified ophthalmic technician using a 'gold standard' device for assessing vision.



Distribution of our vision screening program is focused on schools, through government health programs, pediatricians, and other health care professions. Schools represent the most efficient distribution channel because they offer the single largest venue for reaching children.

EyeSpy 20/20™ System Features


- Easy-to-implement program and will operate on most computers currently in schools.
- Automated, age-appropriate vision screening protocols, for ages 3 years old and up.
- Professional, easy-to-read report documenting screening results in either English or Spanish.




- Integrates with EMR and EHR Systems
- HIPAA Compliant Data Storage
- Results are automatically and immediately available for school health records as well as for the child's guardian.
- Child's screening history is retained in a database, enabling rescreening without record duplication as well as reports that can be reprinted at any time, individually or collectively.
- Can be installed on an unlimited number of computers within a school.
- Offers schools and other organizations a fund-raising opportunity while protecting their students vision.

EyeSpy 20/20™	300400
Annual Renewal License	300401
Testing Goggles.	300402
UBS Mouse 10' Extension Cable. . .	300403
Eye Patches	300404


Distance Vision Screening



Critical Line Testing with Indicator



MassTest Acuity Format



Optotypes with 50% Crowding Bars

VSS: Computerized Vision Screening Software

The Good-Lite® Company has developed a computer-based Vision Screening System (VSS) that provides acuity tests for children as young as 3 years to adults.

The screening software will display LEA Symbols, HOTV, LEA Numbers, and Sloan Letters optotypes ranging from 20/100 to 20/10 (6/30 to 6/3) equivalent, 0.20 to 2.00. Optotypes can be displayed as either single line, single line with 50% spaced rectangles, or single optotypes with 50% crowding bars. The VSS can also quickly display critical lines only of 20/40 (6/12) for 3-5 yr olds and 20/32 (6/9.5) for 6 and older.

The screening software can be quickly calibrated to a testing distance of either 5, 10, or 13 feet.

Features:

- LEA Symbols, LEA Numbers, HOTV, and Sloan Letter optotypes.
- Critical line testing for 3-5 year olds and 6 and older.
- Single optotype with 50% crowding bars.
- Lines with 50% crowding rectangles.
- Threshold acuity testing.

Easily cycle through each acuity line using the included remote control. Run the VSS Screening software on your existing computer equipment. Windows XP and Vista compatible.

This kit includes:

- Vision Screening Software
- USB Remote Control

Requirements

Intel Celeron or AMD Athlon CPU, 512 MB RAM, 40 GB hard drive, Windows XP.

..... **300000**

iPad Vision Apps

for iPads with high resolution displays only



Near/Distance Vision Screening Testing



Good-Lite's Near/Distance Vision Screening and Testing app for the iPad 3 with retina display features evidence-based, proportionally spaced eye chart designs and optotypes for screening and testing both distance and near vision. App also includes Continuous Reading Text.

Near and distance vision threshold charts include LEA Symbols, HOTV, LEA Numbers, and Sloan Letters. Toggle button provides ability to change charts to prevent memorization. Each eye chart and Continuous Reading Text chart has 3 charts choices.

App settings permit screening or testing at different distances and measurements.

Available at iTunes Store



AAPOS Vision Screening App



Quickly and easily screen the vision of children and adults beginning at age 3 years with evidence-based LEA Symbols and Sloan Letters. App is based on the basic distance Vision Screening Kit created by the American Association of Pediatric Ophthalmology and Strabismus (AAPOS).

Charts include both full-chart threshold and critical line formats using either LEA Symbols or Sloan Letters optotypes. Critical line screening is a quick method for screening where individuals are screened using only the line that matches their age.

Each test format can be used for distance screening at 10'/3 m or for near vision screening at 16"/40 cm.

Available at iTunes Store

Distance Acuity Wall Charts



LEA Symbols 10-Line Distance Chart (10 feet/3 meters)

This test is the most commonly used test in the health care of children at the age of 3 to 4 years and in amblyopia screening in many countries.

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/16 (6/60 to 6/4.8) equivalent, 0.25 to 1.25. The 20/160 (6/48) and 20/100 (6/30) lines have been omitted.

The reverse side is black and can be used as the test surface for the color vision test (#260260). Has a hole for hanging the chart on the wall. Includes response key, training cards, and instructions. 16" x 16" (40.6 cm x 40.6 cm). Folds to 16" x 8" (40.6 cm x 20.3 cm).

..... 250250



LEA Symbols 15-Line Distance Chart (10 feet/3 meters)

This test is designed for vision screening at schools and assessment of vision of children with special needs beginning with age of 3 years. Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. The reverse side is black and can be used as the test surface for the color vision test (260260).

Has a hole for hanging the chart on the wall. Includes response key, training cards, and instructions. 17" x 21" (43 cm x 53.4 cm). Folds to 17" x 10.5" (43 cm x 26.7 cm).

..... 250150

Insist On Genuine LEA Products by GOOD-LITE®

lea SYMBOLS®
lea-numbers®



LEA Numbers 15-Line Distance Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Has hole for hanging on the wall. 18" x 20" (45.7 cm x 50.8 cm). Folds to 18" x 10" (45.7 cm x 25.4 cm).

..... 271100



Sloan Letters Folding Chart (10 feet, 3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. 18" x 18" (45.7 cm x 45.7 cm); folds to 18" x 9" (45.7 cm x 22.9 cm).

..... 735000




Cardiff Pediatric Acuity Test

Designed specifically for acuity measurement in children ages 6 months to 3 years. The test may also be helpful in other age groups for people with intellectual impairment, dementia, head injury, or stroke. Cards can be used at either 1 meter or 0.5 meter. This card set includes 36 cards providing 24 acuity levels when used at both distances and Children's Fun Frame Occluders.

Cardiff Pediatric Acuity Test 474000

Cardiff Low Vision Cards 475400

Cardiff Low Contrast Cards 566700



Pacific Acuity Test (PAC Test)

The Pacific Acuity Test is a vanishing optotype test designed with opposing figures to provide a simple forced-choice testing format appropriate for nonverbal children and individuals with developmental disabilities who are unable to respond to standard optotype recognition tests. Each card has two figures; a simple face, and an opposing oval target of the same overall size but with interior details that allow discrimination from the face at the threshold of acuity.

The full test kit consists of 18 cards with two cards at each level. The cards at each acuity level have a different position of the face target compared to the opposing oval figure. Each card can be presented at either 50 cm or 1 meter but the 1 meter test distance is recommended for greater accuracy and comparison to distance visual acuity tests.

The test range from 1.3 logMAR (20/400) down to 0 logMAR (20/20) allows for testing a broad range of visual acuity.



PAC Test comes with protective canvas case and isolation mask.

..... 360000

Distance Acuity Handheld Tests



LEA Symbols Crowded Symbol Book (10 feet/3 meters)

Variation of the line test for children or adults who have difficulty concentrating on a chart. The center symbol is surrounded on four sides by symbols of the same size. The interaction between the center symbol and the surrounding symbols is approximately the same as in a line test. The test size is convenient and easy to administer during screenings, home visits, and unique testing situations. Often, this is the first line test that can measure visual acuity of a young child with amblyopia, because the test produces the effect of crowding. Line sizes range from 20/63 to 20/16 (6/19 to 6/4.8) equivalent, 0.32 to 1.25. Includes response key, training cards, and instructions. 16 offset pages for easy flipping. 5" x 5" (12.7 cm x 12.7 cm).

..... 250700



LEA Symbols 50% Spaced Crowding Bars Book (10 feet/3 meters)

This book features 50% spaced contour interaction bars to provide single optotype acuity with the proper crowding effect. Optotype line sizes range from 20/100 to 20/20 (6/30 to 6/6) equivalent, 0.20 to 1.00. Includes response key, training cards, and instructions. 17 offset pages for easy flipping. 4" x 4" (10.2 cm x 10.2 cm).

..... 259200



LEA Symbols Single Symbol Book (10 feet/3 meters)

The Single Symbol Book is ideal for people who cannot perform when tested with a line test. This book can help determine the best resolution of the visual pathway in a child/adult with amblyopia or impaired vision. Single optotype line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Includes isolating mask, response key, training cards, and instructions. 13 offset pages for easy flipping. 5" x 5" (12.7 cm x 12.7 cm).

..... 250600



LEA Symbols Distance Vision Screener (10 feet/3 meters)

This test is designed for amblyopia testing or screening. The background is gray with the exception of a white rectangle spaced at 50% surrounding each line of symbols. Line sizes range from 20/80 to 20/16 (6/24 to 6/4.8) equivalent, 0.25 to 1.25. Includes response key and training cards. Pages are offset for easy flipping. Instructions are printed on the inside cover. 10 pages, 3.25" x 9" (8.25 cm x 22.9 cm).

..... 259900



LEA Symbols Flash Cards (10 feet/3 meters)

This set of flash cards was designed for testing distance visual acuity in children and adults who have difficulty with visual acuity charts, the Crowded Symbol Book test or the Single Symbol Book test. Line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Includes response key, training cards, and instructions. 30 cards per set, 4" x 5" (10 cm x 12.7 cm) to 5" x 5" (12.7 cm x 12.7 cm).

..... 252700



Single HOTV Letter Book With 50% Spaced Crowding Bars (10 feet/3 meters)

This book features 50% contour interaction bars to provide single optotype acuity with the proper crowding effect. Optotype line sizes range from 20/100 to 20/16 (6/30 to 6/4.8) equivalent, 0.20 to 1.25. Includes response panel, training cards, and instructions. 17 offset pages for easy flipping. 4" x 4" (10.2 cm x 10.2 cm).

..... 600350



LEA Numbers Single Number Book (10 feet/3 meters)

The LEA Numbers Single Number Book is ideal for situations in which the child/adult cannot perform when tested with a line test. This book can help determine the best resolution of the visual pathway in a child/adult with amblyopia or impaired vision. Single optotype line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Includes isolating mask and instructions. 13 offset pages for easy flipping. 5" x 5" (12.7 cm x 12.7 cm).

..... 274000



Magnetic Sloan Letters Chart (20 feet/6 meters)

This durable chart is made of plastic-coated, galvanized steel with 2 sets of 20 magnetic strips (40 total) with 5 characters on each strip. Meets all qualifying standards for military aviators (e.g., officer candidates, pilots, air crewmen, etc). Approved by the Navy Bureau of Medicine and Surgery (BUMED). Chart measures 9.5" x 11" (24.1 cm x 28 cm).

20/20 size 400500
20/40 size (not shown) 400900

Spot Handheld Vision Screener from PediaVision



VISION SCREENING



Screener Love Spot

Whether you're screening 10 or 10,000 kids, any member of your team can conduct vision screenings. Spot is as easy to use as a camera and does the work for you; most operators become proficient within a few practice screenings.


Spot provides a reliable, objective and immediate assessment. Spot can screen up to 60 children per hour. That's enough to surpass the goals of even the most ambitious schools and charities.

Because Spot is fast and simple to use, it's easy to integrate into any charitable, nonprofit or school-based screening program. Spot goes where you go with a lightweight design and battery that powers a typical day of screening. You can do screenings virtually anywhere!

- Simple as taking a photo
- No formal training required
- Can be used by any staff member or volunteer

Screening Solution

With Spot, you get everything you need:

- User guide to help navigate you through Spot's usage
- WiFi enabled for easy printing to optional printer 
- Data can be easily stored and transferred to support various electronic records protocols.



Spot is Thorough

Spot instantly detects the indication of major vision issues, including:

- Anisometropia (unequal refractive power)
- Hyperopia (farsightedness)
- Strabismus (eye misalignment)
- Myopia (nearsightedness)
- Astigmatism (blurred vision)
- Anisocoria (pupil size deviations)

Spot's unique binocular screening for amblyopic precursors (anisometropia/eye misalignment) makes it the ideal tool for early intervention programs targeted at special needs, pre-verbal and pre-school children.

Spot Handheld Vision Screener

Spot handheld vision screener with power cord.

..... 400100



Spot Soft Carrying Case

Protective case to store and transport the Spot Vision Screener.

..... 400102

Replacement Power Cable for Spot

Power cable for the Spot Vision Screener.

..... 400103

1 Year Additional Warranty

Provide your Spot Vision Screener with an additional 1 year of protection.

..... 400104

Contact Good-Lite for information on recommended printer and router for Spot.



Vision Screening Kit



The AAPOS Vision Screening Kit is for screening the vision of individuals beginning at age 3 through adulthood. Spiral-bound, multisided, distance visual acuity eye charts that can be hung on a wall or held in the hands. The distance charts contains both Full-Chart Threshold and Single Critical Line formats with the following optotypes:

- **Sloan Letters** for when children know their letters.
- **LEA Symbols** for preschoolers and for children who do not yet know their letters.

The AAPOS Vision Screening Kit includes Black Handheld Flip Paddle Occluder, Children's Fun Frames Occluder Glasses, 10-foot cord for measuring, Translucent Response Panel, Instructions and DVD with additional information.

..... 120



AAPOS Supplemental Screening Kit

The AAPOS Vision Screening Kit: Supplemental Screening Package is for screening the vision of individuals beginning at age 3 through adulthood.

The near vision chart contains both Full-Chart Threshold and Single Critical Line formats with the following optotypes:

- **Sloan Letters** for when children know their letters.
- **LEA Symbols** for preschoolers and for children who do not yet know their letters.

The AAPOS Supplemental Vision Screening Kit includes Black Handheld Flip Paddle Occluder, Children's Fun Frames Occluder Glasses, 16 Plate Color Vision Test, Random Dot Butterfly stereopsis test, Instructions and DVD with additional information.


..... 130

Intermediate Distance Charts

80 cm and 66 cm Sloan Letters Intermediate Distance Charts

Charts calibrated to test visual acuity at an intermediate distance of 31" (80 cm) or 26" (66 cm). Charts are available in 3 scrambled versions for repetitive testing situations. Sloan Letters optotypes with proportionally spaced lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Notations for Letter Size, Decimal, and logMAR are also on the chart. Chart fits in Quantum Series, Model A+, and Model A illuminated cabinets. Chart size is 9" x 14" (22.9 cm x 35.6 cm). 31" (80 cm) testing distance.

80 cm Intermediate Distance Charts



- A** 80 cm Chart 1 758400
- B** 80 cm Chart 2 758500
- C** 80 cm Chart 3 758600

66 cm Intermediate Distance Charts



- D** 66 cm Chart 1 758410
- E** 66 cm Chart 2 758510
- F** 66 cm Chart 3 758610

Near Vision Charts for ESV1200 & Model A+



LEA Symbols Near Vision Chart

Proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Chart fits in Quantum Series, Model A+, and Model A illuminated cabinet. Chart size is 9" x 14" (22.9 cm x 35.6 cm). Includes response panel and flash cards.

..... 257800



LEA Numbers Near Vision Chart

Proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Chart fits in Quantum Series, Model A+, and Model A illuminated cabinet. Chart size is 9" x 14" (22.9 cm x 35.6 cm).

..... 273000



Sloan Letter Near Vision Charts

Proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Chart fits in Quantum Series, Model A+, and Model A illuminated cabinet. Chart size is 9" x 14" (22.9 cm x 35.6 cm).

- G** 756300
- H** 756400

Near Vision Acuity Cards for 16 inches/40 centimeters



LEA Symbols Near Vision Card (16 inches/40 centimeters)

This test assesses a child's functional vision at near distances. Can also be used to familiarize child with testing procedure before introducing a distance test.

Proportionally spaced (logMAR) lines on one side and more tightly-spaced symbols on the opposite side. Line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Response key printed on test card. Includes training cards and instructions. Card measures 8" x 10" (20.3 cm x 25.4 cm).

..... 250800



LEA Symbols Near Vision Screener (16 inches/40 centimeters)

This test is designed for amblyopia testing or screening. The background is gray with the exception of a white rectangle spaced at 100% surrounding each line of symbols. Line sizes from 20/400 to 20/16 (6/120 to 6/4.8) equivalent, 0.05 to 1.25. Printed on two sides. Response key is printed on the test card. Includes training cards and instructions. Screener measures 8" x 10" (20.3 cm x 25.4 cm).

..... 252000



LEA Symbols & LEA Numbers Near Card (16 inches/40 centimeters)

This card combines a LEA Symbols line test on one side and a LEA Numbers line test on the other. Proportionally spaced (logMAR) lines from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Response key is printed on the test card. Training cards and instructions included. Card measures 8" x 10" (20.3 cm x 25.4 cm).

..... **251000**



Sloan Letters Near Vision Card (16 inches/40 centimeters)

Proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Printed on both sides. 7" x 9" (17.8 cm x 22.9 cm).

..... **729000**



Landolt C and Tumbling E Near Vision Card (16 inches/40 centimeters)

Proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Printed on both sides. 7" x 9" (17.8 cm x 22.9 cm).

..... **728000**



LEA Numbers Near Vision Card (16 inches/40 centimeters)

Proportionally spaced (logMAR) lines on one side and more tightly-spaced numbers on the other. Line sizes from 20/400 to 20/10 (6/120 to 6/3) equivalent, .05 to 1.25. Card measures 8" x 10" (20.3 cm x 25.4 cm).

..... **270900**




Combination Near Vision and Reading Card

Near Vision Card has both continuous text reading examples and ETDRS format Sloan Letters. Sloan Letters line sizes range from 20/400 to 20/10 equivalent with continuous text line sizes ranging from 20/100 to 20/20 equivalent. Notations are indicated in Snellen in Feet/ Meters, logMAR, Decimal, Letter Size and Diopters of Add.

This near vision card is available in different languages.

Near vision card measures 10" x 8" (25.4 cm x 20.32 cm) and comes with attached 16" (40 cm) measuring cord for easy use. Printed on styrene for long lasting durability.

- English **684600**
- Spanish..... **658100**
- Portuguese **658200**
- French..... **658300**



Combination Near Vision and Reading Card with High Contrast and 10% Low Contrast

Near Vision Card has both continuous text reading examples and ETDRS format Sloan Letters. The front side of the card is 100% black optotypes while the reverse side is 10% low contrast level.

Sloan Letters line sizes range from 20/400 to 20/10 equivalent with continuous text line sizes ranging from 20/100 to 20/20 equivalent. Notations are indicated in Snellen in Feet/ Meters, logMAR, Decimal, Letter Size and Diopters of Add.

Near vision card measures 10" x 8" (25.4 cm x 20.32 cm) and comes with attached 16" (40 cm) measuring cord for easy use. Printed on styrene for long lasting durability.

- English **684700**
- Spanish..... **658500**
- Portuguese **658600**
- French..... **658700**



HOTV Near Vision Card (16 inches/40 centimeters)

This double-sided card combines Proportionally spaced (logMAR) lines on one side and 50% spaced rectangle (MassTest™ format) on the other. Line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00 Includes training cards and instructions. 8" x 10" (20.3 cm x 25.4 cm).

..... **731000**

Near Vision and Reading Visual Acuity


Near Vision Pocket Cards for 16"/40 cm



LEA Symbols Pocket-Sized Near Vision Card (16 inches/40 centimeters)

Use this small pocket card to quickly determine whether a child's visual acuity is within the normal range. Proportionally spaced (logMAR) lines on one side and 50% spaced symbols on the opposite side. Line sizes from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00 3.75" x 6.25" (9.5 cm x 15.9 cm).

..... 250950



LEA Symbols Pocket Box Stick (16 inches/40 centimeters)

The LEA Symbols Pocket Box Stick shows 20/125, 20/100, 20/80, 20/60, 20/50, 20/40, 20/30, and 20/20 lines. Each stick is 9/16" (1.4 cm) square and 4.5" (11.5 cm) long with fixation targets on each end.


- A** Proportionally spaced 486200
- B** Linear-spaced 507100



LEA Numbers Pocket-Sized Near Vision Card (16 inches/40 centimeters)

Use this small pocket card to quickly discover whether a child's or an adult's visual acuity is within the normal range. Proportionally spaced (logMAR) lines on one side and more tightly-spaced numbers on the other. Line sizes from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00 3.75" x 6.25" (9.5 cm x 15.9 cm).

..... 271000



Runge Pocket-Sized Near Vision Card (16 inches/40 centimeters)

Sloan Letters, LEA Numbers, or LEA Symbols Pocket Near Cards with line sizes ranging from 20/500 to 20/16 equivalent, 0.04 to 1.25 at 16" (40 cm). Near acuity measurements and their equivalents are listed on reverse side. 3.75" x 6.25" (9.5 cm x 15.9 cm).

- C** Sloan Letters 756000
- D** LEA Numbers 756100
- E** LEA Symbols 756200



Sloan Letters Pocket-Sized Near Vision Card (16 inches/40 centimeters)

Sloan letters with line sizes of 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00 on one side. Continuous text is printed on the reverse side with line sizes 20/80 to 20/16. 3.75" x 6.25" (9.5 cm x 16 cm).

..... 744000



Sloan Letters and Landolt C Near Vision Card

This easy-to-use visual acuity near vision card features Landolt C on one side and Sloan Letters on the other. Both sides have optotypes in proportionally spaced (logMAR) lines with line sizes ranging from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Card size is 3.75" x 6.25" (9.5 cm x 15.9 cm). Testing distance is 16" (40 cm).

..... 758300



LEA Symbols and Landolt C Near Vision Card

This easy-to-use visual acuity near vision card features Landolt C on one side and LEA Symbols pediatric symbols on the other. Both sides have optotypes in Proportionally spaced (logMAR) lines with line sizes ranging from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Card size is 3.75" x 6.25" (9.5 cm x 15.9 cm).

..... 257900



LEA Numbers and Landolt C Near Vision Card

This easy-to-use visual acuity near vision card features Landolt C on one side and LEA Numbers on the other. Both sides have optotypes in proportionally spaced (logMAR) lines with line sizes ranging from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Card size is 3.75" x 6.25" (9.5 cm x 15.9 cm).

..... 273100


Lighthouse Near Vision Cards from Optelec



Landolt C Near Vision Card

This easy-to-use visual acuity near vision card features Landolt C. Both sides have optotypes in proportionally spaced (logMAR) lines with line sizes ranging from 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00. Card size is 3.75" x 6.25" (9.5 cm x 15.9 cm).

..... 257900



Near Acuity Test Chart

This standardized assessment is designed with Sloan letters for testing subnormal low vision. Acuity range 16M to .5M (5/200 to 20/25). Measures 9" x 14" (23 cm x 35.5 cm)

..... 925400



"Bee" Continuous Text Card

This near vision reading card includes simple sentences in proportionally spaced (logMAR) lines with Snellen equivalent notation. Acuity range 8.0M to .4M (20/400 to 20/20). Pediatric card tests children with elementary reading levels (grades 1-3). Two-sided, measures 7" x 9" (18 cm x 23 cm).

..... 926100



Word "Game" Recognition Pocket-Size Test Card

Designed as a rapid test for checking current prescription and accuracy of the add to be prescribed. Line increments in logMAR units in Snellen equivalent notations. Acuity range 8.0M to .4M (20/400 to 20/20). Measures 3.5" x 6.5" (9.5 cm x 16 cm)

..... 925800



Number & Word Recognition Pocket-Size Test Card

Designed as a rapid test for checking current prescription and accuracy of the add to be prescribed. Contains words on one side and numbers on the reverse. Line increments in logMAR units in Snellen equivalent. Acuity range 8.0M to .4M (20/400 to 20/20). Measures 3.5" x 6.5" (9.5 cm x 16 cm)

..... 926000



"Bird" Continuous Text Card

Simple sentences in proportionally spaced (logMAR) lines with Snellen equivalent notations. Acuity range 8.0M to .4M (20/400 to 20/20). Pediatric card tests children with elementary reading levels (grades 1-3). Two-sided, measures 7" x 9" (18 cm x 23 cm).

..... 926400



Number Recognition Pocket-Size Test Card

Designed as a rapid test for checking current prescription and accuracy of the add to be prescribed. Line increments in logMAR units in Snellen equivalent notations. Acuity range 8.0M to .4M (20/400 to 20/20). Measures 3.5" x 6.5" (9.5 cm x 16 cm)

..... 925900

Near Vision and Reading Visual Acuity

Near Vision Reading Cards for 16"/40cm



Crowded HOTV Near Vision Card (16 inches/40 centimeters)

Crowded line sizes from 20/125 to 20/20 (6/38 to 6/6) equivalent, .16 to 1. Instructions printed on reverse. Includes paper isolating mask. 5.5" x 8" (14 cm x 20.3 cm).

..... 726000



Graded Reading Near Vision Booklet

This easy-to-use near vision reading book provides appropriate reading examples for preschool through 6th grade, according to Flesch-Kincaid grade levels. The preschool section provides preschool testing using the clinically tested LEA Symbols together with LEA Numbers, which provide a narrow acuity threshold. Each tab along the side of the book allows the user to quickly go to the grade to be tested.

Each Grade Level has examples for 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, and 20/20. The Reading Booklet also provides double-spaced examples for 2nd and 3rd grade students with suspected amblyopia, an Amsler Grid, an Acuity Equivalent Chart for 7 different optotype measurement nomenclatures, an Accommodation test, and an Astigmatic Chart.

Book has tabbed pages for ease of use and is printed on durable styrene. Coil binding permits the least distracting method of page presentation. Dimensions are 6.25" x 6.25" (16 cm x 16 cm) exclusive of tabs and binding.

..... 571100



Lebensohn Near Vision Card (14 inches/35.6 centimeters)

This 2-sided contemporary reading card, designed by James Lebensohn,

MD, provides a full range of near vision exercises, including words, numbers, and OX symbols. Optotypes include 20/800, 20/400, 20/260, 20/200, 20/100, 20/65, 20/50, 20/40, 20/30, 20/25, and 20/20. Notation equivalents include decimal readings in English and metric, Jaeger, point size, and visual efficiency percentages. Red/green bars employed on 5 pt to 12 pt optotypes. Card size is 5.5" x 8.5" (14 cm x 22 cm).

..... 119010



Fonda-Anderson Near Vision Reading Card (16 inches/40 centimeters)

One side displays continuous text while the other side has LEA

Symbols, LEA Numbers, and Sloan Letters. Card measures 5" x 8" (12.7 cm x 20.3 cm).

..... 727000



Contemporary Near Vision Card (16 inches/40 centimeters)

This 2-sided contemporary reading card provides familiar reading exercises, such as a crossword puzzle, help wanted ads, football and basketball scores, music song bar, calendar, playing cards, recipes, instructions, telephone white page listings, and TV listings. Continuous text examples are equivalent to 20/160, 20/100, 20/80, 20/50, 20/40 and 20/30. Notation equivalents include Jaeger and point size. Card size is 8.5" x 11" (21.5 cm x 28 cm).

..... 119000



Wormington Near Reading Card (16 inches/40 centimeters)

The Wormington Near Vision Card corrects inaccuracies of the Rosenbaum Vision Screener noted by Horton and Jones (1997) and others. For example, the "J" Notation, which lacked associated standards, was eliminated. The Wormington Near Vision card measures a patient's single optotype near vision acuity and reading acuity. Card features Times Roman font and Tumbling Es and provides a metric scale and pupil size array to enhance versatility.

..... 756075



Dr. Good's Near Vision Reading Card (14 inches/35.6 centimeters)

First-grade words. Line sizes from 20/160 to 20/20 (6/48 to 6/6) equivalent, 12 to 1. Card measures 5.5" x 8" (14 cm x 20.3 cm).

..... 724000




Sloan Letters Near Card (16 inches/40 centimeters)

Linear-spaced lines. Line sizes range from 20/400 to 20/16 (6/120 to 6/4.8) equivalent, 0.05 to 1.25. 5.5" x 8" (14 cm x 20.3 cm).

..... 725000

Spanish Near Vision and Reading Cards for 16"/40 cm



9x14 Spanish Near Reading Charts for Illuminated Cabinets

Continuous text reading cards in Spanish designed to fit in the ESV1200 or ESV1500 illuminated light box. Chart set includes 2 charts and available in two versions for repetitive testing.

Line sizes range from 20/320 to 20/250 on one chart and 20/200 to 20/16 equivalent. 16"/40 cm testing distance.

Chart Set 1 **756780**

Chart Set 2 **757670**



Spanish Contemporary Near Vision Card (16 inches/40 centimeters)

This card provides a variety of Spanish language reading examples in order to ensure that one is familiar to your patients. These include a song, train schedule, movie schedule, calendar, obituary, recipe, and soccer schedule. In addition, the chart includes reading examples for acuity levels at 20/30, 20/40, 20/50, 20/80, 20/100, and 20/160. Card is approximately . 8.5" x 10.5" (28 cm x 22 cm) and is laminated for durability.

. **547300**

Spanish Sloan Letters Near Pocket Card (16 inches/40 centimeters)

Sloan Letters with line sizes of 20/400 to 20/10 (6/120 to 6/3) equivalent, 0.05 to 2.00 on one side. Continuous text in Spanish is printed on the reverse side with line sizes 20/80 to 20/16. 3.75" x 6.25" (9.5 cm x 16 cm).


. **745000**



Spanish Wormington Near Reading Card (16 inches/40 centimeters)

This is a Spanish translation of the Wormington card. The known inaccuracies of the Rosenbaum Vision Screener as noted by Horton and Jones (1997) and others have been corrected in the Wormington card. The 'J' notation has been eliminated due to the lack of standards associated with it. This card measures a patient's single letter near acuity and reading acuity. Font used is Times Roman, typical of newspapers and books. Tumbling E series, metric scale, and pupil size array enhance the card's versatility.

. **756085**



Spanish Graded Reading Near Vision Booklet

This easy-to-use Spanish near vision reading book provides appropriate reading examples for preschool through 6th grade, according to Flesch-Kincaid grade levels. The Preschool section provides preschool testing using the clinically tested LEA Symbols together with LEA Numbers. Each tab along the side of the book allows the user to quickly go to the grade to be tested.

Each Grade Level has examples for 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, and 20/20. This reading booklet also provides double-spaced examples for 2nd and 3rd grade students with suspected amblyopia, an Amsler grid, an Acuity Equivalent Chart for 7 different optotype measurement nomenclatures, an Accommodation test, and an Astigmatic Chart.

Book has tabbed pages for ease of use and is printed on durable styrene. Coil binding permits the least distracting method of page presentation. Dimensions are 6.25" x 6.25" (16 cm x 16 cm) exclusive of tabs and binding.

. **592300**



Combination Near Vision and Reading Card

Near Vision Card has both continuous text reading examples and ETDRS format Sloan Letters. Sloan Letters line sizes range from 20/400 to 20/10 equivalent with continuous text line sizes ranging from 20/100 to 20/20 equivalent. Notations are indicated in Snellen in Feet/ Meters, logMAR, Decimal, Letter Size and Diopters of Add.

Near vision card measures 10" x 8" (25.4 cm x 20.32 cm) and comes with attached 16" (40 cm) measuring cord for easy use. Printed on styrene for long lasting durability.

Spanish Language. **658100**
with 10% Contrast Level on Back . . **658500**

Spanish Near Point Card (16 inches/40 centimeters)

This Spanish language continuous text near point reading card includes notation equivalents for Jaeger, letter size, font size, Snellen feet, and Snellen meters. Determines acuity levels at 20/25, 20/30, 20/40, 20/60, 20/100 and 20/200. 5" x 7.5" (13 cm x 18 cm) card is printed on white styrene plastic.

. **119760**



Miscellaneous Near Vision Items



LEA Symbols Playing Cards

These specially designed playing cards make it easy to measure near visual acuity in young children, while teaching concepts of similar/different, big/small, bigger/smaller. The 4 packs contain 16 cards each with symbols of varying sizes: 16M and 10M, 6.3M and 4M, 2.5M and 1.6M, and 1M and 0.63M. Includes training cards and instructions. The paper cards measure 1.6" x 2.75" (4 cm x 7 cm).

..... **252500**




LEA Symbols Domino Cards

Assess the near visual acuity of 2- and 3-year-old children, or use the cards at any age (older than 2+ years) to train an amblyopic eye. Set contains 3 packs of 16 domino-like cards with a large symbol at one end and a small symbol at the opposite end. Symbol sizes are 3.2M and 0.8M, 2.0M and 0.5M, and 1.25M and 0.3M. Each symbol on the card is printed with its M unit, decimal value, 6 meter, and 20 foot equivalency for visual acuity value at 40 cm testing distance. Includes training cards and instructions. Paper cards measure 1.6" x 2.75" (4 cm x 7 cm).

..... **251500**

Phoropter Roto Chart

This versatile near point Roto chart for the phoropter contains a variety of tests for the eye care professional. Fits on a phoropter using both the phoropter rod and clip. Notched wheel for easy chart selection. The Roto Chart comes in 2 versions and contains the following tests:



Phoropter Roto Chart 1 Features:

- Sloan Letters: Line sizes ranging from 20/200 to 12.5/20
- LEA Symbols: Line sizes ranging from 20/160 to 20/12.5
- LEA Numbers: Line sizes ranging from 20/160 to 20/12.5
- HOTV: Line sizes ranging from 20/160 to 20/12.5
- Tumbling E: Line sizes ranging from 20/200 to 20/12.5
- Astigmatic Wheel
- English and Spanish continuous text
- 20/32 and 20/63 Sloan Letters Vertical Line
- Crosshatch Target


..... **400400**



Roto Chart Clip

Used in conjunction with the Phoropter Rod to hold near point roto charts on your phoropter.

..... **560110**



Phoropter Roto Chart 2 Features:

- Sloan Letter; line sizes ranging from 20/200 to 12.5/20
- LEA Symbols; line sizes ranging from 20/160 to 12.5/20
- LEA Numbers; line sizes ranging from 20/160 to 12.5/20
- English and Spanish continuous text
- 20/40 Sloan Letter Vertical and Horizontal Line
- 20/20 Sloan Letter Horizontal Line
- Cross Hatch Target and more

..... **400450**



Phoropter Rod

Durable metal rod contains measurements for inches and diopters on one side and centimeters on the other. Phoropter Rod breaks down into two separate pieces for easy storage. 26.5" (67.31 cm) long in length.

..... **560100**

Low Vision Charts



LEA Symbols Chart for Vision Rehabilitation (40 inches/1 meter)

This low vision test is designed for examining children and adults with severe visual impairments. The largest symbols are 50M in size, corresponding to 1/50 (20/1000, 6/300, 0.02). The smallest symbols are 1M in size, corresponding to 1/1 (20/20, 6/6, 1.00) at a distance of 1 meter. The test comes with attached cord to maintain the proper testing distance. Includes response key, training cards, and instructions. Designed courtesy of August Colenbrander, MD. Chart measures 16" x 21.5" (40.7 cm x 54.6 cm). Folds to 16" x 10.75" (40.7 cm x 27.3 cm).


..... 258000



LEA Numbers Chart for Vision Rehabilitation (40 inches/1 meter)

This low vision test is designed for examining children and adults with severe visual impairments. The largest numbers are 50M in size corresponding to 1/50 (20/1000, 6/300, 0.02) and the smallest symbols are 1M in size corresponding to 1/1 (20/20, 6/6, 1.00) at a distance of 1 meter. The test comes with an attached cord to maintain the proper testing distance. Includes instructions. Designed courtesy of August Colenbrander, MD. 16" x 21.5" (40.7 cm x 54.6 cm). Folds to 16" x 10.75" (40.7 cm x 27.3 cm).

..... 272100



LEA Symbols Cards, 40M and 60M

These five LEA Symbols Cards can be used for presentation as single optotypes or arranged in a row as a line test. 60M cards are approximately 8.6" x 8.6" (21.9 cm x 21.9 cm). 40M cards are approximately 6" x 6" (15.25 cm x 15.25 cm).

..... 256100



LEA Numbers Cards, 40M and 60M

These five LEA Numbers Cards can be used for presentation as single optotypes or arranged in a row as a line test. 60M cards are approximately 6" x 9" (15.25 cm x 22.9 cm). 40M cards are approximately 4" x 6" (10.2 cm x 15.25 cm).

..... 272400



Low Vision Reading Card

Available in either English or German language near vision reading card with 16" (40 cm) attached measuring cord for easy measurement. Chart contains continuous text with line sizes ranging from 20/320 to 20/16 (6/95 to 6/4.8) equivalent. The card also contains Letter Size, Decimal, logMAR, and M Size notations. Chart measures 15" x 10" (25.5 cm x 38.5 cm).

English 777001
German 777000



Effective Magnification Acuity Chart

The EMAG Chart allows for easy determination of the starting dioptric power needed to read standard size print. The chart is used at 10" (25 cm) with a +4.00 add and gives the doctor or technician a uncomplicated approach for incorporating low vision in their practice. Acuity lines range from 20/400 to newspaper size print. The chart comes with instructions and has an attached measuring cord that ensures proper testing distance.

..... 777500



Runge Low Vision Near Card (16 inches/40 centimeters)

This easy-to-use low vision card was developed to mimic reading from left to right. Printed answers on back of card allows the tester to follow the patient's testing results. The back includes Snellen, meters, decimal, logMAR, font size, Jaegar, and M unit notations. A ruler and pupil gage is also included.

Sloan Letters with line sizes ranging from 20/1000 to 20/16 (equivalent). 13" x 7.5" (33 cm x 19 cm).

..... 756900

Low Vision



LEA Symbols Low Vision Book

LEA Symbols low vision test is designed for examining children and adults with severe visual impairments. This complete book includes testing for 10' (3 meters), 5' (1.5 meters) and for near vision at 16" (40 cm). Distance testing line sizes range from 20/1600 to 20/16 equivalent, 0.25 to 1.25. Near vision testing line sizes range from 20/800 to 20/16 equivalent, 0.25 to 1.25.

This book comes bound either on the left side or top. Kit includes Butterfly occluder glasses for children and metal frame occluder glasses for adults are included.

Left-Side Bound

Kit **527400**
 Book Only **575500**

Top Bound

Kit **527500**
 Book Only **575600**



LEA Numbers Low Vision Book

LEA Numbers low vision test is designed for examining children and adults with severe visual impairments. This complete book includes testing for 10' (3 meters), 5' (1.5 meters) and 16" (40 cm) near vision. Distance testing line sizes range from 20/1600 to 20/16 equivalent, 0.25 to 1.25. Near vision testing line sizes range from 20/800 to 20/16 equivalent, 0.25 to 1.25.

This book comes bound either on the left side or top. Kit includes Butterfly occluder glasses for children and metal frame occluder glasses for adults are included.

Left-Side Bound

Kit **513100**
 Book Only **575300**

Top Bound

Kit **527300**
 Book Only **575400**



Twin-Lux LED Stand Magnifier

The Schweizer 4.4X Twin-Lux LED Stand Magnifier is equipped with LEDs for bright, evenly dispersed, high contrast light. The Twin-Lux stand magnifier features the largest field of view at a 4.4x magnification. This is created using 2 aspheric lightweight lenses built-in at a precisely engineered distance. The Twin Lux magnifier enables a comfortable reading posture and utilizes the LEDs low energy consumption and long battery life. Requires 3 AA batteries (included).

Power: 4.4X
 Diopter: 8D + 6D
 Lens Size: 100 x 75mm + 100 x 50mm
 **582100**




531800 - 532300
 532410

Schweizer Sempral Round Hand Magnifiers

Schweizer Sempral Round Handheld Magnifiers are a lightweight and durable visual aid, with a sleek and elegant design. The hand magnifier features a high quality aspheric lens to provide a clear, quality image and a round grip that includes an integrated lanyard holder for ease-of-use. The 100 x 50 mm is the only magnifier offered with a rectangular lens.

Power	Diopter	Lens Size	Part Number
2.5X	6D	100 x 50mm	532410
3X	8D	85mm	532310
4X	12D	70mm	532210
5X	16D	60mm	532100
6X	20D	55mm	532000
8X	28D	35mm	531900
10.75X	39D	35mm	531800



Peripheral Field Loss



Combined Loss



Color Deficient

VisualEyes Vision Simulator Glasses

Simulator glasses are educational tools, developed by The Lighthouse Inc., to demonstrate the impact of different types of vision loss or near and distance activities. Experience what those with low vision see each day. The simulators in this kit allow you to try on a range of vision impairments and gain valuable understanding and insight into the daily challenges of those living with visual impairment.

Glasses are available in 7 versions with each one simulating a different visual impairment.

Combined Loss **200773**
 Central Field Loss **200774**
 Peripheral Field Loss **200775**
 Overall Blur **200776**
 Hemifield Loss. **200777**
 Low Contrast. **200778**
 Color Deficient. **200779**

Low Vision Magnifiers from SCHWEIZER



Schweizer Okolux Plus Handheld Magnifiers

The Schweizer Okolux Plus Mobil Handheld Magnifier features SMD LEDs (white clear light- 4500K light temperature). The SMD LEDs provides a bright, evenly dispersed, high contrast light. The aspheric lightweight magnifying lens is made of acrylic glass that is scratch-resistant. The lens provides a sharply defined, clear image that is further

enhanced by the SMD LEDs. Okolux Plus Mobil line of handheld magnifier's have been engineered for low energy consumption and a long-lasting battery life. The battery compartment was designed for easy battery replacement.

Power	Diopter	Lens Size	Head Type	Part Number
3.5	10D	75 x 50mm	Angled	582000
3.5X	10D	75 x 50mm	Straight	580900
4X	12D	70mm	Straight	580800
5X	16D	60mm	Straight	580700
6X	20D	55mm	Straight	580600
7X	24D	55mm	Straight	580500
8X	28D	35mm	Straight	580400
10.75X	39D	35mm	Straight	580300
13X	48D	35mm	Straight	580200
15X	56D	35mm	Straight	581900



Schweizer 3X/8D Okolux Plus SMD LED Stand Magnifier

The Schweizer Okolux Plus Stand Magnifier features SMD LEDs (white clear light- 4500K light temperature). The SMD LEDs provides a bright, evenly dispersed, high contrast light. Low energy consumption and a long-lasting battery life. The aspheric lightweight magnifying lens is made of acrylic glass that is scratch-resistant. The

lens provides a sharply defined, clear image that is further enhanced by the SMD LEDs. Okolux Plus Stand magnifiers feature easily interchangeable magnification heads, allowing you to create and update your personal magnifier if your future vision needs change or to replace a broken lens.

Power	Diopter	Lens Size	Part Number
3X	8D	100 x 75mm	581500
4X	12D	70mm	581400
5X	16D	60mm	581600
6X	20D	55mm	581300
7X	24D	55mm	581200
8X	28D	35mm	581100
10.75X	39D	35mm	581000
13X	48D	35mm	581800
15X	56D	35mm	581700

Low Vision



Schweizer 2.5X/6D Modular LED Stand Magnifier

This magnifier features three 4-chip SMD LEDs. It is a bright, consistent and high contrast light with the three precisely aligned LEDs. The Aspheric lens is lightweight and mounted on an open and adjustable stand that can be inclined to accommodate

right- or left-handed users and the lower powered magnifiers can be used when writing. This line of magnifiers also features interchangeable magnification heads.

Power	Diopter	Lens Size	Part Number
2.5X	6D	100mm	583110
3X	8D	100mm	583000
4X	12D	70mm	582900
5X	16D	60mm	582800
6X	20D	55mm	582700
8X	28D	35mm	582600



Schweizer Okolux Plus Low Power LED Illuminated Stand 4 Magnifier Kit

This Schweizer kit offers a range of magnifiers with 4 chip SMD LEDs packed in a compact tray box. They provide bright consistent and high-contrast light, the ideal illumination for reading or viewing objects. Each magnifier requires 3 AA batteries, included, and carry a 10-year limited warranty.

Kit includes:

- 581500 8D Stand Magnifier
- 581400 12D Stand Magnifier
- 581600 16D Stand Magnifier
- 581300 20D Stand Magnifier

..... **572600**



Optelec Low Vision Diagnostic Tool

Even if your specialty isn't low vision, this tool will allow you to easily guide more low vision patients to rehabilitation, support and vision aids. Assess, refer and retain patients and improve their quality of life in 10 minutes or less!

Assess - Conduct the test after an eye exam or it can be used in a home evaluation. It will determine the patient's reading ability and need for electronic or optical vision aids.

Refer - Write your recommendation on a Low Vision Passport and give it to your patient. This Educational Guide introduces various types of Optelec low vision aids and an informational end-user website. It helps them connect with local support networks, including Low Vision Rehabilitation



Specialists and Authorized Dealers. Build a partnership with professionals who will refer back to you!

Retain - Expand your patient care, partner for referrals, and distinguish your practice. Low vision patients get instant gratification and a new sense of hope despite their vision loss. This will help guide them towards a clear path for maintaining independence. Patients continue to return for check-ups and further treatment.

This kit includes:

- Diagnostic Tool
- Instruction Manual with Quick Start Guide
- Two-Sided Reading Test Card
- 15 Low Vision Passport Guides

..... **940100**



Schweizer Round Pocket Magnifier

The Schweizer Round Pocket Magnifier is a sleek, easy-to-use magnifier with Schweizer quality optics. The biconvex glass magnifying lens is mounted in transparent plastic. A strong grain leather casing protects the lens when not in use. The compact design allows you to store the pocket magnifier in a handbag, briefcase or pocket to take with you anywhere.

- 3.5X/10D - 60mm Lens**
- Black **544300**
- Red **578200**

- 6X/20D - 30 mm Lens**
- Black **544000**
- Red **543700**



Schweizer Bar Magnifier with Reading Line

The Schweizer bar magnifier or magnifying ruler is an ideal tool for magnifying text and following lines while reading. A handle is included on the reading magnifier to adjust while reading, and a flat bottom surface allows you to lay it flat on paper for hands-free use. An integral red-line guide helps to aid while you read. The bar magnifier is made of pure optical plastic.

- 1.5X Power**
- 3.5 cm x 20 cm Long **530200**
- 3.5 cm x 33 cm Long **530400**
- 2X Power**
- 2 cm x 16.5 cm Long **530000**



Optelec Compact Mini

The Optelec Compact Mini is the smallest and lightest electronic video magnifier offered by Optelec, designed as a simple and effective, high contrast alternative to a traditional optical magnifier. Read medicine labels, airline tickets, restaurant menus or price tags; the Compact Mini delivers convenience in a modern, pocket-sized solution.

The 3.5" full color screen has adjustable magnification of 2X, 5X, 8X, and 11X with multiple high contrast viewing options. Pocket-sized, lightweight at only 4.7 oz. 3-hour rechargeable battery; auto-shutoff to retain battery life.



Protective carrying case and wrist strap (included) with a 2-year standard warranty. **940000**



Optelec 4X/12D PowerDome Dome Magnifier

Optelec PowerDome 4X/12D Dome Magnifier is a great low vision aid for simple reading tasks requiring magnification. The dome magnifier is small and lightweight making them ideal for people having difficulty in holding a magnifier away from the page. The base easily glides over surfaces of reading material.

- Each PowerDome comes with a soft micro-fiber case for cleaning and protection.
- 65 mm. **076400**
- 80 mm. **076500**
- 95 mm. **076200**



PowerScope Handheld Monocular Telescopeonocular

The Optelec PowerScope Handheld Monocular Telescope is a compact and lightweight, achromatic mini monocular telescope. The focal barrel rotation is smooth and fluid, allowing you to easily set your ideal viewing distance.

Available in 3 different magnification levels: 2.5X25, 4X25, and 6X25.

Due to the new design of the optics, the PowerScope 2.5X25 monocular telescope can be used for both distance and near focus viewing.

The PowerScope Clip-On Monocular includes a protective case and detachable neck strap and comes with a Limited Lifetime Warranty.

- 2.5X25 **556300**
- 4X12 **037600**
- 6x12 **037700**



Optelec PowerTV Binocular Glasses

Lightweight, adjustable Optelec Power TV Glasses provide 2.1X magnification ideal for distance viewing. These comfortable glasses are perfect for watching television, movies, sporting events or other activities.

Using a fully adjustable twin-lens design to set a desired viewing distance, the easy-to-use glasses will maximize the size of your television screen and provide complete enjoyment of multiple activities. Optimal focal distance of 10'.

Weighs only 1.7 oz. Includes a micro-fiber cleaning cloth and protective carrying case. Limited One Year Warranty.

- 2.1X **460950**

Low Vision

Low Vision Lamps from



Ultra Slim Magnifying Lamp XR

With bright 28 watts daylight™ light, the unique 7"/17.5 cm XR™ lens and the new Free-Motion head-joint™, this lamp combines the latest

innovations in lighting design. These new features make your lamp extremely easy to operate and reliable. Our XR technology™ makes the lens 50% lighter and easy to clean. And the new head-joint guarantees you precise and smooth positioning of your lamp without constantly having to tighten the wing-nut.

..... 022080




MAGnificent Floor/ Table LED Magnifying Lamp

The MAGnificent Lamp is ideal for all types of hobbies, crafts and reading. By simply removing the pole, it converts from a floor to a table magnifying lamp. It provides 27

high performance, energy efficient daylight™ LEDs which only use 3 Watts of energy but provide 65 watts of traditional bulb light. The large 5.75"/14.6 cm semi-rimless 1.75X magnifying lens allows close-work activity.

..... 025050



Folds Up

Foldi™ LED Lamp

The portable LED Foldi™ Lamp features 30 bright daylight™ LEDs to give you comfortable full-spectrum natural daylight which is kind to your eyes, reducing glare which helps reduce eye strain, headaches and red-eye. The Foldi™ will give up to 8 hours of bright white light on batteries and unlimited power using the USB connection to your laptop, making it the ideal travel companion.

White..... 045000
Black..... 045001



Twist portable lamp

daylight's™ popular portable lamp now comes with an ingenious Easy-Twist™ shade for maximum light angle and maximum work space! The 'Easy Twist' shade has been designed by daylight™ to give

you the ultimate spread of light across the whole of your work area. The lamp itself is lightweight and designed to fold neatly away, making it easy to transport. Now for the technical bit - this top lamp boasts even more practical features, all designed to make sure you craft in comfort. The lamps detachable base makes it extremely stable and with daylight™ technology there is absolutely no strobe effect so eyestrain is reduced and potential headaches are eliminated.

..... 033700



Flexi-vision Floor Lamp

This product is the result of extensive research and collaboration with the RNIB and other low vision specialists. The Flexi-vision™ Floor Lamp regroups all the best features to make it the ultimate reading lamp. The 20 watts daylight™ energy saving bulb (285 watts equiv.) reduces glare and improves contrast

allowing you to read and write for a longer time in optimal comfort. With the new multi-swivel system and flexible arm, you can direct the light exactly where you need it and use this lamp near a sofa or at a table. The Flexi-vision™ Floor Lamp is also one of the safest lamps on the market with clear on/off switch location, anti-trip cable holders and low heat shade.

- Lamp brightness 8800 Lux at 15 cm and 2650 Lux at 30 cm
- Total amount of visible light emitted in Lumens is 1250 lm

..... 031067



The Swan Lamp

The EasyTwist™ shade with extra long flexi-arm neck reach of 14"/ 36 cm allows easy positioning of the powerful 18 watts energy saving light. The lamp also includes a crystal clear 1.75X magnification lens on a new longer removable 18"/46 cm flexi-arm.

Table Lamp 023020
Floor Lamp 023030



Hobby Lamp

The powerful 27 watts Full Spectrum Naturalight™ tube will give you the best color rendering and allow to read or work in comfort. With its flexible neck, you can also direct the light exactly where you are, not the other way around.

Table Lamp 001062
Floor Lamp 001072

Low Contrast Charts




Chart Size: 9" x 14", 22.9 cm x 35.6 cm
Recommended for use only with the Quantum™ Series Cabinets.

Translucent Low Contrast Charts (10 feet/3 meters)

A low contrast test is a special way to measure visual acuity by determining how well eyes function in low light situations and/or the ability to distinguish objects from similarly colored or shaded backgrounds. Testing is identical to the measurement of visual acuity at high contrast levels. The 2.5% test is the most practical test in clinical use. Available in 5 contrast levels: 1.25%, 2.5%, 5%, 10%, and 25%. LEA Symbols charts include response key and training cards. All contrast charts include recording forms and instructions.

A LEA Symbols*

1.25% (not shown)	253700
2.5% (shown)	253800
5% (not shown)	253900
10% (not shown)	254000
25% (not shown)	254100

*LEA Symbols charts come with both response panel and training cards.

B LEA Numbers

1.25% (not shown)	271500
2.5% (shown)	271600
5% (not shown)	271700
10% (not shown)	271800
25% (not shown)	271900

C Sloan Letters

1.25% (not shown)	600754
2.5% (shown)	600753
5% (not shown)	600752
10% (not shown)	600751
25% (not shown)	600750



Low Contrast Flip Charts

A low contrast flip chart is an easy way to measure, record, and detect changes in the transfer of visual information when the change affects only visual acuity at low contrast levels. Chart presents LEA Symbols, LEA Numbers, and Sloan Letters at the following contrast levels: black, 25%, 10%, 5%, 2.5%, and 1.25%. Pages are offset for easy flipping. LEA Symbols contrast flip charts include response key and training cards. All flip charts include recording forms and instructions. 6 pages, 3.25" x 9" (8.25 cm x 22.9 cm).

D LEA Symbols

.	251100
Mixed Order	251200

E LEA Numbers

.	270400
Mixed Order	270500

F Sloan Letters

.	736000
Mixed Order	737000



Hiding Heidi Low Contrast Face Test

Visual communication is the most important communication method during the first year of life. The ability to detect objects of low contrast is an important component of the visual system. For example, facial expressions are mediated by faint shadows and changes in the contours of the mouth and eyes. Determining the levels of contrast that an infant can detect helps provide baseline data

for evaluating future changes. For example, deviations from usual behavior may indicate disorders that leave vision at high-contrast levels unaffected.

Four cards printed on both sides in the following contrast levels: black, 25%, 10%, 5%, 2.5%, and 1.25%. Cards are 9" x 9" (23 cm x 23 cm). Includes instructions.

.	253500
-----------	---------------

Low Contrast Tests



Evans Letter Contrast Test (ELCT)

The Evans Letter Contrast Test (ELCT) chart is a standardized test for letter contrast sensitivity. The ELCT chart is designed to be used with either the ESV3000 or ESC2000 providing consistent lighting for improved testing results. The ELCT chart presents Sloan Letters optotypes line size of 20/630 with a testing distance of 1m. The varying contrast in sets of 3 letters of equal contrast level. All triplets are balanced for letter difficulty and differ in contrast by 0.15 log contrast step.

ELCT Chart 1 contains 12 contrast levels ranging from standard 100% acuity to 1.65 log contrast or 2.24% contrast level. Studies show that 1.65 log contrast is the average letter contrast sensitivity for the normal older population. (Weber contrast calculation).

ELCT Chart 2 presents 4 contrast levels from 1.8 to 2.25 log contrast.

ELCT charts are available as either Chart 1 only or as a set of both Chart 1 and Chart 2. Kit with charts and ESV3000 are also available that include ESV3000 Standardized Cabinet, remote control and floor stand.

- ELCT Charts Only**
- Chart 1 Only **500031**
 - Chart 2 Only **500032**
 - Set with Chart 1 and Chart 2 **500030**

- ESV3000 Kits**
- Chart 1 Only **500631**
 - Set with Chart 1 and Chart 2 **500630**



Adult Near Vision Low Contrast Test

This adult near vision low contrast test features both ETDRS format and Continuous Text in 100%, 25%, 10%, 5%, and 2.5% contrast levels. Continuous text letter sizes range from 20/100 to 20/20 (6/30 to 6/6) equivalent in each contrast level. ETDRS format letter sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent in 17 increments for each contrast level. Each contrast level page is tabbed for easy access with included instructions. Printed on durable long-lasting styrene plastic and clear coated for added protection. The near vision contrast booklet contains 7 pages measuring 6.25" x 6.25" (15.8 cm x 15.8 cm).


Testing distance is 16" (40 cm).
..... **587700**



Pediatric Low Contrast Near

Pediatric Near Low Contrast Test Booklet with contrast sensitivity levels of full contrast (100%), 25%, 10%, 5%, 2.5% and 1.25% proportionally spaced (logMAR) lines ranging from 20/400 to 20/16 using Michelson measures. One side of each page provides LEA Symbols for testing of 3 to 5 year olds. The reverse side shows LEA Numbers for testing 5 years olds and older. Near Contrast testing provides an efficient means to detect or track ocular pathologies, including glaucoma, ARMD, diabetic retinopathy, cataract, and keratoconus. Low contrast acuity loss measurement can be used to demonstrate progress in disease treatment and help children's parents understand the need for continued support.

Score card, recording form and comprehensive instructions are also provided to support plotting of a graph with the child's results indicating the relationship of visual acuity and contrast sensitivity. Normal ratios are also indicated to permit comparison. The booklet is printed on high-quality styrene plastic for long-term durability.
..... **610800**



Low Contrast Translucent ETDRS Charts from Good-Lite

Good-Lite offers a wide selection of low contrast charts to fit your testing needs. Charts are available in different contrast levels and testing distances.

All charts are proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00. 24.5"H x 25.5"W (62 cm x 65 cm).

- 13 feet/4 meters Testing Distance**
LEA Numbers Optotype
- 2.5% **274200**
 - 10% **274300**

- 13 feet/4 meters Testing Distance**
Sloan Letters Optotype
- 1.25% **500022**
 - 2.5% **500021**
 - 5% **500020**
 - 10% **500019**
 - 25% **500023**

- 10 feet/3 meters Testing Distance**
Sloan Letters Optotype
- 1.25% **500030**
 - 2.5% **500029**
 - 5% **500028**
 - 10% **500027**
 - 25% **500026**

Quantum Series Illuminated Cabinet and Charts



ESV1200 Quantum™ Cabinet

LEDs provide uniform illumination across the entire testing surface. Universal power supply automatically adjusts power. Can be powered by an optional battery pack. The ESV1200 has a compartment in back for chart




Only 3.25" (8.3 cm) thick

storage. Charts are not provided with this unit and must be ordered separately. 9.25"W x 15"H x 1.75"D (38.1 cm x 23.5 cm x 4.4 cm). CE and CSA listed.

..... **ESV1200**

Accessories



Stand for Professional, ESV1500, ESV1200, ESV1018, LED Insta-Line, Models A+, A, or M Cabinets

Add mobility to your light box with this deluxe universal 5-leg stand. Available with or without casters. **Note:** Using stand with LED Insta-Line requires LED Insta-Line bracket 914100 found on page 47 .

- A** Casters **705000**
- B** Without casters **704000**



Adjustable Wall Mount for ESV1500 and ESV1200 Illuminated Cabinets

This wall mount can easily adjust to patient's eye level for more accurate acuity testing or screening. Extends less than 4" (10.2 cm) from wall. This wall mount is for ESV1500 and ESV1200 illuminated cabinet only.

..... **706100**



Hinged Wall Bracket for Professional, ESV1500, ESV1200, ESV1018, Model A+, A, or M Cabinets

Allows cabinet to swing from wall or door jamb.

..... **700009**



Canvas Carrying Case

For Quantum Series, Model A+, Model A, LED Insta-Line, and Insta-Line PLUS™. Adjustable shoulder strap provides comfort when carrying an illuminated cabinet.

..... **904020**



ESV1500 Quantum™ Cabinet Calibrated for Research Vision Testing

Similar to the ESV3000 using advanced patented technology, self-calibrating to provided consistent photopic (85 cd/m²) and mesopic (3 cd/m²) illumination on the chart surface. Controlled by infrared wireless remote. The ESV1500 has a compartment in back for chart storage. Charts are not provided with this unit and must be ordered separately. 9.25"W x 15"H x 3.25"D (38.1 cm x 23.5 cm x 8.3 cm). CE and CSA listed.

..... **ESV1500**


Similar to the ESV3000 found on page 40 but Portable!

Illuminated Cabinets & Charts

10 feet/3 meters Distance Acuity Translucent Charts

ESV1500, ESV1200, Model A+, and Model A Translucent Charts

Chart Size: 9" x 14", 22.9 cm x 35.6 cm



A LEA Symbols Set or Chart (10 feet/3 meters)

Linear-spaced lines; line sizes range from 20/100 to 20/16 (6/30 to 6/3) equivalent, 0.20 to 1.25. Set includes response key, training cards, and instructions.

Set 255000
Chart Only 255001

B LEA Symbols Set or Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/8 (6/38 to 6/2.4) equivalent, 0.16 to 2.50. Set includes response key, training cards, and instructions.

Set 250400
Chart Only 250401

10' and 20' notations only (not shown)

Set 250411
Chart Only 250412

C LEA Symbols Set or Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/8 (6/38 to 6/2.4) equivalent, 0.16 to 2.50 with repetitive left eye/right eye groups. Includes response key, training cards, and instructions.

Set 252400
Chart Only 252401



D Crowded HOTV Set or Chart (10 feet/3 meters)

Crowded line sizes from 20/100 to 20/16 equivalent, 0.20 to 1.25. Set includes chart, response key, training cards, and instructions.

Set 600001
Chart Only 600724

E HOTV Set or Chart (10 feet/3 meters)

Linear-spaced lines; line sizes range from 20/100 to 20/16 equivalent, 0.20 to 1.25. Set includes chart, response key, training cards, and instructions.

Set 600017
Chart Only 600717

F HOTV Test Set or Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/8 (6/38 to 6/2.4) equivalent, 0.16 to 2.50. Set includes chart, response key, training cards, and instructions.

Set 600018
Chart Only 600736

G HOTV Test Set or Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/8 (6/38 to 6/2.4) equivalent, 0.16 to 2.50. Set includes chart, response key, training cards, and instructions.

Set 600759
Chart Only 600758

PLEASE NOTE: All charts may be held in the hands, mounted, or hung on a wall. Charts should be used at eye level with the person screened or examined.

Test Sets for LEA Symbols & HOTV

Test Sets for both LEA Symbols and HOTV optotypes include 1 chart, 1 response key with all four optotypes, and 4 training cards with one optotype on each card.

Response Panel and Cards are used to assist in preparing children for testing or help in identifying what optotype the children see on the chart.



LEA Symbols Response Key and Flash Cards



HOTV Response Key and Flash Cards

Replacement Response Panel and Flash Cards



A LEA Symbols Clear Response Key
 Response key with LEA Symbols. Printed on one side. Response Key is now translucent to allowing to easily view which optotype the child is pointing to.
 251701

B LEA Symbols Training Cards
 Training cards with LEA Symbols. Printed on one side.
 251800

LEA Symbols & HOTV Combination Response Key
 Response key with LEA Symbols on one side and HOTV on the other. Printed two sides.
 (not shown). 700030

LEA Symbols & HOTV Combination Training Cards
 Training cards with LEA Symbols on one side and HOTV on the other. Printed two sides.
 (not shown). 700031




C HOTV Response Key
 Response Key with HOTV optotypes. Printed on one side.
 700523

D HOTV Training Cards
 Training cards with HOTV optotypes. Printed on one side.
 700515

10 feet/3 meters Distance Acuity Translucent Charts

ESV1500, ESV1200, Model A+, and Model A Translucent Charts
 Chart Size: 9" x 14", 22.9 cm x 35.6 cm




LEA Numbers Translucent Distance Chart (10 feet/3 meters)
 Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/8 (6/38 to 6/2.4) equivalent, 0.16 to 2.50.
 271200



Tumbling E Chart (10 feet/3 meters)
 Linear-spaced lines; line sizes range from 20/100 to 20/16 (6/30 to 6/4.8) equivalent, 0.25 to 1.25. Lines are labeled with 10' and 20' notations only.
 600738




E Sloan Letters Chart (10 feet/3 meters)
 Linear-spaced line sizes; line sizes range from 20/100 to 20/16 (6/30 to 6/4.8) equivalent, 0.20 to 1.25. Lines are labeled with 10' and 20' notations only.
 This chart is recommended by the Texas Department of Public Health for vision testing and screening.
 600723



Tumbling E Chart (10 feet/3 meters)
 Wide-spaced lines; line sizes range from 20/100 to 20/16 (6/30 to 6/4.8) equivalent, 0.20 to 1.25. Lines are labeled with 10' and 20' notations only.
 This chart is recommended by the Texas Department of Public Health for vision screening.
 600722

F Sloan Letters Chart (10 feet/3 meters)
 Proportionally spaced (logMAR) lines; line sizes range from 20/160 to 20/10 (6/48 to 6/3) equivalent, 0.12 to 2.00.
 600705
 Sloan Letters Chart with 10' and 20' notations only (not shown).
 600711

13 ft/4 m Translucent Distance Charts



LEA Symbols Set or Chart (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/100 to 20/8 (6/30 to 6/2.4) equivalent, 0.20 to 2.50. Set includes response key, training cards, and instructions.

- Set 255100
- Chart Only 255101



A LEA Numbers Translucent Distance Chart (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/100 to 20/8 (6/30 to 6/2.4) equivalent, 0.20 to 2.50.

- 271300

B Sloan Letters Chart (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/125 to 20/10 (6/38 to 6/3) equivalent, 0.16 to 2.00.

- 600739

20 feet/6 meters Translucent Distance Charts



HOTV Test Set or Chart (20 feet/6 meters)

Crowded lines; line sizes ranging from 20/80 to 20/16 (6/24 to 6/4.8), 0.25 to 1.25. Set includes chart, response key, training cards, and instructions.

- Set 600006
- Chart Only 600731



Sloan Letters Chart (20 feet/6 meters)

Linear-spaced lines; line sizes range from 20/80 to 20/16 (6/24 to 6/4.8), 0.25 to 1.25. Lines are labeled with 20' notations only.


- 600725



Sloan Letters Chart (20 feet/6 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/80 to 20/10 (6/24 to 6/4.8), 0.25 to 2.00. Lines are labeled with 20' notations only.

- 600737



Tumbling E Chart (20 feet/6 meters)

Linear-spaced lines; line sizes range from 20/80 to 20/16 (6/24 to 6/4.8), 0.25 to 1.25. Lines are labeled with 20' notations only.

- 600735



Tumbling E Chart (20 feet/6 meters)

Wide-spaced lines; line sizes range from 20/100 to 20/16 (6/24 to 6/4.8), 0.20 to 1.25. Lines are labeled with 20' notations only.

- 600728



Snellen Letters Chart (20 feet/6 meters)

Linear-spaced lines; line sizes range from 20/100 to 20/16 (6/30 to 6/4.8), 0.20 to 1.25. Lines are labeled with 20' notations only.


- 600727

Intermediate Distance Charts

80 cm and 66 cm Sloan Letters Intermediate Distance Charts

Charts calibrated to test visual acuity at an intermediate distance of 31" (80 cm) or 26" (66 cm). Charts are available in 3 different scrambled versions for repetitive testing situations. Sloan Letters optotype with proportionally spaced lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Notations for Letter Size, Decimal, and logMAR are also on the chart. Chart fits in Quantum Series, Model A+, and Model A illuminated cabinets. Chart size is 9" x 14" (22.9 cm x 35.6 cm). 31" (80 cm) testing distance.

80 cm Intermediate Distance Charts



- A** 80 cm Chart 1 **758400**
- B** 80 cm Chart 2 **758500**
- C** 80 cm Chart 3 **758600**

66 cm Intermediate Distance Charts



- D** 66 cm Chart 1 **758410**
- E** 66 cm Chart 2 **758510**
- F** 66 cm Chart 3 **758610**

16 inches/40 centimeters Near Vision Charts



LEA Symbols Near Vision Charts

Proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Chart fits in Quantum Series, Model A+, and Model A illuminated cabinets. Chart size is 9" x 14" (22.9 cm x 35.6 cm). 16" (40 cm) testing distance.


..... **257800**



LEA Numbers Near Vision Charts

Proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Chart fits in Quantum Series, Model A+, and Model A illuminated cabinets. Chart size is 9" x 14" (22.9 cm x 35.6 cm). 16" (40 cm) testing distance.

..... **273000**



Sloan Letters Near Vision Charts

Proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Chart fits in Quantum Series, Model A+, and Model A illuminated cabinets. Chart size is 9" x 14" (22.9 cm x 35.6 cm). 16" (40 cm) testing distance.

- G** **756300**
- H** **756400**

Spanish Near Vision Reading Charts

Spanish Near Vision Reading chart for Illuminated Cabinet

A set of 2 near vision chart with continuous text in Spanish language. Charts fit in either the Model A+ or ESV1200 illuminated cabinet. Three different chart combination versions are available when repetitive testing is involved. Line sizes range from 20/320 to 20/250 (6/60 to 6/4.8) equivalent, 0.10 to 1.25 on one chart and 20/200 to 20/16 (6/60 to 6/4.8) equivalent, 0.10 to 1.25 on the other.

Charts come in 2 different versions incase repetitive testing is involved.

Testing distance is 16" (40 cm).



Chart Set 1 **756780**



Chart Set 2 **756760**

ESV1018 Professional Illuminated Cabinet & Charts

10 feet/3 meters ESV1018 Translucent Distance Charts

Chart Size: 10" x 18", 25.4 cm x 45.7 cm



ESV1018 LED Professional Illuminated Cabinet

The ESV1018 Professional Illuminated Cabinet holds 10" x 18" (25.4 cm x 45.7 cm) charts and uses light emitting diodes (LEDs) to provide uniform illumination for standardized testing. Controlled by included wireless remote. Universal power supply. Measures 18.125"H x 10.125"W x 3"D (46 cm x 25.7 cm x 7.6 cm). CE and CSA Listed. Does not include charts, purchased separately.


.....**ESV1018**
Rechargeable Battery Pack (Optional)
 Rechargeable battery pack for all LED series illuminated cabinets.
 **914003**



Adjustable Wall Mount for ESV1018 Cabinets

This wall mount can easily adjust to patient's eye level for more accurate acuity testing or screening. Extends less than 4" (10.2 cm) from wall. This wall mount is for ESV1018 illuminated cabinet only.

..... **706000**



LEA Symbols Set or Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Set includes response key, training cards, and instructions.

Set **800722**
 Chart Only **800721**



Crowded HOTV Translucent Chart (10 feet/3 meters)

Crowded lines; line sizes range from 20/160 to 20/16 (6/48 to 6/4.8) equivalent, 0.12 to 1.25. Set includes chart, response key, training cards, and instructions.

Set **800004**
 Chart Only **800726**

PLEASE NOTE: All charts may also be mounted or hung on a wall.



LEA Numbers Translucent Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/160 to 20/8 (6/48 to 6/2.4) equivalent, 0.12 to 2.50.


..... **800737**



A Sloan Letters Chart (10 feet/3 meters)

Linear-spaced lines; line sizes range from 20/160 to 20/16 (6/48 to 6/4.8) equivalent, 0.12 to 1.25.


..... **800713**



B Sloan Letters Translucent Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/160 to 20/8 (6/48 to 6/2.4) equivalent, 0.12 to 2.50.

..... **800735**



**Tumbling E Chart
(10 feet/3 meters)**

Linear-spaced lines; line sizes range from 20/160 to 20/16 (6/48 to 6/4.8) equivalent, 0.12 to 1.25.

..... 800714


13 feet/4 meters ESV1018 Translucent Distance Charts



**Sloan Letters Translucent Chart
(13 feet/4 meters)**

Proportionally spaced (logMAR) lines; line sizes range from 20/160 to 20/10 (6/48 to 6/3) equivalent, 0.12 to 2.00.

..... 800738




Low Vision Sloan Letters Translucent Chart (20 feet/6 meters)

Low vision reversed Sloan Letters chart; line sizes range from 20/400 to 20/160 (6/120 to 6/48), 0.05 to 0.125. Available as either standard or reversed letters.

A..... 800741
B (Reversed Letters) 800742

20 feet/6 meters ESV1018 Translucent Distance Charts



**LEA Symbols Set or Chart
(20 feet/6 meters)**

Proportionally spaced (logMAR) lines; line sizes range from 20/80 to 20/8 (6/24 to 6/2.4), 0.25 to 2.50. Set includes response key, training cards, and instructions.

Set 800724


Chart Only 800723



**HOTV Translucent Chart
(20 feet/6 meters)**

Linear-spaced lines; line sizes ranging from 20/100 to 20/16 (6/30 to 6/4.8), 0.20 to 1.25. Response key and training cards are not included.

..... 800704



**LEA Numbers Translucent Chart
(20 feet/6 meters)**

Proportionally spaced (logMAR) lines; line sizes range from 20/80 to 20/8 (6/24 to 6/2.4) equivalent, 0.25 to 2.50.

..... 800701



**C Sloan Letters Chart
(20 feet/6 meters)**

Linear-spaced lines; line sizes range from 20/100 to 20/16 (6/30 to 6/4.8), 0.20 to 1.25.

..... 800733



**D Sloan Letters Translucent Chart
(20 feet/6 meters)**

Proportionally spaced (logMAR) lines; line sizes range from 20/100 to 20/10 (6/30 to 6/3), 0.20 to 2.00.


..... 800740



E Sloan Letters Reversed Translucent Chart (20 feet/6 meters)

Reversed proportionally spaced (logMAR) lines; line sizes range from 20/100 to 20/16 (6/30 to 6/4.8), 0.20 to 1.25.

..... 800739




**F Tumbling E Chart
(20 feet/6 meters)**

Linear-spaced lines; line sizes range from 20/100 to 20/20 (6/30 to 6/6), 0.20 to 1.00.

..... 800729

ETDRS Illuminated Cabinet & Charts



Simplify your ETDRS testing with either the ESV3000 or ESC2000

ESV3000 ETDRS Illuminated Cabinet

Contains patented internal calibration technology, which automatically controls the standardized LEDs (light emitting diodes) output. When this lighting technology is combined with a universal standard power supply, the result is consistent standardized lighting that makes the ESV3000 the most accurate, easy-to-use tester for evaluating ETDRS and logMAR acuity. Simply turn on the device and with a push of a button, the ESV3000 automatically calibrates to a photopic light level of 85 cd/m² or a mesopic level of 3 cd/m². These light levels are

recommended by the National Academy of Sciences' Committee on Vision Testing Standards.

The ESC2000 Illuminator Cabinet

This is our more cost effective ETDRS Cabinet providing a uniform retro-illuminated testing surface, ideal for clinical testing of high contrast and low contrast logMAR tests, including ETDRS, Lea Symbols, Lea Numbers, Landolt C, etc. This makes it the perfect choice when precision calibrated luminance is not required.

Both measure 24.5"H x 25.5"W x 3.75"D (62 cm x 65 cm x 9.75 cm). CE and CSA listed.

Compare Cabinets	ESV3000 with Calibrated LED Lights	ESC2000 with LED Lights
Self Calibrate to 85 cd/m ²	✓ Yes	✗ No
Provide Mesopic Testing	✓ Yes	✓ Yes [†]
Provide Uniform Light	✓ Yes	✓ Yes
Ready for Immediate Use	✓ Yes	✓ Yes
Universal Power Supply	✓ Yes	✓ Yes
Custom Light Levels Available	✓ Yes	✓ Yes
Built-in Glare Light Receptacle	✓ Yes	✗ No
Remote Controlled	✓ Yes	✗ No

ESV3000 ETDRS Standardized Viewer

The only ETDRS viewing cabinet that self-calibrates illumination on the chart surface for both photopic (85 cd/m²) and mesopic (3 cd/m²) light level testing. The unit utilizes the most advanced self-monitoring LED light source patented technology, providing uniform and highly calibrated illumination across the entire testing surface (*Patent 5,078,486*).

ESV3000 Kit With Original Series 4 Meter ETDRS Charts

Includes the ESV3000 viewer, ESV3000 stand, and 3 Original Series 4 meter ETDRS charts.*

..... 500600

ESV3000 Kit With "2000" Series 4 Meter ETDRS Charts

Includes the ESV3000 viewer, ESV3000 stand, and 3 "2000" Series 4 meter ETDRS charts.*

..... 500601

ESV3000 Kit With "2000" Series 3 Meter ETDRS Charts

Includes the ESV3000 viewer, ESV3000 stand, and 3 "2000" Series 3 meter ETDRS charts.*

..... 500602

ESV3000 ETDRS Viewer Individual Items

ESV3000 ETDRS Standardized viewer only.*

.....ESV3000

ESC2000 ETDRS Standardized Viewer

Our more cost effective ETDRS Viewer when calibrated light levels are not required.

ESC2000 Kit With Original Series 4 Meter ETDRS Charts

Includes the ESC2000 viewer, ESC2000 stand, and 3 Original Series 4 meter ETDRS charts.*

..... 500700

ESC2000 Kit With "2000" Series 4 Meter ETDRS Charts

Includes the ESC2000 viewer, ESC2000 stand, and 3 "2000" Series 4 meter ETDRS charts.*

..... 500701

ESC2000 Kit With "2000" Series 3 Meter ETDRS Charts

Includes the ESC2000 viewer, ESC2000 stand, and 3 "2000" Series 3 meter ETDRS charts.*

..... 500702

ESC2000 ETDRS Viewer Individual Items

ESC2000 ETDRS Standardized viewer only.*

..... ESC2000

[†] Manual user adjustment required.

* Power plug adaptor may be needed for outlets outside the United States.

Accessories and Options



ESV3000 Glare Lights

Easy-to-use Glare Lights simply plug in. Turn on and off remotely when used in conjunction with the ESV3000. Set of 2. Can only be used with the ESV3000.

..... 500100




ESV3000 Floor Stand

ESV3000 Floor Stand with casters. The base of the stand measures 27.5" square (69.85 cm). The height is adjustable from a minimum of 33" (83.82 cm) to maximum of 49" high (124.46 cm).

..... 500004

13 feet/4 meters ETDRS Distance Charts


Chart Size: 24.5" x 25.5", 62 cm x 65 cm



"2000" Series ETDRS Translucent Distance Charts (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00. Revised version of the original ETDRS series based on more recent research.


- A Chart "1" 500016
- B Chart "2" 500017
- C Chart "3" 500018



Original Series ETDRS Translucent Distance Charts (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00. Standardized format is used in National Eye Institute-sponsored clinical studies and in low vision evaluations when accurate and repeatable visual acuity measurements are needed.


- D Chart "R" 500013
- E Chart "1" 500014
- F Chart "2" 500015



LEA Symbols ETDRS Translucent Distance Charts (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/160 to 20/8 (6/48 to 6/2.4) equivalent, 0.90 to 2.50. Includes response key, training cards, and instructions.

..... 253333



LEA Numbers ETDRS Translucent Distance Chart (13 feet/4 meters)

Proportionally spaced (logMAR) lines; lines range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Chart is available in 3 scrambled versions for repetitive testing.


- G 270200
- H 270200A
- I 270200B

Note: All charts, except for the low contrast charts, can be hung on a wall.

Illuminated Cabinets & Charts

13 feet/4 meters ETDRS Distance Charts

Chart Size: 24.5" x 25.5", 62 cm x 65 cm



Landolt C ETDRS Translucent Distance Charts (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00.


..... 500037




Tumbling E ETDRS Translucent Distance Charts (13 feet/4 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00.

..... 500044



New Roll-up ETDRS Charts can be found on page 44.




HOTV ETDRS Translucent Distance Chart (13 feet/4 meters)

Proportionally spaced (logMAR) lines; lines range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Chart is available in 2 scrambled versions for repetitive testing.

A 500025

B 500034




Crowded HOTV ETDRS Translucent Distance Chart (13 feet/4 meters)

HOTV letters with 50% contour interaction bars to provide single optotypes with crowding. Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Chart is available in 2 scrambled versions for repetitive testing.

C 500035

D 500036

10 feet/3 meters ETDRS Distance Charts




"2000" Series ETDRS Translucent Distance Charts (10 feet/3 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00. Revised version of the original ETDRS series based on more recent research.

E Chart "1" 500010

F Chart "2" 500011

G Chart "3" 500012



Original Series ETDRS Translucent Distance Charts (10 feet/3 meters)


Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00. Revised version of the original ETDRS series based on more recent research.

H Chart "R" 500046

I Chart "1" 500047

J Chart "2" 500048

Note: All charts, except for the low contrast charts, can be hung on a wall.



LEA Symbols ETDRS Translucent Distance Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; lines range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50. Includes response key, training cards, and instructions. 24.5" x 25.5" (62 cm x 65 cm).

..... 250300



LEA Numbers ETDRS Translucent Distance Chart (10 feet/3 meters)

Proportionally spaced (logMAR) lines; lines range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50.

..... 271400

6.5 feet/2 meters ETDRS Distance Charts



"2000" Series ETDRS Translucent Distance Charts (6.5 feet/2 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00. Revised version of the original ETDRS series based on more recent research.

- K** Chart "1" 500041
- L** Chart "2" 500042
- M** Chart "3" 500043





Original Series ETDRS Translucent Distance Charts (6.5 feet/2 meters)

Proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00. Revised version of the original ETDRS series based on more recent research.

- N** Chart "R" 500038
- O** Chart "1" 500039
- P** Chart "2" 500040

ETDRS Low Contrast Charts



Low Contrast Translucent ETDRS Charts from Good-Lite

Good-Lite offers a wide selection of low charts to fit your testing needs. Charts are available at different contrast levels and testing distance.

Some charts are available in "Scrambled Versions" with optotypes in different order for when repetitive testing is involved.

All charts are proportionally spaced (logMAR) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent, 0.10 to 2.00. 24.5"H x 25.5"W (62 cm x 65 cm).

13 feet/4 meters Testing Distance

Sloan Letters Optotype

- 1.25% 500022
- 1.25% - Scrambled Version 500033
- 2.5% 500021
- 2.5% - Scrambled Version 500032
- 5% 500020
- 10% 500019
- 25% 500023

10 feet/3 meters Testing Distance

Sloan Letters Optotype


- 1.25% 500030
- 2.5% 500029
- 5% 500028
- 10% 500027
- 25% 500026

13 feet/4 meters Testing Distance

LEA Numbers Optotype

- 2.5% 274200
- 10% 274300

ETDRS Rollup Charts



ETDRS Translucent Distance Charts

Our standard ETDRS charts printed on a thinner styrene plastic allowing the chart to be easily rolled up for storage or transporting. Each chart comes with plastic wall hanging clip.

Each chart is proportionally spaced (logMAR) lines; lines range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent, 0.10 to 2.50.

13ft/4m ETDRS Charts

Original Series Chart R **511113**

Original Series Chart 1 **511114**

LEA Numbers Chart **270211**





10ft/3m ETDRS Charts

LEA Symbols ETDRS Chart **250311**

Super Pinhole Cabinet & Charts



now using LED



Super Pinhole Macula Charts (5 feet/1.5 meters)

High contrast charts for the Super Pinhole Macula Cabinet. Optotypes are Proportionally spaced (logMAR) lines; line sizes ranging from 20/100 to 20/20 (6/30 to 6/6) equivalent, 0.25 to 1.25.

F Sloan Letters **300304**

G Sloan Letters Scrambled **300305**

H LEA Numbers **300306**

I LEA Symbols **300307**

Replacement Extension Cord

Replacement extension cord for the Super Pinhole Macula Cabinet. **300320**

Super Pinhole Macula Cabinet

This Cabinet can assist with testing potential visual acuity when considering cataract surgery for patients with retinal disease, such as macular degeneration. This greatly improved Super Pinhole Macula Cabinet provides uniform illumination across the entire chart. LEDs provide uniform illumination. Universal power supply, Sloan Letters Super Pinhole chart (300304), Extension Cord (300320), and Multiple Pinhole Occluder (700010) are included with the cabinet.

. **300300**

Accessories



Multiple Pinhole Spectacle Occluder

Pinhole Spectacle . . . **700010**

Spectacle Only **700012**



Reversible Pinhole Occluding Glasses

One side is occluded with 5 pinholes, the other is fully occluded. Sturdy metal frames with straight temples for easy reversibility.

. **204006**



Small Pinhole Occluder

Miniature occluder with 17 pinholes. 4" x 1.25

. **120400**



J Short-handled black occluder **120100**

K Short-handled single pinhole occluder **120301**

L Short-handled 48 pinhole occluder **120300**

VectorVision CSV-1000 Contrast Sensitivity



Remote Control



The New CSV-1000HGT

The CSV-1000HGT tests a combination of contrast and glare sensitivity providing the most accurate measure to evaluate cataracts, glaucoma, diabetic eye disease, contact lens performance, refractive surgery, etc. Patented photocell circuitry automatically monitors and calibrates instrument testing light levels. The CSV-1000HGT includes the following: CSV-1000HGT instrument, remote control (with battery), score pads, your choice of test face(s), and instruction manual.

CSV-1000HGT Contrast Testing System

CSV-1000HGT Instrument with Glare Lights and 1 Test Face. Stand sold separately.

(shown) **500050**

Accessories for CSV-1000

Mesopic Testing Lens

..... **500053**

Floor Stand with Casters

..... **500004**

CSV-1000 Contrast Testing System

The CSV-1000 is the only self-standardized vision testing instrument. In each row of the device, photocell circuitry continually monitors and calibrates the instrument light level. This unique system senses both the external light striking the face of the instrument and the internal light which emanates from the miniature fluorescent bulbs. By measuring both light sources, the patented circuitry (U.S. Patent Number 5,078,486) determines the total light reflected from the instrument face and provides for a constant test luminance level of 85 cd/m².

CSV-1000 Contrast Testing System comes with CSV-1000 unit and one contrast testing chart. Floor stand sold separately.


..... **500000**



CSV-1000 ETDRS Test Face

ETDRS acuity testing is the worldwide standard for visual acuity testing. The CSV-1000 ETDRS

test face provides standardized testing of ETDRS acuity. Three different ETDRS tests are available: Chart R, Chart 1, and Chart 2. Chart R **500003**



CSV-1000S Test Face

The CSV-1000S test face is widely used for in- or out-of-office screenings and for cataract evaluations.

This test presents standard visual acuity from 20/15 to 20/200, 2 rows of spatial frequencies (6 and 12 cycles/degree), and a real-world driver's scene. The test is useful for cataract documentation because it provides a standard acuity score, a functional acuity score, and a real-world simulation.


..... **500001**



CSV1000CA Contrast Acuity (CA) Chart

Tests high contrast and 10% contrast from 20/100 to 20/10.

..... **500052**



CSV1000E Contrast Chart

This test provides for 4 rows of sine-wave gratings. At the recommended test distance of 8 ft (2.5 meters), these

gratings test the spatial frequencies of 3, 6, 12, and 18 cycles/degree. The CSV-1000E provides a full contrast sensitivity curve, which is useful for the evaluation of ocular disease (particularly cataracts, glaucoma, optic neuritis, diabetes, and macular degeneration), contact lenses, and refractive surgery. For cataract evaluation/documentation, a functional acuity score can be obtained directly from the test scoring sheet.

..... **500051**

LED Insta-Line Quantum Cabinet & Charts



Complete LED Insta-Line Quantum™

All LED Insta-Line Quantum cabinets use light emitting diodes (LEDs) to provide uniform illumination for standardized testing. The increased efficiency of the bulbs and the rechargeable battery (included) permit testing in almost any location. Controlled by radio frequency (RF) wireless remote. A power pack is included for recharging and using the unit with or without the battery.

One person can administer a complete series of vision tests in less than 2 minutes. Ideal for ages 3 years through adulthood. Used for testing visual acuity at 10 feet and for assessing muscle imbalance. Includes LEA Symbols® Test Set, LEA Numbers®, Tumbling E Chart, Sloan Letters Chart, and HOTV Test Set, occluder, +1.75, +2.25, and +2.50 hyperopia glasses. Completely self-contained and weighs only 10 lbs. Universal power supply. 14.5" H x 9.25" W x 5" D. CSA Listed.

..... **915000**

Preschool LED Insta-Line Quantum™

Complete and accurate screening for preliterate children as young as 3 years. Portable, durable, and works with the touch of a finger. Includes LEA Symbols test set, HOTV test set, and occluder. Universal power supply. CSA Listed.

..... **916000**

Basic LED Insta-Line Quantum™

Cabinet only, charts not included. Universal power supply. CSA Listed.

..... **914000**

10 feet/3 meters Testing Distance Charts and Sets for ALL Insta-Line Cabinets



LEA Symbols Test Set

Linear-spaced optotypes in the following size equivalents: 20/60, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, and 20/12.5. Two lines of each size with 5 optotypes on each line. Includes chart, response key, training cards, switch box overlay, and instructions. This test set is included with 915000 Complete LED Insta-Line Quantum™ and 916000 LED Preschool Insta-Line Quantum™. 9" x 14" (22.9 cm x 35.6 cm).

..... **904025**

LEA Symbols Chart Only

Includes switch box overlay. 9" x 14" (22.9 cm x 35.6 cm).

..... **904032**



LEA Numbers Chart

Linear-spaced optotypes in the following size equivalents: 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, and 20/12.5. Two lines of each size with 5 optotypes on each line. Includes chart and switch box overlay. Included with 915000 Complete LED Insta-Line Quantum™. 9" x 14" (22.9 cm x 35.6 cm).

..... **904026**



HOTV Test Set

Linear-spaced optotypes in the following size equivalents: 20/60, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, and 20/12.5. Two lines of each size with 5 optotypes on each line. Includes chart, response key, training cards, switch box overlay, and instructions. This test set is included with 915000 Complete LED Insta-Line Quantum™ and 916000 LED Preschool Insta-Line Quantum™. 9" x 14" (22.9 cm x 35.6 cm).

..... **904030**

HOTV Chart Only

Includes switch box overlay. 9" x 14" (22.9 cm x 35.6 cm).

..... **904031**



Sloan Letters Chart

Linear-spaced optotypes in the following size equivalents: 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, and 20/12.5. Two lines of each size with 5 optotypes on each line. Includes chart and switch box overlay. This chart is included with 915000 Complete LED Insta-Line Quantum™. 9" x 14" (22.9 cm x 35.6 cm).

..... **904029**



LEA Symbols Combination Distance and Near Test Set

Easily screen both near and distance vision with one chart using the Insta-Line illuminated cabinet.

This pediatric eye chart features LEA Symbols pediatric symbols with 2 lines of linear-spaced optotypes in the following size equivalents: 20/40, 20/32, 20/25 and 20/20. Distance screening for these lines is 10 ft, 3 meters.

Near vision screening features LEA Symbols with 2 lines of linear-spaces optotypes in the following size equivalents: 20/200, 20/40 and 20/32. Near vision portion also contains 2 columns of proportionally spaced characters with size equivalents of 20/50 to 20/16.

Chart measures 9" x 14" (22.9 cm x 35.6 cm). Chart set includes overlay, flash cards and response panel.

..... **904034**



Crowded HOTV Test Set

Crowded optotypes in the following size equivalents: 20/50, 20/40, 20/30, and 20/20. Includes chart, response key, training cards, switch box overlay, and instructions. 9" x 14" (22.9 cm x 35.6 cm).

..... **904015**

Crowded HOTV Chart Only

Includes switch box overlay. 9" x 14" (22.9 cm x 35.6 cm).


..... **900707**



LED Insta-Line Mounting Bracket

Add portability to your LED Insta-Line with the Insta-Line Mounting Bracket. Adding this bracket to our deluxe universal stand allows increased mobility. Mounting bracket can only be used with the LED Insta-Line Quantum.

..... **914100**



Stand for Professional, ESV1500, ESV1200, ESV1018, LED Insta-Line, Models A+, A, or M Cabinets

Add mobility to your illuminated cabinet with this deluxe universal 5-leg stand. Available with or without casters. **Note:** Using stand with LED Insta-Line requires LED Insta-Line bracket 914100 found above.

A Casters **705000**

B Without casters **704000**



Hyperopia Glasses (Hard Lens)

+1.75. Recommended for children 10 years or older.

..... **900003**

+2.25. Recommended for children ages 6 to 10 years who have greater accommodation.

..... **900002**

+2.50. Recommended for children ages 6 to 10 years who have greater accommodation.

..... **900005**



Canvas Carrying Case

For Quantum Series, Model A+, Model A, LED Insta-Line, and Insta-Line PLUS™. Adjustable shoulder strap provides comfort when carrying an illuminated cabinet.

..... **904020**

Children's Vision Screening/Testing Fun Frames



Lens Type	Butterfly	Horse	Parrot	Tiger
Opaque Occluder	460800	460900	461000	461100
Frosted Occluder	460811	460911	461011	461111
Polarized	460400	462000	462100	460700
Anaglyph	461200	461300	461400	461500
+1.50 Hyperopia	516000	525900	525700	528900
+2.00 Hyperopia	491000	491100	491200	491300
+2.50 Hyperopia	500300	500400	500500	500611
+3.00 Hyperopia	514900	515100	515000	515200

Stereoacuity Screening and Testing



Combination LEA Symbols and Stereo Fly Acuity Test

Tests designed to rapidly test for amblyopia and strabismus using gross thru fine stereopsis (4800 to 20 sec of arc). Graded circle test from 400 sec now down to 20 sec with no monocular clues. New improved booklet has answer key on back cover and includes both adult and child polarized viewers.

..... 100050



Combination LEA Symbols and Stereo Butterfly Acuity Test

Tests designed to rapidly test for amblyopia and strabismus using gross thru fine stereopsis (2500 to 20 sec of arc). Graded

circle test from 400 sec now down to 20 sec with no monocular clues. New improved booklet has answer key on back cover and includes both adult and child polarized viewers.

..... 101050



Random Dot E Stereoacuity Test

Designed to rapidly test for amblyopia and strabismus in early and non-readers and non-verbal children and adults. Choice of dot size acuity (standard and low vision). Includes adult and new child-friendly pediatric viewers with improved test frames.

..... 101550



Random Dot 2 LEA Symbols Stereoacuity Test

Designed to rapidly test for amblyopia and strabismus in early and non-readers and non-verbal children and adults. Expanded random dot LEA test (500, 250, 125 sec of arc). Graded circle test now down to 12.5 sec with no monocular clues. New improved booklet has answer key on back cover and includes both adult and child polarized viewers.

..... 100750



Random Dot 3 LEA Symbols Stereoacuity Test

Designed to rapidly test for amblyopia and strabismus in early and non-readers and non-verbal children and adults. Expanded Random Dot 3 with LEA Symbols tests gross disparity at 900, 650, 500 and 400 seconds of arc. Test children at 400, 200, 100 seconds of arc. Graded circle test now down to 12.5 seconds with no monocular clues. New improved booklet has answer key on back cover and includes polarized viewers.

..... 100900



Stereo Vision Viewers

Polarized stereoacuity test viewers.

A Adult 103050

B Intermediate. 103175

C Children 103150



PASS Test™ 1

Preschool Assessment of Stereopsis with a Smile (PASS) - Stereoacuity test for basic screening of stereopsis in children and non-readers. Includes: 3 Cards (720, Demo, and Blank), 1 Pair Pediatric Polarized Viewers, and 1 Pair Intermediate Polarized Viewers.

..... 101100



PASS Test™ 2

Preschool Assessment of Stereopsis with a Smile (PASS) - Stereoacuity test for school assessment of stereopsis in children and non-readers. Includes: 4 Cards (480, 240, Demo, and Blank) with Carrying Case, 1 Pair Pediatric Polarized Viewers, and 1 Pair Intermediate Polarized Viewers.

This test is also part of the VIP Kit developed by the Vision in Preschoolers Study Group for screening the vision of preschoolers.

..... 101200



PASS Test™ 3

Preschool Assessment of Stereopsis with a Smile (PASS) - Stereoacuity test for quantitative assessment of stereopsis in children and non-readers. Includes: 6 Cards (480, 240, 120, 60, Demo, and Blank) with Carrying Case, 1 Pair Pediatric Polarized Viewers, and 1 Pair Intermediate Polarized Viewers.

..... 101300

Polarized Variable Vectographs



Vectograph Kit

Polarized Variable Vectographs

Polarized Variable Vectograph Vision Therapy System for Vision Training. Improve your patients' vision therapy outcomes using the new home reinforcement Vision Therapy System.

All Vectographs include therapy binder, doctor and patient instruction manuals, patient vision therapy recording form, and standard polarized viewers.



Vortex Polarized Variable Vectograph

Stereopsis depth on Vortex image of 6100 seconds of arc. Letters' orient to the Vortex image above, below, and even with the image at +/-600 seconds of arc. 16-diopter range for base-in (divergence) training and 24-diopter range for base-out (convergence) training for a total accommodation range of 40 diopters.

..... 105550



VFD Variable Fixation Disparity™ Polarized Variable Vectograph

Two-dimensional Fixation Disparity target without central Fusion Lock and without peripheral Gem image. 16-diopter range for Base-In (Divergence) Training and 24-diopter range for Base-Out (Convergence) Training with total accommodation range of 40-diopters.

..... 106590



Small Saccadic Polarized Variable Vectograph

Treatment Vectograph comes with five rows of saccadic Fixation Disparity Cross stimuli presented at a 20/63 acuity level to monitor visual, binocular and accommodative functions more reflective of reading demands at higher grade levels (smaller print books). Two-dimensional Fixation Disparity Target stimuli with Number Central Fusion Lock help relate binocular alignment to accommodation (clarity of the target) with a total of accommodation range of 40 diopters.

..... 107101



Gem Polarized Variable Vectograph

Stereopsis depth of 700 seconds of arc, 16-diopter range for base-in (divergence) training and 24-diopter range for base-out (convergence) training for a total accommodation range of 40 diopters.

..... 106055



NFL Gem PL-NFL™ Polarized Variable Vectograph

Stereopsis depth of 700 seconds of arc and Fixation Disparity target without central Fusion Lock. 16-diopter range for Base-In (Divergence) Training and 24-diopter range for Base-Out (Convergence) Training with total accommodation range of 40-diopters.

..... 106595



Chain Polarized Variable Vectograph Vision Therapy System

Thress-dimensional Chain image to provide an introduction to binocular experiences. 16-diopter range for Base-In (Divergence) Training and 24-diopter range for Base-Out (Convergence) Training with total accommodation range of 40-diopters.

..... 105600



Frog Polarized Variable Vectograph Vision Therapy System

Vectograph has stereopsis depth of 5300 seconds of arc. Frog image has stereopsis depth of 2500 seconds of arc. 16-diopter range for base-in (divergence) training and 24-diopter range for base-out (convergence) training for a total accommodation range of 40 diopters.

..... 105800



Gem PL Polarized Variable Vectograph

Stereopsis depth of 700 seconds of arc including fixation disparity target with central fusion lock. 16-diopter range for base-in (divergence) training and 24-diopter range for base-out (convergence) training for a total accommodation range of 40 diopters.

..... 106050

Large Saccadic Polarized Variable Vectograph

Treatment Vectograph contains four rows of saccadic Fixation Disparity Cross stimuli presented at a 20/100 acuity level, reflective of the visual demands of larger print books. Two-dimensional Fixation Disparity Target stimuli with Number Central Fusion Lock help relate binocular alignment to accommodation (clarity of the target). Total accommodation range is 40 diopters.

..... 107100

Disparity Polarized Targets



Far Polarized Fixation Disparity Target

Vectographic 2-dimensional Fixation Disparity Target for assessing Far Point of Fixation Disparity and Associated Vergence Measures at Distance. Can be administered in free space or behind a refractor and be presented at a variety of working distances. Includes Doctor Manual and Standard Polarized Viewers.


..... 106575



Near Polarized Fixation Disparity Target

Vectographic 2-dimensional Fixation Disparity Target for assessing Near Point of Fixation Disparity and Associated Vergence Measures at Near. Can be administered in free space or behind a refractor and be presented at a variety of working distances. Includes Doctor Manual and Standard Polarized Viewers.

..... 106550



Randomized Fixation Disparity Polarized Target

Assessment tool for measuring binocular/ accommodative function during a saccadic eye movement task using a sequence of saccadic identification targets. Vectographic silvered card with four rows of saccadic Fixation Disparity Cross stimuli presented with a 20/100 target size, reflective of the visual demands of larger print books. Two-dimensional Fixation Disparity Cross stimuli with Number Central Fusion Lock help relate binocular alignment to accommodation (clarity of the target)

Large Targets 107103

Small Targets..... 107104



Fixation Disparity Polarized Target

Fixation Disparity Polarized Target Vectographic 2-dimensional card with 2 Fixation Disparity Crosses positioned vertical to one another to facilitate the comparison of binocular performance when using a lower versus higher spatial frequency or acuity demand. Top Fixation Disparity Cross at 20/95 acuity level mimics reading demands of larger print books.

The Small Target bottom Fixation Disparity Cross at 20/63 acuity level mimics reading demands at higher grade levels (smaller print books) while the Large Target bottom cross is at 20/100 acuity level.

Fixation Disparity Polarized Target Vectograph includes manual and standard polarized viewer.

Large Targets 107102

Small Targets..... 107105

Kits



PL Binocular Vision Dysfunction Diagnostic and Treatment System™

Quickly and accurately diagnoses large to small fixation disparities with this all inclusive kit, which contains the following items:

- 106575 - Far Fixation Disparity Target™
- 106550 - Near Fixation Disparity Target™
- 106055 - Gem™ Polarized Variable Vectograph
- 106050 - Gem PL™ Polarized Variable Vectograph
- 106595 - GemPL NFL™ Polarized Variable Vectograph

Therapy Binder with Pen, Instruction Manual, and Standard Polarized Viewers.

..... 107000



Binocular Accommodative Saccadic Series

BASS-PL (Binocular Accommodative Saccadic Series) Vision Therapy System for improving patient visual performance and reading efficiency, consisting of diagnostic Fixation Disparity Targets and treatment Vectographs with saccadic identification targets at two acuity levels.

- Fixation Disparity Polarized Targets (2 cross targets) (Large and Small)
- Randomized Fixation Disparity Polarized Targets (Large and Small)
- Saccadic Polarized Variable Vectographs (Large and Small)

..... 107500



Vectograph Viewers

Greatly simplifies vision therapy by holding vectographs at the proper angle and allows sliding the vectographic panels to the "breaking point" and back to the "recovery point". These vectograph viewers support the full extension of the panels to their maximum base-in and base-out capabilities. A white background behind the vectograph supports base-out vergence work (convergence). The white background can be removed to reveal the clear background. This permits use of the vectographs for fusional divergence work, which requires that the patient view beyond the plane of the vectograph.

Available as either a single or double Vectograph Viewing stand.

- A** Single. **585600**
- B** Double. **587800**



Illuminated Vectograph Viewer

The Illuminated Vectograph Viewer is dual height to permit viewing of two VA vectographs at the same time. Provides even illumination and supports the full extension of vectographs.

The Vectograph Viewer is specifically designed to accommodate vectographs that measure 4.5" H x 6.5" W (11.5 cm x 16.5 cm) allowing for 2 rows top and bottom as shown. Sits at the proper viewing angle and includes a handle on top for easy portability.

The Illuminated Vectograph Viewer measures 10" x 16" x 2.5" (26 cm x 41 cm x 7 cm). Uses LED illumination with universal power supply (50-60 Hz, 110/240 AC) ensures long life.

. **621700**



Lichtenstein Fixation Box

The Lichtenstein Fixation Box gives you the added ability to test fusion, as well as to maintain fixation for your younger patients. The box contains a Worth 4-Dot examination, a white fixation light for adults and older children, and a fixation light that changes colors for strabismus measurements or refraction. LEDs provide increased bulb life. Universal power supply. 7.25" H x 6" W x 1.5" D (18.4 cm H x 15.3 cm W x 3.8 cm D). CE and CSA Listed.

. **995000**



Worth 4-Dot Test Attachment

The Worth 4-Dot Attachment fits on your Welch Allyn or Keeler transilluminator and supports Worth 4-Dot testing using red/green goggles. Compared to conventional flashlight versions, the Worth 4-Dot Test Attachment does not require cumbersome flashlight units or D batteries. Use with red/green anaglyph adult glasses 955100 or one of our Children's Fun Frame anaglyph glasses found on page 56.

Standard Worth 4-Dot **475900**

LEA Symbols Worth 4-Dot. **499600**



Combination Polarized and Anaglyph Flipper

Combination polarized and red/green stereo lenses into one flipper. Use the red/green for the Worth 4-Dot or other anaglyph test; then, flip to polarized lenses for stereoacuity test.

. **488400**



Berens 3 Character Test

The Berens 3 Character Test is designed for patients who are unable to respond to the Worth 4-Dot Test, such as children who cannot name colors and individuals with low color perception. The 3 characters are: a red girl, a green elephant and a white ball. All are reproduced on a special insert that goes over the flashlight lens. The test set consists of a high-quality, three-cell flashlight with a pair of universal spectacles together and complete instructions. 3 D-cell batteries included.

. **153090**



C Worth 4-Dot Flashlight with Anaglyph Glasses.

Batteries not included.

Standard **950100**

LEA Symbols. **520600**

D Reversible Red/Green Anaglyph Glasses

Durable metal frame glasses with spring-loaded temples. Storage case included.

. **955100**



Red & Green Wraparound Glasses

Wraparound anaglyphic (red/green) glasses offer a functional advantage for vision testing and therapy; they prevent side light regardless of where the patient moves his or her eyes or head.

E Adult **153160**

F Children **153170**

Examination Tools




Vergence Facility Prism

Quick screen for binocular vision. Use to get the patient to cycle between the 2 prisms up to 15 cycles per minute. Use at 16" (40 cm) with 20/30 near point line.

Prisms are 12D base-out and 3D base-in and 1.5 inch high (38 mm). Comes either mounted on a stainless stick holder or loose.

- Mounted (Shown) **477700**
- Prism Only **477701**



Stainless Steel Stick Prism Set

Prisms are mounted on a stainless steel stick for easier and more sanitary presentation. Set includes: Red lens plus diopters of 0.5, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25, 30, 35, 40, 45, and 50 in a wooden storage box. Prisms are

- 37 mm square.
- 16-piece set **467100**



Loose Prism Sets

These high quality prisms are used to measure the deviation in strabismus. It is important to conduct polarized testing through prisms that are known to cause the least interference with polarization. Prisms are available loose or with a magnetic attachable handle.

Set of 16 includes: 0.5, 1, 2, 4, 6, 8, 10, 15, 20, 25, 30, 35, 40, 45, 50 diopter prisms, and red lens.

Set of 22 includes: 0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50 diopter prisms, and red lens.

Each measures 37 mm square and includes case.

Loose Prism Sets

- 16-piece set **444900**
- 22-piece set **444500**

Loose Prism Sets with Magnetic Handle

- 16-piece set **444910**
- 22-piece set **444510**



Prism Bars

Made from durable, lightweight acrylic. 16 segments that provide 1, 2, 4, 6, 8, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, and 45 PD. Vertical prism bars have 15 segments that provide 1, 2, 3, 4, 5, 6, 8, 10, 12, 14, 16, 18, 20, 25, and 30 PD. Includes protective case.

- A Horizontal **444600**
- B Vertical **444700**
- A and B Combined **444800**



Retinoscopy Paddles with Axis Compass

This two-paddle version has one paddle with plus lenses and the other for minus. Optics on both are 0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 8.0, and 12. Includes 180° axis compass printed directly along the top of each paddle. Lenses are 16 mm diameter, recessed to help prevent scratching. Lenses are glass for greater accuracy and durability. Includes protective case.

- **506600**



Risley Rotary Prism

Large aperture (1 1/23 mm) rotating prisms permit more flexibility compared to your Phoropter. Two versions are offered: +/-20 PD or +/-30 PD with handles.

The prisms are controlled with a convenient knurled rotating knob. Accurately calibrated dial shows prism diopter selected and rotates through 30 PD right-exo to 30 PD left-eso.

- 20 PD **467300**
- 30 PD **467400**



Red Maddox Insert

- Red maddox insert for Risley Rotary Prism. **207002**



Prism Flipper Set

Five Prism Flipper Set includes 2, 3, 4, 6 and 8 diopters with base in and base out plus one occluding Flipper and Stand.

Lens Types are:

- 2 Diopter (2^+2^)= 4^ total
- 3 Diopter (3^+3^)= 6^ total
- 4 Diopter (4^+4^)= 8^ total
- 6 Diopter (6^+6^)= 12^ total
- 8 Diopter (8^+8^)= 16^ total

- **537900**



Individual Prism Flipper Lens

Handheld individual prism flippers. Orientation of the prism flippers are Base In and Out.

Individual Prism Flipper Lens

- 1 Diopter (1^+1^)= 2^ total **137001**
- 2 Diopter (2^+2^)= 4^ total **137002**
- 3 Diopter (3^+3^)= 6^ total **137003**
- 4 Diopter (4^+4^)= 8^ total **137004**
- 6 Diopter (6^+6^)= 12^ total **137006**
- 8 Diopter (8^+8^)= 16^ total **137008**



38mm & 48mm Corrected Curve Prisms

Corrected Curve Prisms are used with our Yoked Prism Frames. The lenses are available individually in base powers from 0.25D to 12D. Each lens is clearly marked with the diopter power and a 'line' indicating the base. The 38 mm are recommended for use with Pediatric Patents.

Visit www.good-lite.com for more



Yoked Prisms Frames

Good-Lite offers a variety of Yoked Prism Frames for your practice. 38 mm for Pediatrics, 45 mm for VTE (Italy) lenses and 48 mm for adults with different options.

Visit www.good-lite.com for more



Color-Coded Flipper Set

Flipper Set with 6 flippers, plus and minus each of 0.25, 0.50, 0.75, 1.00, 1.50, and 2.00 diopters. Each is marked with diopter power and handles are color-coded. The lenses are durable and will last for years. Includes attractive holder.

- 6 Color-Coded Flipper Set** (shown) **499700**
- 12 Color-Coded Flipper Set** **517900**



Individual Flipper Lens

Individual flippers available in various lens configurations. Durable with a soft handle.

- Individual Flipper Lens**
- 127025** +/- 0.25 Lens
 - 127050** +/- 0.50 Lens
 - 127075** +/- 0.75 Lens
 - 127100** +/- 1.00 Lens
 - 127125** +/- 1.25 Lens
 - 127150** +/- 1.50 Lens
 - 127175** +/- 1.75 Lens
 - 127200** +/- 2.00 Lens
 - 127225** +/- 2.25 Lens
 - 127250** +/- 2.50 Lens
 - 127275** +/- 2.75 Lens
 - 127300** +/- 3.00 Lens
 - 127400** +/- 4.00 Lens



Flipper in Canvas Pouch

Set of 6 Flippers in Canvas Pouch with elastic closure. Black. Includes 6 Flippers: with plus on one side and minus on the other of 0.25, 0.50, 0.75, 1.00, 1.50, and 2.00 diopters.

- **552500**

Thorington Card

The Thorington Card quickly measures heterophoria. The near card is used with a Maddox Rod. Scoring reveals if the patient has normal alignment when fusion is broken or if the patient has abnormal alignment.

The card helps diagnose Uncrossed or Crossed Diplopia and RhyppoP, LhyperP, RhyperP, and LhypoP. Includes a 1 mm hole in center for the fixation light. With the near adult card, each optotype represents 1 prism diopter, and measures +/-12 diopters. The card measures 5.75 x 4.75 inches (147 x 121 mm) and is available in 2 versions: Penlight/Transilluminator attachment or Touch LED Light.



Adult Thorington Card

Uses a combination of both numbers and letters. Available in 2 versions.

- Transilluminator/Penlight **489300**
- Touch LED. **555800**



Pediatric Thorington Card

- Available in 2 versions.
- Transilluminator/Penlight **539000**
 - Touch LED. **555900**

Thorington Options



Transilluminator



Penlight



LED Touch Light



Pocket Neutral Density Bar

Used to measure Afferent Pupillary Defects (APD) which is defined as diminished pupil reaction to light levels. It uses 6 standard ND levels for quantifying APD. The log units quantify amount of light transmitted through the filter. Includes 6 filters rated as 0.3, 0.6, 0.9, 1.2, 1.5, and 1.8 log unit optical density values. The 0.3 ND filter is most useful for uncovering a borderline pupil defect.

The filter bar can be used in conjunction with the "Swinging Flashlight" test to detect and quantify the depth of any afferent conduction defect.

Fits in shirt pocket. Includes leather slip cover. Overall bar is 6.5" x 1" (16.5 cm long and 3 cm wide). Each filter measures 7/8" by 1" (2.5 cm x 3 cm).

- **470550**



Neutral Density Filter Bar

Used to measure afferent pupillary defects (APD) also called "diminished pupil reaction to light levels". Provides the filter values: 0.3, 0.6, 0.9, 1.2, 1.5, and 1.8 density log units. Each filter is 1.5" (38 mm) square with an overall bar size of 10" by 1.5" (25.3 cm by 3.8 cm). Comes with a handle for easy use and case for storage.

- **470200**



Neutral Density Bar Slide

- Will add 1.8 density to filter bar.
- **488800**



Fixed PD Trial Frames

These axis adjustable Fixed PD Trial Frames are comfortable because they are lightweight and sized to the patient's needs and adjustable for temple length. Available in sets or individually. The set includes a PD rule to make accurate frame selection quick and easy. Each frame holds 3 trial lenses in front and a 4th lens behind on each side. All have a knob to rotate the axis of the front lenses.

The light blue fixed trial frames come in smaller PD sizes of 52, 54, 56, 58, and 60 mm. The 52, 54, and 56 mm PD trial frames will hold 3 lenses while the 58 and 60 mm light blue frames will both hold 4.

The light gray fixed trial frames come in larger PD sizes of 62, 64, 66, 68, and 70 mm and each will hold 4 lenses.

Light Blue Frames

- 52 mm PD. **530600**
- 54 mm PD. **530700**
- 56 mm PD. **530800**
- 58 mm PD. **530900**
- 60 mm PD. **531000**

Gray Frames

- 62 mm PD. **531100**
- 64 mm PD. **531200**
- 66 mm PD. **531300**
- 68 mm PD. **531400**
- 70 mm PD. **531500**

Adult Axis Adjustable Trial Frame Set

Includes sizes of 62, 64, 66, 68, and 70 mm, a storage box that holds 5 frames and a PD rule to make accurate frame selection quick and easy.

..... **532200**

Pediatric Axis Adjustable Trial Frame Set

Includes the sizes of 52, 54, 56, 58, and 60 mm, a storage box that holds 5 frames, and a PD rule to make accurate frame selection quick and easy.

..... **531700**



Adjustable Adult Trial Frames

Easy-to-use and lightweight adjustable adult trial frames. Lens holders on right/left side available for inserting 5 standard trial lenses simultaneously.

- Range of PD scale: 50~80 mm
- Single side (right or left).
- Scale: 25~40 mm;
- Minimum scale Value: 1 mm
- Axis scale step: 5°
- Range of nose bridge adjustment:
- Height: 0~23 mm
- Angle: 0°~360°

These light-weight frames weight only 61g.
..... **455500**



Universal Adult Frames

The following adjustments can be made: temple angle, temple length 4.125" to 5.5" (98 mm to 135 mm), nose bridge height 0 mm to 14 mm, nose angle 0° to 30°, PD adjustable to 24 mm to 40 mm on each side (48 mm to 80 mm total both eyes).

- Axial calibration on each lens:
 - left - 120° to 0° to 135°
 - right - 45° to 180° to 60°
 - Lens rotation is full 360°

Accommodates three 38 mm lenses in front and 1 behind. Clear, precisely printed scales. Foam-covered nose bridge for comfort.

..... **455400**



Raff Level Trial Flipper

Trial lens flipper with built-in level helps speed the lens trial and prescribing process creating a more accurate presentation while eliminating inadvertent off-axis lens rotation. The higher accuracy provides a more positive patient experience and the higher-tech appearance helps build patient confidence.

..... **151560**



Pediatric Trial Frames

Pediatric spring-hinged trial frames are both lightweight and easy to use. PD can be adjusted from 54 mm to 70 mm and accommodates up to four 38 mm trial lenses on each side. Also contains axis of astigmatism and one touch adjustment on nose pad.

- Blue. **469800**
- Charcoal **469900**



Clip-On Trial Frames

Attaches to existing Rx glasses for a quick opportunity to over refract. Use of the clips reduces vertex distance error compared to phoropter. Once the patient's distance Rx is correct, clips can help check presbyopic addition. Axis markings provided for use with cylinders. Set of 2 Clip-On Frames

..... **301100**



8-Well Trial Lens Holder

Easily use your existing trial lenses with this handheld 8-Well Trial Lens Holder. Contains 2 lens holders, 2 lens holders with axis scale, and nose bridge with a handle. Made of sturdy black plastic for long lasting durability.

..... **532300**



10-Well Trial Lens Holder

Easily use your existing trial lenses with this handheld 10-Well Trial Lens Holder. Contains 5 wells with a handle. Made of sturdy black plastic for durability.


..... **151460**

Bagolini Striated Lenses

Finely striated, clear lenses used for evaluating retinal correspondence in patients with eye deviations. Consists of 2 circular lenses, indexed and mounted in 38 mm trial lens metal rings with tabs for use by the examiner. In testing suppression, the lenses are placed in a trial frame with the striations at an angle of 135° and 45° for each eye respectively.



..... **346700**



Maddox Lenses for Trial Frames

Maddox lenses mount in trial rings oriented to the handle. Available as either Red, Green, or White.

- A** Red **473000**
- B** Green..... **473100**
- C** White..... **473200**



Trial Lens Set

These are packaged in an attractive lockable briefcase and rims are plated for quick identification. The Trial Lenses kit colored rims with Red = Plus and Black = Minus or chrome-rimmed lenses. Also includes case, Prisms and accessory lenses. Space is provided for your choice of trial frame.

- Reversed Red/Black Rims **592000**
- Chrome Rims **449300**



Occluders for Trial Frames

Various occluders mounted in trial frames.

- D** Opaque Occluder **532400**
- E** 1 mm Pinhole Occluder **532500**
- F** Frosted..... **533100**



Trial Lenses for Your Practice

Good-Lite has a variety of individual trial lenses for your Practice available in both Cylinder and Sphere. Individual lenses have a chrome rim and come in standard sizes. Contact us with the diopter increment that you need for your set.

Cylinder Trial Lens Range

- ±20.0D to ±8.0D — 1.0 diopter increment
- ±7.50D to ±5.50D — 0.5 diopter increment
- ±5.25D to ±0.12D — 0.25 diopter increment

Sphere Trial Lens Range

- ±20.0D to ±8.0D — 1.0 diopter increment
- ±7.50D to ±5.50D — 0.5 diopter increment
- ±5.25D to ±0.12D — 0.25 diopter increment

Note: ±0.37D and ±0.12D are available for both Cylinder and Sphere Trial Lenses.

- 128000** Cylinder - Specify Diopter Increment
- 129000** Sphere - Specify Diopter Increment

Occluders



Multiple Pinhole Spectacle Occluder

- Pinhole Spectacle **700010**
- Spectacle Only **700012**



Multiple Pinhole Glasses

Each lens contains 49 pinholes in circular formation. Lightweight and easy to use.

..... **204007**



Reversible Pinhole Occluding Glasses

One side is occluded with 5 pinholes, the other is fully occluded. Sturdy metal frames with straight temples for easy reversibility.

..... **204006**



Reversible Occluding Glasses

One side is occluded, the other is open. Sturdy metal frames with straight temples for easy reversibility.

..... **452300**




Mini Handheld Occluder for Trial Frames

These easy-to-use and portable Mini Occluders are only 4 inches (10.1 cm) long. Available in 4 different versions; full occlusion, frosted, single pinhole or 17 pinhole occluder. Handle in printed with half moon to the edge with 8 sizes: 2 mm to 9 mm and 6 cm ruler and **will fit in Trial Frames.**

- Opaque **605700**
- Frosted **618600**
- Single Pinhole **605800**
- 17 Pinholes **605900**
- Red Maddox **606000**
- White Maddox..... **607400**
- Green Maddox **607300**

Examination Tools



Children's Vision Screening/Testing Fun Frames

Who said serious tests had to be done with boring test glasses? Add fun to your pediatric tests with our animal motif frames. Use these fun frames as occluders to increase testability, as pinhole occluders or hyperopia glasses, and during tests requiring anaglyph or polarized glasses,

such as stereoacuity tests. Occluders come in left and right pairs with the same frame style. These covered frames are more colorful than shown here and they are surprisingly durable! Not for possession by children below age 3.

Lens Type	Butterfly	Horse	Parrot	Tiger
Opaque Occluder	460800	460900	461000	461100
Frosted Occluder	460811	460911	461011	461111
Polarized	460400	462000	462100	460700
Anaglyph	461200	461300	461400	461500
+1.50 Hyperopia	516000	525900	525700	528900
+2.00 Hyperopia	491000	491100	491200	491300
+2.50 Hyperopia	500300	500400	500500	500611
+3.00 Hyperopia	514900	515100	515000	515200



Sunflower Occluder Glasses

These fun sturdy occluding glasses have flip-up sunflowers to permit occlusion of each eye. Supports 55 mm to 65 mm PD and recommended for children 3 years to 5 years old.

..... 548500



Happy Flip-up Occluders

Fun colorful Happy Froggy or Pink Blossom pediatric occluding glasses have flip-up lenses to permit occlusion of each eye. Supports 55 mm to 65 mm PD and recommended for children 3 years to 5 years old.

Pink Blossom 587100
Happy Froggy 587200



LEA Face Disposable Occluder

These easy-to-use disposable occluders are a perfect addition to any vision testing or screening program. Printed on sturdy card stock with LEA Smiling Face for fixation. Each package contains 50 occluders glued along the edge for easy storage.

..... 253200



Adhesive Eye Occluder

Lightweight, breathable elastic cloth for comfort and adhesive for a secure fit and easy removal.

Individually packed for added protection. Dimensions: 1.875" x 2.5". Recommended by eye care professionals.

..... 760001



Opaque Plastic Occluders

Set of 6 in either black or gray.

Black 751000
Gray 750000

Available individually.

Black 700545
Gray 700544



Clip-On Occluder

Lightweight black plastic occluder clips onto any existing frame using metal spring clip.

..... 121002



Multiple Pinhole Clip-On Occluder

Lightweight black metal multiple pinhole occluder securely clips onto any existing frame.

..... 121100



Frosted Clip-On Occluder

Lightweight frosted occluder securely clips onto any existing frame.

..... 121001

Opaque and Frosted Occluders

approximately 9.5" (24 cm) in length.



- A** Single-ended black occluder.
..... **204004**
- B** Single-ended single pinhole occluder.
..... **479500**
- C** Single pinhole and black occluder.
..... **120101**
- D** 17-pinhole and black occluder.
..... **481800**
- E** 49-pinhole and black occluder.
..... **120302**
- F** Frosted and black occluder.
..... **681000**
- G** Frosted "See the Eye" occluder.
..... **760000**
- H** White opaque occluder.
..... **204005**



Multiple Pinhole Spectacle Occluder

- Pinhole Spectacle **700010**
- Spectacle Only **700012**



Mirrored Opaque Occluder


- Long handle single ended occluder with acrylic mirror.
..... **618700**

Maddox Lenses



- I** Red Maddox with 45° lines from center.
..... **205001**
- J** Green Maddox with 45° lines from center.
..... **205000**
- K** Clear Maddox with 45° lines from center.
..... **205002**

Combination Maddox Lenses




- L** Red Maddox with 45° and vertical lines from center.
..... **211000**
- M** Red Maddox with 45° lines and red filter.
..... **209000**
- N** Red Maddox with 45° lines and occluder.
..... **207000**

Combination Maddox Lenses and Pinhole Occluder



- O** Red Maddox with 45° lines and single pinhole occluder.
..... **543500**
- P** Clear Maddox with 45° lines and single pinhole occluder.
..... **623800**
- Q** Green Maddox with 45° lines and single pinhole occluder.
..... **623900**

Frosted Occluder and Maddox



- Combination frosted occluder and Maddox with 45° lines from center.
- R** Clear **207001**
- S** Red **583100**
- T** Green **595300**

Filter Lenses



- U** Red Filter lens.
..... **206001**
- V** Red Filter lens and opaque occluder.
..... **208000**

Examination Tools



Leaf Style Clip-on Occluders

Available in 44 mm, 46 mm, 48 mm or 50 mm in height; Black or Beige color.

Set of 6 occluders

Beige		Black	
481144	44 mm	481044	44 mm
481146	46 mm	481046	46 mm
481148	48 mm	481048	48 mm
481150	50 mm	481050	50 mm



Round Style Clip-on Occluders

Available in 44 mm, 46 mm, 48 mm or 50 mm in height; Black or Beige color.

Set of 6 occluders

Beige		Black	
481244	44 mm	480144	44 mm
481246	46 mm	480146	46 mm
481248	48 mm	480148	48 mm
481250	50 mm	480150	50 mm



Medical Eye Shield

This plastic shield provides temporary protection for the eye. Comes with elastic string in packages of 24 shields. Available in either Adult or Pediatric sizes.

Adult Eye Shield – 77 mm x 53 mm

Black	131000
Beige	133000

Pediatric Eye Shield – 63 mm x 44 mm

Black	571600
Beige	571700



Foam-Edged Clear Moisture Chamber

Foam-Edged Clear Moisture Chamber with clear lens permits post surgical observation of the condition of the eye. Also permits the wearer to see through the patch which makes this item appropriate for Bells palsy patients and for use as a moisture chamber. For Bells palsy patients, this patch prevents air flow on the eye. Uses non-latex foam, fits over either eye with wide elastic band ensuring a secure fit.

Available in two sizes; Large chambers are 2.8125" x 1" (7.14 cm x 5 cm) and Small chambers are 2.375" x 1.6875" (6 cm x 4.5 cm).

Large Size: 2.8125" x 1" (7.14 cm x 5 cm)	
Each	675801
Set of 2	675802
Small Size: 2.375" x 1.6875" (6 cm x 4.5 cm)	
Each	675901
Set of 2	675902



Visual Field Fatigue Diffuser

Designed by John Karickhoff, MD, to block vision by diffusion rather than occlusion. The diffusion level was selected to block viewing of the largest targets and fixation spot of all the standard parameters but allows 97% of the light of the bowl to enter the covered eye while keeping air away from the eye. This eliminates the continuous ocular rivalry caused by standard opaque occluders. Clear Visual Field Fatigue Diffuser with white elastic straps to fit all sizes.

Sold in a package of 3, Size: 3" x 2" (75 mm x 50 mm)
 **166000**



Moisture Chamber Goggles

Large Dry Eye Moisture Chamber Goggles provide a chamber to ease the pain for dry eye sufferers and those with Bell's palsy providing protection for both eyes. Dry Eye Moisture Chamber Goggles will protect eyes from air movement, dust, smoke, sprays and other particulates. Also, the Dry Eye Moisture Chamber Goggles retain moisture helping to keep the eye lubricated. Dry Eye Moisture Chamber Goggles come with a comfortable wide elastic to avoid "pinching" and can be worn under glasses.

The durable clear material is made of PETG and is optically "clean". Dry Eye Moisture Chamber Goggles is latex free and can be easily cleaned with a damp cloth.

Large Size: 2.8125" x 1.9375" (7.2 cm x 5 cm)	
Each	681700
Set of 6	681706
Small Size: 2.375" x 1.6875" (6 cm x 4.5 cm)	
Each	681800
Set of 6	681806



Foam-Edged Patches

Foam edging ensures comfort and excellent occlusion. Durable; can be used repeatedly with careful cleaning. Adult patches are 2.8125" x 1" (7.2 cm x 5 cm). Children's patches are 2.375" x 1.6875" (6 cm x 4.5 cm). White Children's Eye Patch with Elastic

A	487900
Black Adult Eye Patch with Elastic	
B	487700
White Adult Eye Patch with Elastic	
C	487800

Note: All foam-edged patches contain latex.

Fixation Products




Small and Large Fixation Paddles

The small 2" (5 cm) diameter "smiling face" fixation paddle was developed based on studies by Robert Fanz. Features include grating on reverse side for preferential looking, small pictures of a bunny and a giraffe on end of stick for fixation, and marks for measuring pupil size and corneal diameter measurement on the handle.

The medium 5" (13 cm) diameter fixation paddle is for assessing children with severe disabilities to determine the distance at which children respond to a high-contrast face figure.

- A** Medium Target **253100**
- B** Small Target **253000**



Toucan Fixation PD Rule

Toucan fixation and PD rule with near optotypes, including LEA Symbols, LEA Numbers, and HOTV optotypes for 20/32, 20/63, 20/80, and 20/100. Includes inch/mm rules and diopter scale plus pupilometer in half mm increments. 6" (15.2 cm) long.

- **497000**



Flashing Fixation Stars or Happy Face

For mounting on a retinoscope or ophthalmoscope; includes 3 Flashing Fixation Stars or 2 Flashing Happy Faces. Both have 9 mm magnetic disks in box and double-sided tape for mounting. Fits most retinoscopes.

- Flashing Stars **489000**
- Flashing Happy Faces **533900**



FixiStix™ Wiggle Pictures

Wiggle Pictures are back, but with a brand new style! This classic product features new designs recreated especially for the pediatric market. Set includes 3 plastic spatula sticks with wiggle pictures on the front and back.

- **935800**



Blue Hippo Triple Fixation Set

Set of 3 double-sided fixation sticks with 12 colorful images that will not rub off. Images include Birthday Cake, Toucan, Blue Hippo, Sun, Sunflower, Tree, Snowflake, Lion, Present, Maracas, Dog, and Cat. Size: 5.75" x 1.325" (14.6 cm x 3.4 cm).

- **543200**



Detailed Fixation Sticks

For detailed accommodation stimulation targets. For use alone or together. Precisely calibrated. For saccade testing. Durable printing.

- **542055**



Mini Spinning Fixation Globes

The Mini-Spinner is only 4" (10 cm) tall with a 1.5" (37 mm) globe where the blue, red and yellow LEDs spin. Light pattern changes as ball spins. The

button must be held down for the lights to light and the ball to spin. Runs for a long time on three LR44 hearing aid batteries. Comes 2 to a package and includes instructions for battery replacement.

Note: Colors may vary.

- **591500**



Pediatric Fixation Sticks

Pediatric Fixation Sticks provide a series of detailed accommodation stimulation targets for the young patient. The targets are calibrated for near. Included are (as examples) a Clown Face, a Dog, a Beach Ball, etc. Each target is in brilliant color for patient ease.

The set consists of 3 durable plastic sticks with targets on both sides.

Single sticks can be used for performing various tests: unilateral cover/uncover test at near, alternating cover test at near, pursuit testing, near point of convergence (NPC), motilities (ductions and versions), vergence facility at near, fixation stability at near in diagnostic positions of gaze, amplitude of accommodation, and accommodative facility. Dual sticks with different optotypes exposed to the patient can be used for saccade testing.

- **542075**



Face Fixation Sticks

Provides a series of detailed accommodation stimulation targets for the young patient. The set consists of 3 durable plastic sticks with fixation faces smiling on one side and frowning on the other side. Single sticks can be used for performing various tests: unilateral cover/uncover test at near, alternating cover test at near, pursuit testing, near point of convergence (NPC), motilities (ductions and versions), vergence facility at near, fixation stability at near in diagnostic positions of gaze, amplitude of accommodation, and accommodative facility.

- **542065**

Examination Tools




DXDT's Quantitative Retina Test Grid at-Home Retina Health Monitoring System

The system is used to monitor for changes in retina health for patients under medication that may affect retina health and individuals who want to get the earliest possible warning if retina issues emerge. Free analysis tools will estimate the rate of improvement or degradation in retina health or to test for the effectiveness of treatment interventions.

The system, preloaded on a USB drive, operates on any Microsoft Windows-based personal computer. The system provides grids of selectable sizes and color schemes.

..... 300026



Amsler Grid Book

This is an expanded version of our popular Amsler Grid Book and provides 7 Amsler Grids and a score sheet master. In addition to the 4 grids we've always included, this expanded booklet adds the White Dots, Parallel Line, and Parallel Line with Central Line Grids. In addition to the Standard Grid, the Red on Black is used to diagnose a color scotoma. The White Dots Grid reveals a scotoma without distortion. The two Parallel Line Grids are used with patients needing a more simplified Amsler Grid and are presented both horizontally and vertically. Each of the small squares on the Grid with Small Squares approximates 0.5 degree and is used for patients with small, central scotomata.

..... 730026



Damato Campimeter

This handheld campimeter enables practitioners to assess the central 30-degree visual field when other methods are impractical. It consists of a double-sided card with a test grid on each side. Each test grid displays a central stimulus surrounded by a series of numbered fixation targets. The patient looks at each numbered fixation target in turn, so that the stimulus moves through 30 relevant areas of the visual field. While the patient gazes at a selected target, the examiner dials a disc so that a black-on-white stimulus briefly appears in a central window on the test card. The normal blind spot is used as a control. An occluder for the eye is attached to the side of the test card by a rigid arm to ensure that the patient is at the correct working distance. This test is fully portable and does not require a power supply. The "perimetry by numbers" technique makes the examination relatively easy to learn and administer. Includes recording forms.

..... 293000



Amsler Grid Viewer

Provides an amsler grid in a simple key chain viewer ensuring monocular viewing. Priced low enough to give to patients for home use. Available in assorted colors and color shipped can not be guaranteed. Purchase either individually or in a package of 5 viewers.

Each 563300

5-Pack 563400

The Automated Confrontation Visual Field Tester (ACV)

The ACV is both cost-effective and easy-to-use where both eyes can be tested in less than a minute. It is also more accurate and sensitive than using finger confrontation methods.

The ACV tester consists of a black unit measuring 12 cm by 17 cm. Four 1.0 mm red LEDs are positioned 11.25 cm apart, which results in an angular separation of 15° when the device is 60 cm from the patient.

The tester is powered by 2 standard AA batteries (included).


..... 476550



LEA Flicker Wand

The LEA Flicker Wand is 26.5" (67.3 cm) long with a 20" (50.8 cm) flexible wand with a diode at the end that can be used as either a flickering or non-flickering stimulus at 4, 40, and 400 cd/m². The curved wand is moved from behind the person forward in all 4 quadrants of the visual field to detect losses in the periphery of the visual field.

..... 280000



Tangent Screen

Tangent Screens provide a way to chart the central and Para-central field of vision for diagnosing pathological conditions. Testing can be conducted at 3' to 6' (1 m to 2 m). Markings are stitched in semi-visible black thread showing concentric circles, radial lines, and blind spots. Markings ensure ease of testing, scoring, and transcribing the mapped fields for visual field analysis. A white fixation button is attached. Special orders can be filled for white or red thread. Tangent Screens are available using a variety of mountings.

1 m 15 deg **931000**



Tangent Screen Wand

Test object is inserted into the fabric leaving only the object visible. Half-turn rotation obscures the test object behind the wand. Alternating the rotation, objects can be made to disappear and reappear for confirmation response.

18" (45.72 cm) Wand. **920000**

39" (1 m) Wand **922000**



Tangent Screen Test Objects

Flat test objects of various colors and sizes are mounted on non-reflective

washable vinyl with a black felt background and pin shank for use with tangent screens. **918000**



Optokinetic Drum

Tests for nystagmus when rotated at 8-10 rpm. Drums are available in three sizes: 10", 8" or 4" high by 6" diameter (25.4, 20.32, or 10.16 cm x 15.2 cm). Lightweight drum of professional quality with ball bearings for effortless spin mechanism. Standard black and white striped drum for use 24" to 30" (61 cm to 76 cm) from patient.


10" Drum. **547500**
 8" Drum. **576300**
 4" Drum. **616700**



Optokinetic Flag

This flag is used to produce optokinetic nystagmus (OKN) physiologic jerk in a patient's eyes. This is elicited when the patient attempts to fixate on moving, alternative, contrasting squares that are moved in the field of vision. Includes instructions. 37.25" (94.62 cm) long.

Red & White with 1.5" segments . . . **547400**
 Black & White with .75" segments . . **623400**
 Black & White with 1.5" segments . . **607700**
 Black & White with 3.0" segments . . **623500**



Optokinetic Banner


This Optokinetic (OKN) Banner is large enough to fill a patient's visual field. OKN Banner is used to produce involuntary Optokinetic Nystagmus (OKN) physiologic jerk in a patient's eyes. The Red and White stripes on the OKN side are 1.5" (4 cm) wide. Dowels are sewn into each end to provide vertical stability and ease of use. Measures 27" x 10" (68 cm x 25 cm) high.

Red and White. **588600**
 Black and White **607800**

Wolff Fixation Wands

Wolff Wands provide specialized fixation targets used in various procedures such as retinoscopy, Cover Test, and Ocular Motility. Various techniques employ the wands to evaluate accommodative lag, anisometropia, near point convergence, tracking, peripheral awareness and saccadic eye movements.

The rods are covered in comfortable black plastic coating material for durability.



Acrylic Wolff Wands

Available with either a 1" (25 mm) or 1.5" (38 mm) polished acrylic ball on the end.

Small. **657700**
 Large. **657800**



Brass Ended Wolff Wand Set

The basic Wolff Wand set contains a 0.5" (13 mm) polished brass ball and a silvery nickel-plated ball on each end of the 10.5" (28 cm) steel rods.

. **657500**



Double-Ended Wolff Wand

Double-Ended Wolff Wand has a 0.5" (13 mm) polished brass ball on one end and a silvery nickel-plated ball on the opposite end of a 10.5" (28 cm) steel rod.

. **657600**

Wolff Wand Kit With Stand

The Wolff Wand Kit contains a 0.5" (13 mm) polished brass ball, a nickel-plated ball, and both a large and small Acrylic Fixation Wand. The Wolff Wand Stand is also included in this kit.

The rods are covered in comfortable black plastic coating material for durability.

Kit with Stand **660700**
 Stand Only **657900**



Wolff Wand Clear Storage Case

Clear Plastic Box (top and bottom) storage case for the 657500 Basic Wolff Wand Set. Wolff Wands are sold separately.

. **443200**

Examination Tools



Brock String

The Brock String is used during treatment of convergence insufficiency and other binocular vision anomalies. Used to develop convergence skills and disrupt suppression of one of the eyes. Available in varying lengths and in quantities of 1 or in a package of 12.

20' Brock String with 3 Beads

Each **603100**
 12-Pack **593700**

12' Brock String with 5 Beads

Each **603300**
 12-Pack **602400**

10' Brock String with 3 Beads

Each **603200**
 12-Pack **593800**

6' Brock String with 5 Beads

Each **603400**
 12-Pack **593900**



Visual Edge Slant Board

Designed to be lightweight, durable, transportable, multifunctional, and affordable. Angle is 22 degrees. A Dry Erase and Magnetic board surface for accomplishing direct "on board" tasks. Paper clip at the top holds work firmly to the board. A removable "book catch" at the bottom allows for ease of reading any book or for unobstructed writing. Folding legs lock into extended position for stability; folds flat for compact storage. Fits into 2-inch wide space. 12" x 18" (30.48 cm x 45.72 cm). Weighs 1.5 lbs. (.68 k).


. **300700**



Clear Pupilometer

Pupil sizes from 1.5 mm to 10 mm in 0.5 mm increments are easily measured. Can be used vertically or horizontally. The reverse side has a mm scale for measurement of PD. Ruler is printed on clear semi-rigid plastic.

. **566000**



Toucan Fixation PD Rule

Toucan Fixation and PD Rule with near optotypes, including LEA Symbols, LEA Numbers, and HOTV optotypes for 20/32, 20/63, 20/80, and 20/100. Includes inch/mm rules and diopter scale plus pupilometer in half mm increments. 6" (15.2 mm) long.

. **497000**



White to White Corneal Measure

This simple White to White gauge can be used to measure the horizontal cornea diameter for IOL and implantable contact lenses. Gauge provides comparison rules for corneal measurement in 0.5 mm increments from 8.5 mm to 14 mm. Application is similar to Holladay-Godwin Corneal Gauge.

As an option, screw fastener can be removed to use gauges separately. Clean with alcohol. Overall dimensions when closed are 3.5" x 1.4375" (8.9 cm x 3.7 cm).

. **599700**

Marsden Ball

A basic vision therapy tool no vision therapy practice should be without. A 3" hollow ball made of durable vinyl with applied black letters and numbers on a 10' cord. Improved Marsden Balls now have a better weight which offers a better swing.

. **108100**



Hertel Exophthalmometer

This instrument is used for precise exophthalmos measurements (e.g., exophthalmic goiter). The basis of the measurement for exophthalmos determination using the Hertel version is the outer orbital rim (orbital wall) and the apex of the cornea. The integral refraction mirrors with millimeter scales for the left-hand and right-hand measuring halves of the exophthalmometer are calibrated so that the zero mark on the scale will be located in the plane of the resting points of the supports. For exophthalmus measurements the supports are placed against the two temporal orbital walls so that the orbital rim contacts the deepest point of the supports.

. **521500**



Luedde Exophthalmometer


Fastest, easiest, and least expensive method for measuring exophthalmos. By sighting at right angles, this device eliminates parallax and allows the clinician to measure the amount of protrusion or bulging eye from its orbit.

Prisms are clear plastic with easy-to-read millimeter lines on both sides. Eye is viewed through lower half of scale and the scale is seen through the upper half, eliminating parallax and aiding rapid determination of protrusion.

The Luedde Exophthalmometer is made of lightweight durable acrylic. Includes quality light transmission properties and is imprinted with quality screened score marks. Edges are polished by hand to ensure safety. Length is 3" (7.6 cm).

Single Luedde **781000**


Double Luedde with Case **161030**



Hendrickson Lifesaver Card

Base-in Lifesaver vision therapy cards are used with patients who need to increase their ability to diverge without effort. These cards are often sent home with the patient for therapy training. Available in either white or clear/translucent card.


White..... **108500**
Clear..... **108600**



Eccentric Circles

Paired cards are useful in positive and negative fusional range development.

..... **108000**



Transparent Target Forms

Bullseye targets based on procedures developed by Lawrence W. MacDonald, OD and Baxter Swartwout, OD. Each comes in a pack of 8 translucent forms.

Form A is a set of 8 clear bullseye targets measuring 2" x 3.5"

Form B is a set of 8 clear targets with five lines of words appearing above a bullseye, each line decreasing in font size, and two lines of words below the bullseye. Each target measures 3" x 3.5"

Form C is a set of 8 clear targets with three lines of letters below the bullseye. Each target measures 1.5" x 3.5".

Form A..... **108200**
Form B..... **108300**
Form C..... **108400**



JLC Funduscope Grading Card

The Funduscope Grading Card is used to help grade various funduscope challenges. The primary tool will assist the provider to better determine

C/D ratios for glaucoma patients. The JLC Funduscope Grading Card can also be used to select proper A/V (artery/vein) and ALR (arteriolar light reflex) ratios. Card is 2-sided styrene approximately 4.25" x 6" (15.5 cm x 10.5 cm).

..... **490100**



Anti-Suppression Goggles

Red/Green Anti-Suppression Goggles are intended to be used with various games and exercises for anti-suppression training. The Red/Green Anti-Suppression Goggles can also be used in exercises for ocular motor difficulties, accommodation development and visual perception.

Large Size: 2.8125" x 1.9375" (7.2 cm x 5 cm)
Each..... **681900**
6 Pack..... **681906**
Small Size: 2.375" x 1.9375" (6 cm x 4.5 cm)
Each..... **682000**
6 Pack..... **682006**



Anti-Suppression Patches

Red or Green Anti-Suppression Goggles are intended to be used with various games and exercises in anti-suppression training. Sold in a package of two patches.

Large Size: 2.8125" x 1.9375" (7.2 cm x 5 cm)
Red..... **676202**
Green..... **682102**
Small Size: 2.375 x 1.6875" (6 cm x 4.5 cm)
Red..... **676302**
Green..... **682202**



Alger Brush

Used by ophthalmologists, ED physicians, and optometrists trained and licensed to remove foreign bodies and rust rings from patients' eyes. The rotary instrument uses a very low torque motor powered by a single AA battery (included) to "brush" rather than drill the rust ring leaving a smooth surface on the cornea or sclera. Brushing promotes faster health time.

Carbide burrs (for corneal rust ring removal) are available in 1 mm or .5 mm diameter (the .5 mm size is most commonly used). The diamond burrs for assistance in pterygectomies come in the following grits and configurations: 2.5 mm round fine grit, 3.5 mm round medium grit, 4.0 mm wheel fine grit and 5.0 mm wheel medium grit.

..... **222000**
..... **222001**
..... **222100**
..... **222005**
..... **222001**
..... **222100**
..... **222005**
..... **222001**
..... **222005**

Alger Brush II Replacement Burrs

Unsterilized - 5 per box
0.5 mm Burr..... **222002**
1 mm Burr..... **222003**
Pre-sterilized - 5 per package
0.5 mm Burr..... **222004**
1 mm Burr..... **222005**

Chuck and Burr with Carbide Tip

0.5 mm and 1 mm replacement chuck and burr.
0.5 mm..... **222006**
1 mm..... **222007**

Replacement Motor

..... **222008**

Replacement Chuck


..... **222009**

Replacement Housing

..... **222010**

Algerbrush Tip Removal Kit

..... **222300**



Small Eye Model

Oversized model is 5 times the size of the human eye and shows inner anatomy including optic nerve, disc, macula, retina, central retinal artery, and vein. Lens and cornea are removable. Size 5" x 3" x 4" (13 cm x 7.5 cm x 10 cm). Information card is 6.5" x 5.25" (16.5 cm x 13 cm).

..... 493900



Cataract Eye Model

Oversized model includes interchangeable lenses that show different types of cataract including subcapsular, capsular, mature, cortical, and nuclear. Card highlights visual effects of cataracts. Size 5" x 3" x 4" (12.7 cm x 7.5 cm x 10 cm).

..... 493800



Model of the Eye

Plastic molded eye is 6 times the size of the human eye. Disassembles into 7 parts including top, vitreous body, bottom section, 2-part

lens, iris, and cornea. Mounted on a stand. Detailed coloring. Size is 5" x 5" x 6" (12.7 cm x 12.7 cm x 15.2 cm).

..... 512500



Lasik Eye Model

The LASIK Eye Model is made of hard, durable plastic and measures 5.75" (14.6 cm) in diameter. It is realistically colored and textured, and is designed primarily for patient education in LASIK surgery.

The realistic "flap" in the front of the model simulates the incision by the surgeon during the LASIK surgery, whereby the flap pivots back and forth.

In addition to reassuring the patient, the model also provides the opportunity to demonstrate the surgical effect on the focal length of the eye, and how acuity is improved through surgery. Quick and easy-to-use, this model will save the doctor's or technician's time.

..... 151200



Intraocular Eye Model

The IOL Eye Model demonstrates how vision is regained after cataract removal. Six interchangeable lenses simulate IOL and toric implants of clear posterior capsule, cloudy capsule, and post YAG posterior capsulotomy.

The model also shows a cross section of the eye in 3 layers; the retina, choroid, and sclera with veins and arteries, including central retinal. Measures 3.75" (9.5 cm) in diameter, and is mounted on a white, opaque plastic base. It is realistically colored and textured and the cornea, optic nerve, and insertion of muscle are shown.

..... 151490



Corneal Eye Model

Includes 5 interchangeable corneas that have the following conditions: Bullous Keratopathy, Fuch's Endothelial Dystrophy, Keratoconus, Recurrent Corneal Erosion, and normal.

The model is made of hard, durable plastic and measures 3.75" (9.5 cm) inches in diameter. Mounted on a plastic base and the lenses are 2" (5 cm) in diameter. The model is realistically colored and textured and is designed primarily for patient education.

..... 151820



Astigmatism Eye Model

This model of the eye is for patient education in astigmatism, myopia, and hyperopia - or refractive error. The Astigmatism/Myopia/Hyperopia Model demonstrates how the focal length of the eye is changed by lens correction for better acuity. The two chords realistically simulate both a myopic or hyperopic foci.

..... 151230



Visit online at
www.good-lite.com

to see the complete line of Eye Models available from Good-Lite.

Slit Lamp and Lenses



All Glass Lenses and Made in the USA



LEGEND Series Slit Lamp Model YZ-5F1

Precisely engineered German-imported Schott optics reduce optical distortion, and increase clarity and depth of focus. Single joy-stick design easily controls both horizontal focus and vertical elevation. Tilt

mechanism for ease of angled examination. Eyepiece diopter adjustment is 12.5X (-5D ~ +3D). Contains five-step drum magnification at 6X, 10X, 16X, 25X, and 40X with the Sclera vein observable at 40X. Halogen light source for acute illumination.

All optical lenses are moisture-proof, mildew-proof, and anti-reflective.

..... 995300



LEGEND Series Slit Lamp Model YZ-5F

Precisely engineered German-imported Schott optics reduce optical distortion, and increase clarity and depth of focus. Single joy-stick design easily controls both horizontal focus and vertical elevation. Tilt mechanism for ease of

angled examination.

Eyepiece diopter adjustment to 10X (±8D) and 16x (±10D). Contains magnification at 10X (φ18 mm), 16X (φ14.5 mm), 16X (φ11.25 mm), and 25.6X (φ9 mm). Halogen light source for acute illumination.

All optical lenses are moisture-proof, mildew-proof, and anti-reflective.

..... 995400

Slit Lamp Table

..... 995310



ezView 20D Slit Lamp Lens

There's a 20D in every physician's pocket. The computer optimized lens surfaces eliminate optical aberrations.

Viewing Range: 46°-60°
Image Mag.: 3.08X
Laser Spot Mag.: .32X

..... 995020



ezView 28D Slit Lamp Lens

The 25D is ideal for use in patients with smaller palpebral fissures and offers median field of view and magnification.

Viewing Range: 53°-69°
Image Mag.: 2.27X
Laser Spot Mag.: .44X

..... 995028



ezView 78D Slit Lamp Lens

The ion VISION ezView 78D non-contact slit lamp lens is an ideal general diagnostic and treatment lens.

Viewing Range: 97°
Image Mag.: .93X
Laser Spot Mag.: 1.08x

..... 995078



ezView 90D Slit Lamp Lens

The ion VISION ezView 90D slit lamp lens delivers superior clarity, stereopsis and outstanding small pupil capabilities.

Viewing Range: 89°
Image Mag.: .76X
Laser Spot Mag.: 1.32X

..... 995090



ezView 125 Slit Lamp Lens

It is an excellent general use slit lamp lens balancing high magnification viewing, wide field of view and superior small pupil capability.

Viewing Range: 124°
Image Mag.: .57X
Laser Spot Mag.: 1.75X

..... 995125



OmniView 165 Slit Lamp Lens

The ONLY ALL GLASS Wide Field Indirect Laser Lens on the market!

Ultra-wide 165° field imaging to the ora serrata. The OmniView 165 provides a clear view of the retina from the posterior pole to the periphery without the need of shifting or tilting the lens. Perfect for PRP or any mid-to-wide field applications.

Viewing Range: 165°
Image Mag.: .5X
Laser Spot Mag.: 2.0X

..... 995165



DirectView 3 Mirror Slit Lamp Lens

ion VISION's 3 mirror glass lenses are of a traditional design, incorporating 3 mirrors angled at 60°, 66° and 78°. Magnification and laser spot factor is an exacting 1.0X.

The ALL GLASS contact element eliminates the need for messy fluid interface and features a "flanged" 14.5 mm contact surface that is ideal for laser treatment. No fluid design due to ALL GLASS contacting element - an ion VISION exclusive feature!

Viewing Range: 60°/66°/78°
Image Mag.: 1.0X
Laser Spot Mag.: 1.0X

..... 995103



Visit online at
www.good-lite.com

to see the complete line of ion lenses available for the slit lamp.



GOOD-LITE®

The Quality Always Shines Through

1155 Jansen Farm Drive

Elgin, IL 60123

Phone: 800-362-3860 Fax: 888-362-2576

Phone: 847-841-1145 Fax: 847-841-1149

www.good-lite.com



The Original and Still the Best Value in Vision Testing and Screening

Since 1930, Good-Lite™ has been a leading name in the medical supply industry. Good-Lite™ provides eye care professionals, primary care physicians, nurse practitioners, physician assistants, medical practice nurses, medical assistants, school health nurses, and Head Start and early childhood program personnel with the best value and quality in vision testing and screening equipment, eye charts, devices for instrument-based screening, and accessories.