Auto Kerato-Refractometer KR-800S





Auto Kerato-Refractometer with objective and subjective testing



5 IN 1 AUTO KERATO-REFRACTOMETER KR-800S

KR-800S AUTO KERATO-REFRACTOMETER

Topcon is proud to introduce the new KR-800S Auto Kerato-Refractometer. The KR-800S is unique because it features not only objective autorefraction and keratometry but it also performs subjective far and near testing as well as 3 function tests. These 5-in-1 Functions assure quick and accurate results and enhance your test workflow.





MORE THAN A SIMPLE AUTO KERATO-REFRACTOMETER

All measured data can be observed from the wide 8.5 inch color touch screen panel, allowing the user to quickly see each data point and explain the results to the patient. Moreover, Topcon's ability to engineer a weight reduction of approximately 23%, as compared to older Topcon autorefractors, as well as the new Auto-Vertical mode have contributed to a smooth control of the unit during the measuring process. The KR-800S is more than a simple Auto Kerato-Refractometer, it will perfectly match your needs.

OBJECTIVE AND SUBJECTIVE DATA

Both right and left eye information appear on a single 8.5inch wide color touch screen with all obtained data together, which makes it extremely easy to compare: **Objective (SCA), Subjective (SCA & ADD & VA), CL (SCA & ADD & VA), Glare/Grid/Contrast VA.** Utilizing this single display of all the data, the operator can easily understand the current refraction SCA of both eyes, whether the patient has presbyopia, as compared to the current SCA prescription. Moreover, the KR-800S can perform several function tests such as Glare/Grid/Contrast without the need to prepare any other special devices.



OBJECTIVE TEST



The unique Rotary Prism[™] Technology, exclusive to Topcon, allows for unparalleled precision and reliability. Quick measurement becomes available by decentering and rotating the measurement ring projected on the retina rapidly. Moreover, it decreases the influence of an uneven reflection on the eye or a cataract eye.







SUBJECTIVE VA TEST*



The test results of all objective and subjective measurements can be shown on the monitor. Therefore, it is very easy to compare the VA difference between the objective and subjective tests. If a computerized lensmeter is connected, it can also test and show the patient's VA result with their current eyeglasses. Since it is easy to compare VA with the patient's current eyeglasses and the BCVA result, if necessary, new eyeglasses can be introduced.



*Cylinder power and axis cannot be changed for the subjective test. Instead, refractometer data will be used.

*For a precise prescription of eyeglasses, we recommend that you perform the binocular test.

GLARE TEST



The Glare test is useful for screening cataract patients by understanding the patient's VA with/without Glare. The test can also be used to understand the patient's VA of pre & post LASIK surgery or to know the patient's QOV.

* The Glare test can only be performed with the subjective far distance test.



Glare Vision





Normal Vision



The Grid test can be performed for screening AMD and other macular related diseases. Patients with macular disease may see wavy lines or missing lines. The test can be helpful in detecting early signs of abnormality in the eye so a treatment plan can be started.

* The Grid test can only be performed with the subjective near distance test.



Sample image of patient with AMD







The Contrast test is an ideal test to determine the patient's QOV. The number of contrast steps can be changed accordingly; 2.5%, 5%, 10%, 12.5%, 25%, 50%, 100%.

* The Contrast test can only be performed with the subjective far distance test.

* Only the contrast of the background changes.

5 IN 1 AUTO KERATO-REFRACTOMETER KR-800S

PRINTOUT SAMPLE

	-KR 010001- OID : NAME
	2013_12_24 AM 10 : 00 N o . 0 0 0 1 01 SN : SBJ. DATA(REF)
Subjective Measurement	(R) S C A VA -5.00 - 2.00 75 0.9 (L) S C A VA -0.25 -1.00 90 1.2
¹ Subjective refraction Far VA value	NEAR TEST(REF) <r> DIST. ADD VA</r>
² Subjective refraction Near VA value	2 40 cm + 2.50 0.7 ⟨L⟩ DIST. ADD VA 33 cm +2.25 1.2
3 Grid test result	GRID CHART(REF) ⟨R⟩ ⟨L⟩ TS: NG NS: NG NS: 0K TS: 0K
4 Glare test result	C: NG C: 0K TI: 0K NI: 0K NI: 0K TI: NG
5 Contrast test result	$4 \qquad \qquad$
6 Far VA for lensmeter	5 CONTRAST TEST(REF) 5 <r> VA LVL.</r>
7 Near VA for lensmeter	0.8 50% 1.0 50% SBJ. DATA(CL) (R) S C A VA
⁸ Glare test VA for lensmeter	$\begin{array}{c ccccc} 6 & -2.00 & -1.00 & 95 & 0.6 \\ \langle L \rangle & S & C & A & VA \end{array}$
9 Contrast test VA for lensmeter	−0. 25 −1. 00 100 1. 2 NEAR TEST(CL) ⟨R⟩ DIST. ADD VA
	7 40 cm + 1.00 0.5 <l> DIST. ADD VA</l>
	33 cm +1.00 0.7 GLARE TEST(CL)

 $\langle L \rangle$

 $\langle L \rangle$

VA

VA 0.6

LVL. 25%

0.7

TOPCON

 $\langle R \rangle$

VA

 $\langle R \rangle$

VA

LVL.

9

0.3 CONTRAST TEST(CL)

0.5

50 %

Specifications

Refractive Power Measurement	
Spherical refractive power	-25D to +22D (0.12D/0.25D steps)*
Cylindrical refractive power	0D to ±10D (0.12D/0.25D steps)*
Astigmatic axial angle	0°to 180° (in 1°or 5°steps)
Minimal measurable pupil diameter	Ø 2 mm
	0 2 100
Corneal Curvature Measurement	
Corneal curvature radius	5.00 to 10.00mm (0.01mm step)
Corneal refractive power	67.50D to 33.75D (0.12D/0.25D steps) (where, corneal refractive power =1.3375)
Corneal astigmatic refractive power	OD to ±10D (0.12D/0.25 D steps)
Corneal astigmatic axial angle	0°to 180° (1°/ 5°steps)
Range of Subjective refractive check	
pherical refractive power:	Spherical refractive power: -18D to +18D (0.25D steps)
Test chart:	Eyesight test chart of 0.1 to 1.2 or 20/200 to 20/15, Grid display
Chart display:	Overall, Horizontal series, Contrast change
Test items:	Far-sightedness, Near-sightedness, Glare test
PD Measurement Range	20mm to 85mm (0.5mm step)
Data Transport Terminal	USB (Import) /RS-232C (Import/Export) / LAN (Export)
Dimensions	317-341mm (W) × 521-538mm (D) × 447-477 mm (H)
Weight	15 kg
Power Supply	100-240V AC, 50-60Hz, 70VA

*-25D \leq spherical refractive power + cylindrical refractive power or spherical refractive power + cylindrical refractive power \leq +22D

System chart



Topcon Europe Medical B.V. Essebaan 11; 2908 LJ Capelle a/d IJssel; PO. Box 145; 2900 AC Capelle a/d IJssel; The Netherlands Phone: +31-(0)10-4585077; Fax: +31-(0)10-4585045 E-mail: medical@topcon.eu; www.topcon-medical.eu

Topcon Danmark Praestemarksvej 25; 4000 Roskilde, Danmark Phone: +45-46-327500; Fax: +45-46-327555 E-mail: info@topcondanmark.dk www.topcondanmark.dk

 Topcon Scandinavia A.B.

 Neongatan 2; P.O. Box 25; 43151 Mölndal, Sweden

 Phone: +46-(0)31-7109200; Fax: +46-(0)31-7109249

 E-mail: medical@topcon.se; www.topcon.se

Topcon España S.A. HEAD OFFICE; Frederic Mompou, 4; 08960 Sant Just Desverr; Barcelona, Spain Phone: +34-93-4734057; Fax: +34-93-4733932 E-mail: medica@topcones; www.topcones

Topcon Italy Viale dell' Industria 60; 20037 Paderno Dugnano, (MI) Italy Phone: +39-02-9186671; Fax: +39-02-91081091 E-mail: info@topcon.it www.topcon.it

Topcon France BAT A1; 3 route de la révolte, 93206 Saint Denis Cedex Phone: +33-(0)1-49212323; Fax: +33-(0)1-49212324 E-mail: topcon@topcon.fr; www.topcon.fr Topcon Deutschland GmbH Hanns-Martin-Schleyer Strasse 41; D-47877 Willich, Germany Phone: (+49) 2154-885-0; Fax: (+49) 2154-885-177 E-mail: med@topcon.de; www.topcon.de

Topcon Polska Sp. z o.o. ul. Warszawska 23; 42-470 Siewierz; Poland Phone: +48-(0)32-670-50-45; Fax: +48-(0)32-671-34-05 www.topcon-polska.pl

Topcon (Great Britain) Ltd.

Topcon House; Kennet Side; Bone Lane; Newbury Berkshire RG14 5PX; United Kingdom Phone: +44-(0)1635-551120; Fax: +44-(0)1635-551170 E-mail: medical@topcon.co.uk, www.topcon.co.uk Topcon Ireland Unit 276, Blanchardstown; Corporate Park 2 Ballycoolin; Dublin 15, Ireland Phone: +353-18975900; Fax: +353-18293915 E-mail: medical@topcon.ie; www.topcon.ie



TOPCON CORPORATION

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, Japan. Phone: 3-3558-2523/2522, Fax: 3-3960-4214, www.topcon.co.jp