

## Reusable Bimanual Irrigation Aspiration System



EN ISO 13485 / Class IIa / CE 0408

### Product Features:

- Lightweight, colorcoded irrigation and aspiration Handpieces
- Compatible with your existing Phaco or I/A System
- Rough tips available for capsule polishing
- Swiss Design / Swiss Made


### Aspiration:




**BA 100/T**  
Aspiration Handpiece, round cannula Ø 0.80 mm, polished tip, Asp. port 0.30 mm



**BA 102/T**  
Aspiration Handpiece, round cannula Ø 0.80 mm, polished tip, Asp. port 0.40 mm



**BA 104/T**  
Aspiration Handpiece, round cannula Ø 0.80 mm, polished tip, Asp. port 0.50 mm for Viscoelastic



**BA 105/T**  
Aspiration Handpiece, round cannula Ø 0.80 mm, rough tip, Asp. port 0.30 mm

### Irrigation:



**BI 120/T**  
Irrigation Handpiece, round cannula Ø 0.80 mm, polished tip, dual 0.50 mm side ports



**BI 125/T**  
Irrigation Handpiece, round cannula Ø 0.80 mm, rough tip, dual 0.50 mm side ports



**BI 140/T**  
Irrigation Handpiece, round cannula Ø 0.80 mm, polished tip, open end Tip 45°



**BI 121/T**  
Irrigation Handpiece, round cannula Ø 0.90 mm, polished tip, dual 0.65 mm side ports high flow

./.



**BA 108/T**  
Aspiration Handpiece, round  
cannula, Ø 0.90 mm, rough  
tip, Asp. port 0.30 mm



**BI 126/T**  
Irrigation Handpiece, round  
cannula Ø 0.90 mm, polished  
tip, high flow open end Tip 45°



**BAO 105/T**  
Aspiration Handpiece, oval  
cannula 0.55 mm x 0.95 mm,  
mate finish, rough tip, Asp.  
port 0.30 mm



**BIO 120/T**  
Irrigation Handpiece, oval  
cannula 0.55 mm x 0.95 mm,  
mate finish, irrigation port 0.55  
mm



**BAO 102/T**  
Aspiration Handpiece, oval  
cannula 0.70 mm x 1.10 mm,  
mate finish, Asp. port 0.35 mm



**BIO 127/T**  
Irrigation Handpiece, oval  
cannula 0.70 mm x 1.10 mm,  
mate finish, high flow open  
end Tip



**BI 130**  
**Irrigation Handle**  
with Luer adapter for  
irrigation cannulas