


A legacy of innovation and quality

Visitec, a market leader of ophthalmic cannulas, offers a full range of cannulas for several ophthalmic procedures, including cataract, refractive, vitreo-retinal, and oculoplastics. Our cannulas are produced at manufacturing locations in the United States and the United Kingdom where a high level of consistency and product quality is upheld.

Visitec has been focused on providing innovative cannula designs with the highest quality since 1979. A leading brand of Beaver-Visitec International (BVI). Visitec has been driving the ophthalmic market towards single-use cannulas thereby improving the safety of your patients.

Single-use cannulas are highly recommended for surgical use by the ASCRS (American Society for Cataract & Refractive Surgeons) and ESCRS (European Society for Cataract & Refractive Surgeons) because instruments with a small lumen such as a cannula or tubing can be extremely difficult to fully clean. When a cannula is not cleaned and decontaminated effectively, this can lead to the occurrence of Toxic Anterior Segment Syndrome (TASS) and other potential infections. TASS is an acute post-operative inflammatory reaction in which a non-infectious substance enters the anterior segment and induces toxic damage to the intraocular tissue. The utilisation of single-use cannula may help prevent the occurrence of TASS, thereby reducing the risk to the patient.




References:
<http://ascrs.org/sites/default/files/TASS%20Guidelines.pdf>
http://www.es CRS.org/vienna2011/programme/handouts/IC-100/IC-100_Barry_Handout.pdf
<http://emedicine.medscape.com/article/1190343-overview>

NEW Cannulas from BVI



As ophthalmic surgical techniques continue to evolve at a fast pace, BVI is right there to provide new single-use cannulas to meet your diverse needs.

Now with even more options, below are the newest additions to our high-quality line of cannulas.



For ordering information call: +44 (0)1865 601256 or visit www.beaver-visitec.com



Beaver-Visitec International, Sales Limited
 85c Park Drive, Milton Park,
 Abingdon, Oxfordshire, OX14 4RY, UK
 Tel: +44 (0)1865 601256 Fax: +44 (0)1865 595761
www.beaver-visitec.com

BVI, BVI Logo and all other trademarks are property of Beaver-Visitec International (BVI) © 2015 BVI

Visitec® Cannulas

A legacy of innovation and quality








- A market leader of ophthalmic cannulas for over 35 years
- Full range of cannulas encompassing all ophthalmic surgical procedures
- A history of innovation to meet your ever-changing surgical needs
- Manufactured in the United States and United Kingdom for quality assurance you can rely on
- Focused on single-use cannulas to improve patient safety















590361-0715



















A sampling of cannulas from our product catalogue





















Anesthesia	Description	Item #
Retrobulbar (Sharp Point) 10 per box.	  .50 x 38 mm (25G x 1½ in)	581275
Retrobulbar/Peribulbar [Atkinson] 10 per box.	.40 x 32 mm (27G x 1¼ in)	585109
	.50 x 32 mm (25G x 1¼ in)	585108
	.60 x 32 mm (23G x 1¼ in)	585001
Retrobulbar [Atkinson] 10 per box.	.50 x 38 mm (25G x 1½ in)	581637
	.60 x 38 mm (23G x 1½ in)	581638
Sub-Tenon's Anesthesia Cannula 10 per box.	  1.10 x 25 mm (19G x 1 in)	585176







Irrigation	Description	Item #
Anterior Chamber Cannula Angled 45 degrees, 4 mm from end. 10 per box.	  .30 x 22 mm (30G x 7⁄8 in)	581273
	.40 x 22 mm (27G x 7⁄8 in)	581280
	.50 x 22 mm (25G x 7⁄8 in)	581639
	.60 x 22 mm (23G x 7⁄8 in)	585025
	.80 x 22 mm (21G x 7⁄8 in)	581281
Anterior Chamber Cannula Angled 45 degrees, 6 mm from end. 10 per box.	  .30 x 22 mm (30G x 7⁄8 in)	581285 NEW
	.40 x 22 mm (27G x 7⁄8 in)	581284 NEW
	.50 x 22 mm (25G x 7⁄8 in)	581642 NEW
Viscoflow[™] Cannula Angled 45 degrees, 9 mm from end. 10 per box.	  .30 x 22 mm (30G x 7⁄8 in)	585046 NEW
	.40 x 22 mm (27G x 7⁄8 in)	585026 NEW
	.50 x 22 mm (25G x 7⁄8 in)	585045 NEW
.60 x 22 mm (23G x 7⁄8 in)	585009 NEW	









Capsule Polishers	Description	Item #
Flexible Tip Cannula Angled 45 degrees, 7 mm from end, 10 mm sleeve. 10 per box.	  .40 x 22 mm (27G x 7⁄8 in)	585005
	.50 x 22 mm (25G x 7⁄8 in)	585017 NEW
Olive Tip Polisher .33 mm end port, 30 degrees, 10 mm from end. 10 per box.	  .60 x 25 mm (23G x 1 in)	585143
Capsule Polisher [Kratz] .30 mm top port, 35 degrees, 9 mm from end. 10 per box.	  .60 x 22 mm (23G x 7⁄8 in)	581605

Hydrodissection	Description	Item #
Nucleus Hydrodissector [Pearce] Angled 35 degrees, 8 mm from end. 10 per box.	  .40 x 22 mm (27G x 7⁄8 in)	585036 NEW
	.50 x 22 mm (25G x 7⁄8 in)	585037 NEW
Nucleus Hydrodissector [Sauter] Tip is flattened on the vertical plane. 10 per box.	  .40 x 22 mm (27G x 7⁄8 in)	585099
Cortical Cleaving Hydrodissector (Flat End) 45 degrees, 11 mm from end. 10 per box.	  .40 x 22 mm (27G x 7⁄8 in)	585158
Nucleus Hydrodissector Tip is flattened on the horizontal plane. 10 per box.	  .40 x 22 mm (27G x 7⁄8 in)	585155
Micro J Shape Cannula .25 mm open end port. 10 per box.	  .40 x 22 mm (27G x 7⁄8 in)	585018 NEW
	.50 x 22 mm (25G x 7⁄8 in)	585003 NEW
Nucleus Hydrodissector [Chang] 90 degree bend, 1 mm flattened tip. 10 per box.	  .40 x 18 mm (27G x 7⁄8 in) bevel angle 33 degrees.	585294
	.40 x 18 mm (27G x 7⁄8 in) bevel angle 47 degrees.	585295
	  .40 x 22 mm (27G x 7⁄8 in)	585198 NEW
Nucleus Hydrodissector [Akahoshi Style] Tapered tip, thin wall. 10 per box.	  .40 x 22 mm (27G x 7⁄8 in)	585198 NEW
Nucleus Hydrodissector [Buratto] Formed, 1mm bent tip. 10 per box.	  .50 x 22 mm (25G x 7/8 in)	585199 NEW

Cystotomes	Description	Item #
Irrigating Cystotome Formed cystotomes are bent to conform with the convexity of the crystalline lens. 10 per box.	  .40 x 16 mm (27G x 5⁄8 in)	581618
	.50 x 16 mm (25G x 5⁄8 in) [Pearce]	581610
Irrigating Cystotome Formed 8.5 mm from bend to tip. 10 per box.	  .40 x 16 mm (27G x 5⁄8 in)	585008
	.50 x 16 mm (25G x 5⁄8 in)	585004
Irrigating Cystotome [Blumenthal] Formed 8.5 mm from bend to tip. 10 per box.	  .40 x 16 mm (27G x 5⁄8 in)	585042
Irrigating Cystotome [Simcoe] 45 degrees, 10.5 mm from tip. 10 per box.	  .40 x 12.5 mm (27G x 1/2 in)	585093
Irrigating Cystotome [Joyce] Tip .50 mm long. 10 per box.	  .40 x 16 mm (27G x 5⁄8 in)	585039
Irrigating Cystotome 10 per box.	  .50 x 16 mm (25G x 5⁄8 in)	581617

Irrigation & Aspiration	Description	Item #
I and A Cannula [Simcoe] Thin wall cannula. 5 per box	  .60 mm/.60 mm (23G/23G)	585047
Bimanual Aspiration Cannula Curved 45 degrees, .40 mm top aspiration port. 10 per box.	  .60 x 17 mm (23G x 5⁄8 in)	585222
Bimanual Irrigation Cannula Thin wall cannula, curved 45 degrees, .51 mm open end irrigation port. 10 per box.	  .80 x 17 mm (21G x 5⁄8 in)	585223

Lacrimal	Description	Item #
Lacrimal Cannula Curved, 28 mm length. 10 per box.	  .45 x 28 mm (26G x 1½ in)	581276
Lacrimal Cannula Curved, 35 mm length. 10 per box.	  .80 x 35 mm (21G x 1 3⁄8 in)	581615
Lacrimal Cannula Straight, 13 mm length. 10 per box.	  .50 x 13 mm (25G x 1/2 in)	585068 NEW

Refractive	Description	Item #
LASIK Irrigating Cannula 25 degrees, formed 7 mm from end. 10 per box.	  .40 x 22 mm (27G x 7⁄8 in)	585278
LASIK Irrigating Cannula [Solomon] Thin wall cannula. Angled tip features three .25 mm ports on 8 mm extension. 10 per box.	  .50 x 22 mm (25G x 7⁄8 in)	585280
LASIK Irrigating Cannula [Lindstrom] Features three .25 mm ports located on 7 mm extension. 10 per box.	  .50 x 22 mm (25G x 7⁄8 in)	585250
LASIK Flap Spatula/Irrigator [Hersh] 45 degrees, 9 mm from spatulated tip. 10 per box.	  .60 x 22 mm (23G x 7⁄8 in)	585252

For a complete list of our cannulas, please refer to our full product catalogue.