


# Single-Use Instruments





# Moria Quality

*For nearly 200 years, Moria has been recognized for its know-how in the manufacture of ophthalmic micro-surgery instruments.*

***We have put all of our expertise toward the development of a range of sterile disposable instruments, designed and manufactured in France.***


*Our manufacturing processes and our products comply with ISO 13485, FDA (Food And Drug Administration) standards. All our medical devices are in accordance with the European Directive 93/42/EEC.*






# Table of contents

	» SPECULUMS	6
	» FORCEPS	10
	» HOOKS AND RINGS	16
	» TORIC IOL POSITIONING	21
	» CANNULAS AND MICRO-SPONGES	25
	» NEONATOLOGY	28
	» CONSULTATION	32
	» EASYSET	38
	» CORNEAL TRANSPLANT	42
	Instruments	42
	Punches	44
	Trephines	46
	Artificial Chambers	47
	Microkeratome accessories	48








# SPECULUMS

## ▶ Neonatal Wire Speculum


6 and 9 mm blades

Box of 10 units



17258X10






# SPECULUMS

## ▶ Kratz Wire Speculum

15 mm open blades, 1 mm wire

Box of 10 units




17227X10

## ▶ Adjustable Speculum

15 mm blades

Can be used in nasal and temporal sides

Box of 10 units




17508X10

## ▶ Adjustable Speculum


Solid, 15 mm blades

Can be used with a femtocataract laser

Box of 10 units



17231X10







# FORCEPS

## ▶ Capsulorhexis Forceps - Straight, 90° pointed tips

1.8 mm incision

Box of 10 units

one steel



17502X10

## ▶ Capsulorhexis Forceps - Curved, 90° pointed tips

1.8 mm incision

Box of 10 units

one steel



17503X10

## ▶ Capsulorhexis Forceps - Curved, 90° pointed tips

1.8 mm incision

Box of 10 units

one



17237X10




# FORCEPS

## ▶ Capsulorhexis Forceps - Curved, 45° pointed tips

1.8 mm incision

Box of 10 units



one®




17252X10

## ▶ Bonn Forceps - Straight

0.12 mm teeth, 5 mm platforms

Box of 10 units



one® steel




17504X10

## ▶ Bonn Forceps - Angled

0.12 mm teeth, 5 mm platforms

Box of 10 units



one® steel



17505X10




# FORCEPS

## ▶ Bonn Forceps - Straight

0.12 mm teeth, 3.7 mm platforms

Box of 10 units



one®




17221X10


## ▶ Tying Forceps - Angled, Mc Pherson type

8 mm platforms

Box of 10 units



one®




17233X10

## ▶ Tying Forceps - Angled, Mc Pherson type

5 mm platforms

Box of 10 units



one® steel



17500X10



# FORCEPS

## ▶ Tying Forceps - Curved

5 mm platforms

Box of 10 units

one steel




17501X10

## ▶ Tying Forceps - Curved

4 mm platforms

Box of 10 units

one



17225X10

3.50 mm platforms

Box of 10 units

17247X10

## ▶ Tying Forceps - Straight


4 mm platforms

Box of 10 units

one



17222X10








# HOOKS AND RINGS

## ▶ Sinskey IOL Manipulator

Angled, 0.2 mm tip

Box of 10 units

one®



17223X10

## ▶ Lester IOL Manipulator

Angled, collar button tip

Box of 10 units

one®



17228X10

## ▶ Lens Manipulator

Conical, straight, blunt

Box of 10 units

one®




17230X10

# HOOKS AND RINGS

## ▶ Lens Manipulator phaco chopper

1.5 mm tip, inner cutting edge

Box of 10 units



one®




17229X10


## ▶ Drysdale Lens Manipulator

Blunt paddle design

Box of 10 units



one®




17242X10

## ▶ Thornton-Fine Fixation Ring

Blunt teeth, 1.3 mm internal diameter


Box of 10 units



one®




17232X10









# TORIC IOL POSITIONING

## ▶ Bubble Level

Preoperative marking of the 0-180° axis  
Atraumatic and thin marking surface:  
0.4 mm x 3 mm

Box of 10 units



one®

17249X10

## ▶ Mendez Ring

- 12 mm inside diameter
- Indicators at increments of 5°
- Four notches on the internal edge at cardinal positions
- Handle positioned at 0° to hold temporarily
- Matte surface for maximum visibility under the microscope

Box of 10 units




one®

17250X10

## ▶ Bores Marker


- Used inside the Mendez Ring to make two radial marks along the desired axis for precise alignment
- Precise 0.4 mm fine markings
- Blade design follows curvature of the cornea
- Matte surface for maximum visibility under the microscope

Box of 10 units



one®

17251X10







» CANNULAS AND MICRO-SPONGES

---

25





# CANNULAS AND MICRO-SPONGES

## ▶ Retrobulbar Needle - 25 G (0.50 mm)

35 mm length, sterile

Box of 10 units

one®




143BIS

## ▶ Straight Cannula, for children and adults - 25 G (0.50 mm)

Non traumatic and slightly blunt,  
12 mm long, outer diameter 25 G,  
inner 31 G, sterile

Box of 30 units



one®



13300

## ▶ Micro-sponges

Absorbent pads with sleeves,  
sterile, blunt tip

Box of 20 bags of 10 units

Distributed by:

Moria



30601

## ▶ Sponcals

Absorbent pads, non-sterile,  
cuneiform

Box of 100 units

Distributed by:

Moria



30761





» NEONATOLOGY

---

28



# NEONATOLOGY


## ▶ ROP Kit (Speculum and double-ended Depressor)

Neonatal speculum and  
double-ended depressor in  
sterile pack

Box of 20 units



one®



17210


## ▶ ROP Kit (Speculum and loop Depressor)

Neonatal speculum and loop  
depressor in sterile pack

Box of 20 units



one®



17255X20



# NEONATOLOGY

## ▶ ROP Kit (Speculum and double-ended Depressor)

Neonatal speculum and loop depressor in sterile pack

Box of 10 units

one®



17257X10

## ▶ Double-ended Depressor

Box of 20 units

one®




17253X20

## ▶ Neonatal Wire Speculum

6 and 9 mm blades

Box of 10 units



17258X10







» CONSULTATION

---


32

## ▶ Tying Forceps - Angled, Mc Pherson type

5 mm platforms

Box of 10 units

one steel





17500X10

## ▶ Tying Forceps - Angled, Mc Pherson type

8 mm platforms

Box of 10 units

one



17233X10

## ▶ Tying Forceps - Curved

5 mm platforms

Recommended for removing  
10-0 thread or eyelashes

Box of 10 units

one steel



17501X10





# CONSULTATION

## ▶ Tying Forceps - Curved

3.5 mm platforms

Recommended for removing  
10-0 thread or eyelashes

Box of 10 units



17247X10



one®

## ▶ Tying Forceps - Curved

4 mm platforms

Recommended for removing  
10-0 thread or eyelashes

Box of 10 units




17225X10

one®

## ▶ Tying Forceps - Straight

4 mm platforms

Box of 10 units



17222X10


one®

## ▶ Bonn Forceps - Straight

0.12 mm teeth, 5 mm platforms

Recommended for removing 6/0 or 7/0 thread, grasping compresses or filling cavities with compresses

Box of 10 units



one steel




17504X10

## ▶ Bonn Forceps - Straight

0.12 mm teeth, 3.7 mm platforms

Recommended for removing 6/0 or 7/0 thread, grasping compresses or filling cavities with compresses

Box of 10 units



one




17221X10

## ▶ Bonn Forceps - Angled

0.12 mm teeth, 5 mm platforms

Box of 10 units



one steel



17505X10



# CONSULTATION

## ▶ Scissors - Sharp

Recommended for cutting operative field or removing 6/0 to 10/0 thread

Box of 5 units

one steel



17506X5

## ▶ Scissors - Blunt

Recommended for cutting operative field or removing 6/0 to 10/0 thread

Box of 5 units

one steel




17507X5

## ▶ Straight Cannula, for children and adults - 25 G (0.50 mm)


Non traumatic and slightly blunt, 12 mm long, outer diameter 25 G, inner 31 G, sterile

Box of 30 units

one



13300






» EASYSET

---

38



# EASYSET


SPECULUMS	17227	one®
	17508	one® steel
CAPSULORHEXIS FORCEPS	17237	one®
	17252	one®
	17503	one® steel
	17502	one® steel
BONN FORCEPS	17221	one®
	17504	one® steel
	17505	one® steel
TYING FORCEPS	17222	one®
	17225	one®
	17233	one®
	17247	one®
	17501	one® steel
	17500	one® steel
IOL MANIPULATORS	17223	one®
	17228	one®
LENS MANIPULATORS	17229	one®
	17230	one®
	17242	one®
SCISSORS	17506	one® steel
	17507	one® steel

Cataract sets delivered in sterile blister packs containing a range of 3 to 7 instruments.

Easy Set packs are available in customized sets of both one® and one® steel instruments, containing up to:

- 3 Forceps
- 2 Manipulators
- 1 Speculum
- 1 pr. Scissors





# EASYSET



## Easy

- One-finger blister pack opening
- Once the blister pack is opened, the instruments are immediately ready for use

## Fast

- Fast preparation of the instruments

## Convenient

- Compact and ergonomic tray
- Instruments held in their dedicated locations





## » CORNEAL TRANSPLANT

---

<b>Instruments</b>	<b>42</b>
<b>Punches</b>	<b>44</b>
<b>Trephines</b>	<b>46</b>
<b>Artificial Chambers</b>	<b>47</b>
<b>Microkeratome accessories</b>	<b>48</b>




# CORNEAL TRANSPLANT Instruments

## ▶ Price Hook

Reverse Sinskey

Recommended in DMEK, DSAEK for  
descemetorhexis

Box of 5 units



one®



17302X5

## ▶ 90° Spatula

Recommended in DMEK, DSAEK for  
Descemet membrane peeling

Box of 5 units



one®



17303X5




# CORNEAL TRANSPLANT Instruments

one®

## ▶ Busin Spatula

Recommended in DSAEK for  
insertion of the donor graft,  
incision 5.5 mm

Box of 5 units



17300X5

one®

## ▶ Busin Forceps 23G

Recommended in DSAEK for  
insertion of the donor graft

Box of 1 or 5 units



17301X1

17301X5



# CORNEAL TRANSPLANT Punches

## ► Busin Punch

With aspiration  
Straight-walled cut (symmetrical double-bevel blade) throughout the thickness of the graft

Recommended for DMEK, DSAEK, DALK, PK


Available in:

- 6.00 mm
- 6.50 – 9.00 mm (increment of 0.25 mm)
- 9.50 mm

Box of 1 unit



one®



Blade shape

17200DXXX

## ► Muraine Punch

With aspiration  
Straight-walled cut (symmetrical double-bevel blade)  
2 opposed hinges


Recommended for DMEK

Available in 8.00 mm

Box of 1 unit



one®



Blade shape

17209D800



# CORNEAL TRANSPLANT Punches

## ▶ Guarded Punch

With aspiration  
Straight-walled cut (symmetrical double-  
bevel blade)

Recommended for DMEK

Available in:  
7.50, 7.75, 8, 8.50, 9.50, 10 mm

Box of 1 unit

one®



17207DXXX



# CORNEAL TRANSPLANT Trepines

## ▶ Adjustable Trepine

For donor and recipient  
Pre-adjustment of the desired depth,

Straight-walled cut: symmetrical double-bevel blade  
Centering guide in contact with the cornea  
Each turn corresponds to a depth of 250  $\mu\text{m}$   
Can be used with artificial chamber #17204

Available in:  
6, 6.50, 6.75, 7, 7.25, 7.50, 7.75, 8, 8.25, 8.50, 8.75, 9, 9.50,  
10 mm

Box of 1 unit

one®



17202DXXX

## ▶ Simple Trepine


For donor and recipient

Straight-walled cut: symmetrical double-bevel blade  
Centering guide in contact with the cornea  
Each turn corresponds to a depth of 250  $\mu\text{m}$   
Can be used with artificial chamber #17204

Available in:  
7, 7.25, 7.50, 7.75, 8, 8.25, 8.50, 8.75, 9, 9.50 mm

Box of 1 unit

one®



17201DXXX





# CORNEAL TRANSPLANT

## Artificial Chambers

### ▶ Artificial Chamber

To be used with with trephines  
#17201DXXX and #17202DXXX (as a  
replacement for punch #17200DXXX)

Box of 1 unit

one®




17204

### ▶ Artificial Chamber

To be used with microkeratomes One  
Use-Plus, One Use and CBm DSAEK #19175,  
#19155 and #19303 and Muraine Punch  
#17209D800

Box of 1 unit

one®



19182



# CORNEAL TRANSPLANT


## Microkeratome accessories

### ▶ Guide-Ring for One Use-Plus

Guide-Ring for motor One Use-Plus  
#19175

one®

Box of 5 units




19187

### ▶ Guide-Ring for One Use

Guide-Ring for turbine One Use  
#19155

one®

Box of 5 units




19186

### ▶ Guide-Ring for CBm

Guide-Ring for turbine CBm  
#19303

one®

Box of 5 units



19180



# CORNEAL TRANSPLANT

## Microkeratome accessories


### ▶ Heads One Use LC for One Use-Plus and One Use

Heads for motor One Use-Plus #19175 and turbine One Use #19155

Available in: 110, 130, 150, 200, 250, 300, 350, 400, 450, 500, 550 microns

Box of 5 units

one®



19184/XXX


### ▶ Heads CBSU for CBm

Heads for turbine CBm #19303

Available in: 50, 90, 110, 130, 200, 250, 300, 350 microns

Box of 5 units

one®



19178/XXX



# CORNEAL TRANSPLANT


## Microkeratome accessories

### ▶ Tubing Sets

Tubing Set for Artificial Chamber

Box of 10 units

one®



19181

Tubing Set for ACP

Box of 10 units




19192









### Manufacturer

#### MORIA SA

15, rue Georges Besse  
92160 Antony  
FRANCE

Phone: +33 (0) 1 46 74 46 74  
Fax: +33 (0) 1 46 74 46 70

[www.moria-surgical.com](http://www.moria-surgical.com)

#### MORIA Inc

1050 Cross Keys Drive  
Doylestown, PA 18902  
USA

Phone: (800) 441 1314  
Fax: + 1 (215) 230 7670

[www.moria-surgical.com](http://www.moria-surgical.com)

**Moria**  
Ophthalmic Instruments

#### MORIA JAPAN KK

Inami Bldg. 6F  
3-24-2 Hongo Bunkyo-ku  
Tokyo 113-0033  
JAPAN

Phone: + 81 3 5842 6761  
Fax: + 81 3 5842 6762

[www.moriajapan.com](http://www.moriajapan.com)

#### Moria Shanghai Office

RM4, 10 FL., Kai Li Building  
N°432 Huai Hai West Road  
200052 Shanghai  
CHINA

Phone: +86 21 52585278  
Fax: +86 21 52585278

[www.moria-surgical.com](http://www.moria-surgical.com)