



**Medical Workshop**

P.O.Box 461  
9700 AL GRONINGEN  
The Netherlands  
Phone +31 50 527 69 99  
Fax +31 50 527 69 58  
E-mail: [info@medicalworkshop.nl](mailto:info@medicalworkshop.nl)  
[www.medicalworkshop.nl](http://www.medicalworkshop.nl)

**ErgAdap**

Ergonomic 20° microscope adapter  
for slitlamps



MedicalWorkshop

## Medical Workshop **ErgAdap** 20° slitlamp adapters

The Medical Workshop **ErgAdap** eyepiece adapters incline the eyepieces 20° upward, thus assuring a more ergonomic working position for the researcher.

The adapters are constructed in such a way that the height of the eyepieces remains almost the same.



Topcon slitlamp without **ErgAdap**



Topcon slitlamp with **ErgAdap**

The **ErgAdap** fits almost any slitlamp with a 'dove-tail' coupling between the eyepieces and the magnification changer. Three types are available:

### **MW-8000 ErgAdap Black for:**

- Topcon SL-D7 slitlamp
- Topcon SL-D8Z slitlamp (\*)
- Topcon SL-7F slitlamp (\*)
- Topcon SL-8Z slitlamp (\*)
- Topcon SL-D4 slitlamp (\*)
- Topcon SL-D2 slitlamp (\*)
- Nidek SL-1600 slitlamp
- CSO slitlamps
- Zeiss 30 SL(M) slitlamp (\*)
- Zeiss 20 SL slitlamp
- Coherent LDS-10 laser slitlamp
- Coherent LDS-20 laser slitlamp
- Lumenis laser slitlamp

### **MW-8130 ErgAdap for:**

- Zeiss SL-130 (\*)

### **MW-8150 ErgAdap for:**

- Rodenstock RO-3000 (\*)
- Rodenstock RO-4000 (\*)
- Rodenstock RO-5000 (\*)

### **MW-8170 ErgAdap for:**

- Ellex slitlamp (\*)
- Reichert slitlamp (\*)



Rodenstock slitlamp with **ErgAdap**

(\*) These combinations are tested and approved